# SALIVA ELISAS

## Speichel Analytik

# DRG

Testosterone  
Progesterone  
Estradiol  
17alpha-OH Progesterone  
DHEA  
Estriol  
Cortisol  
Lactate


# DRG

# DRG ELISAS

## SALIVA Immunoassays Speichel Analytik

**Sensitive ELISAS**  
for detection of free steroids  
and other parameters in saliva  
samples

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### Benefits:

- Non invasive sampling, can be done also by patients
- Measurement of the free (unbound) fraction of steroids
- Much higher diagnostic relevance with reproducible and reliable results, compared to serum analytics
- Simple and patient-friendly measurement of hormone profiles (Diurnal-, monthly **profiles**)
- Easy sampling in babies and kids
- The measurement of steroids in **veterinary** saliva samples is possible without any further pre-treatment

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- Cat. No.: SLV-4188
- Sensitivity: 0.4 pg/mL
  - Standard range: 1 - 100 pg/mL

## Estrone

- Cat. No.: SLV-5228
- High Sensitivity
  - Short incubation time

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### Sampling device

As saliva contains only tiny amounts of protein non-polar analytes such as steroids tend to adsorb to various plastic surfaces. Glass would tend to be the best material for saliva sampling and for the most part Polyethylene, although it is still in use in many laboratories worldwide, should be avoided. However there is an inexpensive polypropylene device available from DRG (Sali-Tubes). It does not show any interference in saliva even in Progesterone testing, which is the most critical steroid.

### Blood in saliva

Traces of blood might interfere with the results of saliva testing. It is therefore important that the salivary samples should not show any red colour. In cases of visible blood traces in saliva the sample should be discarded. After rinsing the device several times with tap water a new saliva sample can be taken but only after waiting for approximately 15 minutes.