Instructions for Use

17-OH Progesterone ELISA

IVD

CE

REF EIA-1292

Σ 96

DRG

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.

Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées / Modificações introduzidas

The following changes have been made in comparison to the previous version:
 Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen:
 Rispetto alla versione precedente, sono state apportate le seguenti modifiche:
 Se han introducido los siguientes cambios en comparación con la versión anterior:
 Les modifications suivantes ont été apportées par rapport à la version précédente :

Detailed editorial revision. Changed wording in several chapters.
Ausführliche redaktionelle Überarbeitung. Geänderter Wortlaut in mehreren Kapiteln.
Revisione editoriale dettagliata. Modificato il testo in diversi capitoli.
Revisión editorial detallada. Se ha cambiado la redacción de algunos capítulos.
Révision éditoriale détaillée. Modification de la formulation dans plusieurs chapitres.

2 SCIENTIFIC VALIDITY REPORT	Addition of new information about screening incidence
8 REFERENCE VALUES	Improved acceptance range for the children by measuring the additional samples.
10 PERFORMANCE CHARACTERISTICS	Included reference to English version for all other languages.
13 LITERATURE	New references #13,14,15 added

ENGLISH	2
DEUTSCH	13
ITALIANO	21
ESPAÑOL	28
FRANÇAIS	36
LITERATURE	43
SYMBOLS USED	44

1 INTENDED USE

The DRG 17-OH Progesterone ELISA is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of 17- α -OH Progesterone (17- α -OHP) in serum or plasma (K₂-EDTA, K₃- EDTA, lithium heparin or citrate plasma 3.2 %).

For in vitro diagnostic use only. For laboratory professional use.

2 SCIENTIFIC VALIDITY REPORT

The steroid 17- α -Hydroxyprogesterone (17- α -OHP) is produced by both the adrenal cortex and gonads. Even though 17- α -OHP has relatively low pregestational activity, it is of intense clinical interest because it is the immediate precursor to 11-desoxycortisol (Cpd-S). Because Cpd-S is produced by 21-hydroxylation of 17- α -OHP, measurement of 17- α -OHP is a useful indirect indicator of 21-hydroxylase activity. In congenital 21-hydroxylase deficiency, the most common variety of congenital adrenal hyperplasia (CAH), 17- α -OHP is secreted in abundant excess. It is moderately elevated in the 11- β -hydroxylase deficiency as well. Measurement of 17- α -OHP is therefore valuable in the initial diagnosis of CAH.

2.1 Clinical Physiology

Adult non-pregnant women:

In adult non-pregnant women in the childbearing age group, 17- α -OHP concentrations vary over the menstrual cycle with luteal phase concentrations being higher than follicular phase concentrations. This is because 17- α -OHP is secreted parallel with progesterone from maturing follicles or from the corpus luteum. There is also a diurnal variation of 17- α -OHP concentrations.

This rhythm is parallel with adrenal cortisol secretion such that maximum 17- α -OHP concentrations are measured in samples obtained between midnight and 8:00 am.

Adult males:

There is little information available on the systematic variability of 17- α -OHP concentration in adult males.

Pregnant women and newborn children:

The steroid 17- α -OHP is produced in large amounts by the fetus and the adrenals. It is secreted in abundance into both the fetal and maternal circulation. The maternal concentrations of 17- α -OHP increase very sharply after 32 weeks gestational age to about 4-fold above basal concentrations at term.

2.2 Clinical Applications

Congenital adrenal hyperplasia:

The principal application of the 17- α -OHP ELISA is in the diagnosis of CAH in newborns with ambiguous genitalia and in virilized adolescent girls. Since 17- α -OHP is the immediate precursor to 11-desoxycortisol, basal 17- α -OHP concentrations are sharply elevated in patients with 21-hydroxylase deficiency and to a lesser degree in patients with 11-hydroxylase deficiency.

Because 17- α -OHP concentrations are so markedly elevated in newborns and adolescent girls afflicted with CAH, a single basal measurement is all that is normally required to make the diagnosis.

The reported screening incidence for classical CAH in Europe generally varies from 1:10,000 to 1:14,000, similar to the reported incidence results in North America of 1:15,000 to 1:16,000. Higher incidences were noted in Brazil (1:7,500), La Reunion (1:4,000), and Alaska (1:288 in Yupik Eskimos, and 1:800 in Native Alaskans). (13,14). A very high concentration of 17-hydroxyprogesterone at 3 days in full-term infant is diagnostic of classic 21-hydroxylase deficiency (15). Typically, salt-losing patients have higher 17-hydroxyprogesterone concentrations than non-saltlosers. False-positive results from neonatal screening are common with premature infants, and many screening programmes have established reference ranges that are based on weight and gestational age. A corticotropin stimulation test (ACTH) can be used to assess borderline cases. Genetic analysis can be helpful to confirm the diagnosis. Randomly measured 17-hydroxyprogesterone concentrations can be normal in patients with nonclassical CAH. Thus, the gold standard for diagnosis of the non-classic form is a corticotropin stimulation test, with measurement of 17-hydroxyprogesterone at 60 min.

Late onset adrenal hyperplasia:

More recently, 17- α -OHP concentrations have been utilized in the evaluation of androgenized women where late onset of 21-hydroxylase deficiency is suspected. This condition is clinically very subtle and since the presentation is the same as classical polycystic ovarian disease, basal plasma 17- α -OHP concentrations, unlike classical congenital adrenal hyperplasia, are normal. The diagnosis is made by administration of an ACTH stimulation test.

Other applications:

Measurement of 17- α -OHP concentrations is also utilized in evaluation of both men and women with acne vulgaris, male pattern baldness and in some subtle forms of infertility. Experiences with these applications are very limited.

3 PRINCIPLE OF THE TEST

The DRG 17-OH Progesterone ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the **principle of competitive binding**.

The microtiter wells are coated with a polyclonal [rabbit] antibody directed towards antigenic sites of the 17- α -OHP molecule.

Samples are pre-incubated in the coated wells.

During the second incubation, 17- α -OHP in the added sample competes with the added enzyme conjugate, which is 17- α -OHP conjugated to horseradish peroxidase, for binding to the coated antibody.

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of colour is inversely proportional to the concentration of the analyte in the sample.

A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

4 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only. For laboratory professional use only.
- Before starting the assay, read the instructions for use completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Do not reuse microtiter wells.
- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact DRG.
- Microbial contamination of reagents or samples may give false results.
- Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay.
- All indicated volumes must be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into original vials as reagent contamination may occur.

General Precautions

- Follow laboratory quality assurance and laboratory safety guidelines.
- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.

Biohazard Information

- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. However, no known test method can offer total assurance that no infectious agent is present.
- The device contains material of animal origin, which is certified apparently free of infectious or contagious diseases and injurious parasites.
- Bovine components originate from countries where BSE (Bovine spongiform encephalopathy) has not been reported.
- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation. Waste must be discarded according to local rules and regulations.

Information to Chemical Hazards and Hazard Classification

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Substrate Solution contains an ingredient in non-declarable concentrations which causes serious eye irritation. In case of possible contact with eyes, rinse immediately carefully and thoroughly with eye wash or water. After contact with skin, wash with plenty of water. Take-off contaminated clothing and wash it before reuse.
- Avoid contact with Stop Solution containing < 5 % H₂SO₄. It may cause skin irritation and burns.
- Chemicals and prepared or used reagents must be treated as hazardous waste according to the national safety guideline or regulation.
- This product does not contain substances which have carcinogenic, mutagenic or toxic for reproduction (CMR) properties.

The following kit components are classified as hazardous: Standard 0 – 6, Control low, Control high, Enzyme Conjugate, Wash Solution

Warning	Hazard statement(s): H317 - May cause an allergic skin reaction. EUH071 - Corrosive to the respiratory tract
	Precautionary statement(s): P261 - Avoid breathing dust/fume/gas/mist/vapours/ spray. P280 - Wear protective gloves. P333 + P313 - If skin irritation or rash occurs: Get medical advice/attention. P362 + P364 - Take off contaminated clothing and wash it before reuse. P501 - Dispose of contents/ container to an approved waste disposal plant

For detailed information, please refer to the Safety Data Sheet, which is available upon request directly from DRG.

5 MATERIALS

5.1 Materials Provided with the Kit

Symbol	Quantity	Description	Preparation
<i>Microtiterwells</i>	12 x 8 wells (break apart)	Microtiter plate Coated with anti-17- α -OHP antibody (polyclonal).	Ready to use
<i>Standard (Standard 0-6)</i>	7 vials x 1 mL	Standards * Concentrations: 0; 0.15; 0.5; 1.5; 3; 7.5; 20 ng/mL 0; 0.45; 1.5; 4.5; 9.1; 22.7; 60.6 nmol/L Conversion: ng/mL \times 3.03 = nmol/L. <i>Calibrated against the following reference material:</i> <i>Certified Reference Material Cerilliant H-085</i>	Ready to use
<i>Control Low & Control High</i>	2 x 1 mL	Controls * <i>For control values and ranges please refer to vial label or Certificate of Analysis.</i>	Ready to use
<i>Enzyme Conjugate</i>	1 x 25 mL	Enzyme Conjugate * 17- α -OHP conjugated to horseradish peroxidase; Colored red.	Ready to use
<i>Substrate Solution</i>	1 x 25 mL	Substrate Solution Contains 3,3',5,5'-tetramethylbenzidine (TMB). <i>Keep away from direct sun light.</i>	Ready to use
<i>Stop Solution</i>	1 x 14 mL	Stop Solution Contains < 5 % H ₂ SO ₄ . <i>Avoid contact with the stop solution. It may cause skin irritations and burns.</i>	Ready to use
<i>Wash Solution</i>	1 x 30 mL	Wash Solution, 40X concentrate *	See “Reagent Preparation”.
	1 x	Instructions for Use	
	1 x	Certificate of Analysis (CoA)	

* Contain(s) < 0.0015% CMIT/ MIT (3:1)

♦ Contain(s) 0.0108 % CMIT/ MIT (3:1)

Abbreviations:

CMIT: 5-chloro-2-methyl-4-isothiazolin-3-one

MIT: 2-methylisothiazol-3(2H)-one

5.2 Materials Required But Not Provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Manual or automatic equipment for rinsing microtiter plate wells
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

5.3 Storage and Stability of the Kit

Unopened kits and reagents as well as **opened reagents** must be stored at 2 °C to 8 °C.

The microplate must always be stored in the resealable aluminum pouch containing a desiccant. Do not open the pouch until it has reached room temperature. The microtiter plate consists of 12 individual strips. Each strip can be divided into 8 individual wells. Unused wells must be immediately returned to the aluminum pouch with the desiccant and stored again tightly resealed at 2 °C to 8 °C. Once opened, reagent vials must be closed tightly again.

	Storage Temperature	Stability
Unopened kits and unopened reagents	2 °C to 8 °C	Until the expiration date printed on the label. Do not use reagents beyond this date!
Opened kit	2 °C to 8 °C	8 weeks

5.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

Wash Solution

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

Stability after dilution:	at 20 °C to 26 °C	1 week
---------------------------	-------------------	--------

5.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

5.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components must not be used for a test run. They have to be stored until a final solution has been found. After this, they must be disposed of according to the official regulations.

6 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

Human serum or plasma (K₂ EDTA plasma, K₃ EDTA plasma, lithium heparin plasma or citrate plasma 3.2 %)

Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use hemolytic, icteric, or lipemic samples. For further information refer to chapter "*Interfering Substances*".

6.1 Sample Collection

Serum: Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma: Whole blood should be collected into centrifuge tubes containing anticoagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

Whole blood should not be frozen before centrifugation.

Stability of whole blood (10,11)	at 20 °C to 26 °C	up to 7 days
----------------------------------	-------------------	--------------

6.2 Samples Storage

Samples must be stored tightly capped prior to performing the assay. If stored frozen, freeze only once. Thawed samples must be inverted several times prior to testing.

Stability	at 2 °C to 8 °C	7 days
	at -20 °C (in aliquots)	up to 12 months

6.3 Sample Preparation

Samples can be assayed without additional preparation.

7 ASSAY PROCEDURE

7.1 Procedural Notes

- All reagents and samples must be allowed to come to room temperature (20 °C to 26 °C) before use.
- All reagents must be mixed without foaming.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- Use new disposal plastic pipette tips for each standard, control, or sample in order to avoid carry-over.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Mix the contents of the microtiter plate wells thoroughly to ensure good test results.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Once the test has been started, all steps must be completed without interruption and in the same sequence for each step.
- The enzymatic reaction is linearly proportional to time and temperature.
- Optical density is a function of the incubation time and temperature. Respect the incubations times and temperatures as given in chapter "Test Procedure".
- Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- **Important note to wash procedure:**
Washing is critical. Improperly washed wells will give erroneous results. The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
- **Test performance using fully automated analysis devices:**
Automated test performance using fully automated, open-system analysis devices is possible. However, the combination must be validated by the user.

7.2 Test Procedure

Each run must include a standard curve.

The controls serve as internal controls for the reliability of the test procedure. They must be assayed with each test run.

The given test procedure describes manual processing.

Important note: The accuracy of this assay is markedly influenced by the correct incubation temperature.

1. Secure the desired number of microtiter wells in the frame holder.
2. Pipette **25 µL** of each **Standard, Control, and sample** with new disposable tips into appropriate wells.
3. Incubate for 5 minutes at room temperature
4. Add **200 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
5. Incubate for **60 minutes** at room temperature.
6. Wash the wells as follows:
If the wash step is performed manually:
Briskly shake out the contents of the wells.
Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well.
If an automated plate washer is used:
Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well.
At the end of the washing step, always strike the wells sharply on absorbent paper to remove residual droplets!
7. Pipette **200 µL** of **Substrate Solution** to each well.
8. Incubate for **30 minutes** at room temperature.
9. Stop the enzymatic reaction by adding **100 µL** of **Stop Solution** to each well.
10. Measure the optical density (OD) of the solution in each well at **450 nm (measurement wavelength)** and at **620 nm or 630 nm (reference wavelength for recommended background subtraction)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

7.3 Calculation of Results

1. The concentration of the samples can be read directly from the standard curve.
2. For duplicate determinations, the mean of the two optical density (OD) values for each standard, control, and patient sample must be taken. If the two values deviate substantially from one another, DRG recommends retesting the samples.
3. Samples with concentrations exceeding the highest standard can be further diluted with *Standard 0* and re-assayed as described in "Test Procedure", or must be reported as > 20 ng/mL. For the calculation of the concentrations, this dilution factor must be considered.
(Example: dilution 1:10: 10 µL sample + 90 µL Standard 0)
4. Automated method:
The results in the instructions for use have been calculated automatically using a four-parameter logistic (4PL) curve fit.
(4PL Rodbard or 4PL Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. Manual method:
Using semi-logarithmic graph paper, construct a standard curve by plotting the (mean) OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.
Determine the corresponding sample concentration from the standard curve by using the (mean) OD value for each sample.

7.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Density (450 nm)
Standard 0 (0 ng/mL)	2.15
Standard 1 (0.15 ng/mL)	1.77
Standard 2 (0.5 ng/mL)	1.28
Standard 3 (1.5 ng/mL)	0.77
Standard 4 (3.0 ng/mL)	0.49
Standard 5 (7.5 ng/mL)	0.25
Standard 6 (20 ng/mL)	0.12

8 REFERENCE VALUES

The values are only for user's guidance.

It is strongly recommended for each laboratory to establish its own specific values that take into consideration a population indigenous to the area where the laboratory is located.

Values above or below the reference range should be considered as suspicious and require additional testing.

The results alone must not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

In a study conducted with newborns and children, using the DRG 17-OH Progesterone ELISA, the following values are observed:

Population	n		Mean (ng/mL)	Median (ng/mL)	2.5 th - 97.5 th Percentile (ng/mL)	Range (min - max) (ng/mL)
Newborns (boys and girls)	26	1 st month after birth	7.17	6.65	1.00 - 16.97	0.00 - 17.26
	43	2 nd month after birth	4.90	4.57	1.63 - 9.76	0.32 - 13.65
	21	3 rd month after birth	2.34	2.34	0.50 - 4.09	0.06 - 4.15
	12	4 th month after birth	2.13	2.32	0.20 - 4.32	0.20 - 4.60

Population	n	Age (years)	Mean (ng/mL)	Median (ng/mL)	2.5 th - 97.5 th Percentile (ng/mL)	Range (min - max) (ng/mL)
Children	73	1 - 12	0.95	0.83	0.08 – 1.93	0.03 - 2.13
Middle Adolescent males	11	13 - 18	1.17	1.23	0.42 – 2.24	0.41 – 2.35
Middle Adolescent females	-	13 - 18	Refer to the reference values for "Adult females"			

In a study conducted with apparently normal healthy donors, using the DRG 17-OH Progesterone ELISA, the following values are observed:

Population	n		Mean (ng/mL)	Median (ng/mL)	2.5 th – 97.5 th Percentile (ng/mL)	Range (min-max) (ng/mL)
Adult females	95	Follicular phase	1.06	1.10	0.44 – 1.60	0.30 – 1.70
	98	Luteal phase	2.28	2.40	0.90 – 4.06	0.70 – 4.60
	39	Ovulation	1.48	1.50	0.49 – 2.32	0.20 – 2.60
	26	Postmenopausal	0.38	0.32	0.17 – 0.83	0.16 – 1.04
Adult males	50	-	0.99	0.90	0.24 – 2.24	0.16 – 2.54

9 QUALITY CONTROL

Good quality assurance in the laboratory requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the Quality Control Laboratory are stated in the Certificate of Analyses (CoA) added to the kit. The values and ranges stated on the CoA always refer to the current kit lot and must be used for direct comparison of the results.

If available, it is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above-mentioned items without finding any error contact your distributor or DRG directly.

10 PERFORMANCE CHARACTERISTICS

10.1 Specificity of Antibodies (Cross-Reactivity)

The substances listed below were tested for cross-reactivity of the assay. Bias must be < 10 %.

Substance	Conc. Range (ng/mL)	Mean Bias %
17-Benzoate Estradiol	2 - 2000	-2.00
17-Cypionate Estradiol	2 - 2000	-4.75
17-Valerate Estradiol	2 - 2000	-5.32
Aldosterone	2 - 2000	-0.15
Androstenedione	2 - 2000	-6.12
Corticosterone	2 - 200	-0.92
Cortisol	2 - 200	-5.68
Cortisone	2 - 200	-1.01
DHEA	2 - 2000	-6.02
DHEA-S	2 - 2000	-1.21
Estradiol	2 - 2000	-4.39
Estriol	2 - 2000	-8.09
Estrone	2 - 2000	-5.94
Progesterone	2 - 20	1.31
Testosterone	2 - 2000	-9.39

No substantial cross-reactivity of the assay to structurally related substances is detected (<±10 %).

A bias > + 10% was found for Progesterone at concentrations ≥ 200 ng/mL.

A bias > + 10% was found for Corticosterone and Cortisol at concentrations ≥ 2000 ng/mL.

10.2 Detection Capability

Calculated according to CLSI guideline EP17-A2:2012.

Limit of Blank (LoB)	0.047 ng/mL
Limit of Detection (LoD)	0.083 ng/mL
Limit of Quantification (LoQ)	0.156 ng/mL
Lower Limit of Linear Interval (LLLI)	0.156 ng/mL
Measuring range	0.083 – 20 ng/mL
Linear range	0.156 – 20 ng/mL

10.3 Repeatability and Reproducibility

Designed on the basis of CLSI guideline EP5-A3:2014.

10.3.1 Repeatability

The repeatability was determined with 4 patient samples covering the complete measuring range in 1 run with 10 replicates. CV was calculated as mean CV of 10 replicates.

CV must be < 10 %.

Sample	n	Mean (ng/mL)	CV (%)
1 (Li-Heparin plasma)	10	0.23	4.6
2 (Serum)	10	2.10	4.0
3 (Serum)	10	7.74	3.0
4 (Citrate plasma 3.2%)	10	12.11	4.2

10.3.2 Reproducibility (Between-run Precision)

The between-run precision was determined for 4 patient samples covering the measuring range in 3 independent runs on 3 days with 10 determinations. CV was calculated from 30 determinations. CV must be < 15 %.

Sample	n	Mean (ng/mL)	CV (%)
1 (Li-Heparin plasma)	30	0.21	8.7
2 (Serum)	30	1.99	6.3
3 (Serum)	30	7.41	6.3
4 (Citrate plasma 3.2%)	30	11.53	6.2

10.3.3 Reproducibility (Between-lot Precision)

The between-lot variation was determined by 6 measurements of different samples with 3 different kit lots. CV must be < 15 %.

Sample	n	Mean (ng/mL)	CV (%)
1 (K ₃ EDTA Plasma)	18	0.30	9.4
2 (K ₃ EDTA Plasma)	18	0.86	2.5
3 (Serum)	18	1.88	5.7
4 (Serum)	18	7.22	6.7

10.4 Recovery

Recovery was determined by adding increasing amounts of the analyte to different patient samples containing different amounts of endogenous analyte. The percentage recoveries were determined by comparing expected and measured values of the samples.

	Sample 1 Citrate plasma 3.2 %	Sample 2 Li-Heparin plasma	Sample 3 K ₃ EDTA plasma	Sample 4 Serum
Highest concentration added (ng/mL)	1.86	1.86	1.86	1.86
Concentration (ng/mL)	0.17	0.46	0.72	1.92
Average Recovery (%)	97.9	96.0	95.7	93.8
Range of Recovery (%)	from	92.6	89.2	93.3
	to	103.8	103.2	98.6
				96.5

10.5 Linearity

Samples containing different amounts of analyte were serially diluted with *Standard 0*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

	Sample 1 Citrate plasma 3.2 %	Sample 2 Li-Heparin plasma	Sample 3 K ₃ EDTA plasma	Sample 4 Serum	
Highest Dilution	1:2	1:128	1:16	1:16	
Concentration (ng/mL)	1.31	39.96	10.90	15.20	
Average Recovery (%)	104.1	103.1	108.0	105.7	
Range of Recovery (%)	from to	93.9 111.4	86.9 110.2	100.0 113.8	102.4 111.1

11 LIMITATIONS OF THE PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the instructions for use and in compliance with the laboratory quality assurance guidelines.

Any improper handling of samples or modification of this test might influence the results.

11.1 Interfering Substances

11.1.1 Matrix Interference

No interference (bias < ± 20 %) was found for addition of interferent up to concentration stated in the table below.

Bilirubin unconjugated	0.30 mg/mL
Bilirubin conjugated	0.10 mg/mL
Hemoglobin	4.00 mg/mL
Triglyceride	7.50 mg/mL
Cholesterol	4.00 mg/mL
Ethanol	4.00 mg/mL
Glucose	12.50 mg/mL

11.1.2 Heterophilic Antibody Interference

The assay format of the DRG 17-OH Progesterone is not sensible to heterophilic antibodies and it is not affected by poly-specific human anti-mouse antibodies (HAMA) (12). Therefore, there is no need to test these interferences. Nevertheless, complete suppression of their effects cannot be guaranteed (9).

11.1.3 Autoantibody Interference

Autoantibodies including Rheumatoid Factors (RFs) react similar to heterophilic antibodies. The assay format of the DRG 17-OH Progesterone is not affected by heterophilic antibodies (12). Therefore, there is no need to test these interferences. Nevertheless, complete suppression of autoantibody effects cannot be guaranteed (9).

11.1.4 Drug Interferences

The following drugs were tested. Bias must be < 10 %.

Substance	Concentration range	Mean Bias
	ng/mL	%
Coumestrol	5 – 500	-1.91
Daidzein	0 – 5	8.48
Ethisterone	23.4 – 2340	-9.08
Fulvestrant	5 – 500	2.43
Genistein	5 – 500	4.10
Levonorgestrel	0.2 – 20	6.54
Mifepristone	23.4 – 2340	-9.47
Prednisolone	0.2 – 20	-1.31
Prednisone	0.2 – 2	-9.31
Secoisolariciresinol	5 – 500	8.13

A bias > + 10% was found for Daidzein at concentrations ≥ 50 ng/mL.

A bias > - 10% was found for Prednisone at concentrations ≥ 20 ng/mL.

A bias > - 10% was found for Mifepristone at concentrations ≥ 2340 ng/mL.

11.2 High-Dose Hook Effect

"High-Dose Hook Effect" is not detected up to 640 ng/mL of 17-OH Progesterone.

11.3 Trueness (Bias)

The assay is calibrated in the range of the reference material (Cerilliant; H-085). The difference from the expected value was below $\pm 20\%$.

Concentration ng/mL	Bias	Uncertainty of the reference material
0.78	-4.0%	15.1 %
5.25	-16.0 %	7.8 %
12.5	2.8 %	9.6 %


11.3.1 Method Comparison

A comparison of DRG 17-OH Progesterone (EIA-1292) (y) and the reference method (Beckman 17-OH Progesterone RIA) (x) using clinical samples gave the following correlation:

$$n = 32$$

$$r = 0.969$$

$$y = 1.029x + 0.005$$

12 LEGAL ASPECTS

12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the laboratory quality assurance guidelines and applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact DRG.

12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12.4 Reporting of Serious Incident

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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12.4 Notification des incidents graves

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1. Chrousos GP et al. Late onset 21-hydroxylase deficiency mimicking idiopathic hirsutism or polycystic ovarian disease. Annals Intern Med. 1982; 96,143-148.
2. Choi J-H, and Yoo H-W. Management issues of congenital adrenal hyperplasia during transition from pediatric to adult care. Korean J. Pediatr. 2017; 60(2), 31-37.
3. Eisenhoefer G et al. Reference intervals for plasma concentrations of adrenal steroids measured by LC-MS/MS: Impact of gender, age, oral contraceptives, body mass index and blood pressure status. Clin. Chim Acta. 2017; 470, 115-124.
4. Honour JW. 17-Hydroxyprogesterone in children, adolescents and adults. Annals of Clinical Biochemistry, 2014; 51(4), 424-440.
5. Hughes IA, Riad-Fahmy D, and Griffith K. Plasma 17OH-progesterone concentrations in newborn infants. Archives of Disease in Childhood. 1979; 54, 347-349.
6. Ishii T et al. Guidelines for diagnosis and treatment of 21-hydroxylase deficiency (2014 revision). Clin. Pediatr. Endocrinol. 2015; 24(3), 77-105.
7. Kang MJ et al. Relationships of basal level of serum 17-hydroxyprogesterone with that of serum androstenedione and their stimulated responses to a low dose of ACTH in young adult patients with congenital adrenal hyperplasia due to 21-hydroxylase deficiency. J. Korean Med. Sci. 2011; 26(11), 1454-60.
8. Maas KH et al. Relationship between 17-Hydroxyprogesterone responses to human chorionic gonadotropin and markers of ovarian follicle morphology in women with polycystic ovary syndrome. J. Clin. Endocrinol. Metab. 2015; 100(1), 293-300.
9. Ghazal et al. Hormone Immunoassay Interference: A 2021 Update. Ann Lab Med. (2022) 42(1): 3-23.
10. Guder W.G. The Quality of Diagnostic Samples; Recommendation of the working group on preanalytical quality of the German Society for Clinical Chemistry and the German Society for Laboratory Medicine, 2010, 3rd edition, GIT Verlag GmbH.
11. Diver MJ, Hughes JG, Hutton JL, West CR, Hipkin LJ. The long-term stability in whole blood of 14 commonly-requested hormone analytes. Ann Clin Biochem 1994; 31: 561-5.
12. Dimeski G. Interference Testing. Clin Biochem Rev 2008; 29 Suppl, 43-48.
13. Therrell BL. Newborn screening for congenital adrenal hyperplasia. Endocrinology&Metabolism Clinics. 2001, 30(1):15-30.
14. Berglund A et al. Epidemiology and diagnostic trends of congenital adrenal hyperplasia in Denmark: a retrospective, population-based study. The Lancet Regional Health – Europe. 2023;28: 100598.
15. Merke DP and Bornstein SR. Congenital adreal hyperplasia. Lancet 2005; 365: 2125–36.

SYMBOLS USED / VERWENDETE SYMBOLE / SIMBOLI UTILIZZATI / SÍMBOLOS UTILIZADOS / SYMBOLES UTILISES

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use	Gebrauchsanweisung beachten	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	In vitro diagnostic medical device	In-vitro-Diagnostikum	Diagnostica in vitro	Diagnóstico in vitro	Diagnostic in vitro
	Catalogue number	Katalognummer	No. di Cat.	No de catálogo	Référence
	Batch code	Chargen-bezeichnung	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests	Ausreichend für <n> Prüfungen	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit	Temperaturgrenzwerte	Temperatura di conservazione	Temperatura de conservacion	Temperature de conservation
	Use-by date	Verwendbar bis	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer	Hersteller	Fabbricante	Fabricante	Fabricant
	Distributor	Vertriebspartner	Distributore	Distribuidor	Distributeur
	Date of manufacture	Herstellungsdatum	Data di produzione	Fecha de fabricación	Date de production
	Biological risks	Biologische Risiken	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution	Achtung	Attenzione	Precaución	Attention
	Unique device Identifier	eindeutige Produktidentifizierung	Identificativo unico del dispositivo	Identificación exclusiva del dispositivo	Identifiant de dispositif unique
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Piastra per microtitolazione	Placa de microtítulo	Microplaque
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control Low</i>	Control Low	Kontrolle Niedrig	Controllo Bassa	Control Bajo	Contrôle Bas
<i>Control High</i>	Control High	Kontrolle Hoch	Controllo Alto	Control Alto	Contrôle Haut
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Coniugato enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution de substrat
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d'arresto	Solución de parada	Solution d'arrêt
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage