Instructions for Use

Strongyloides IgG/IgM ELISA

REF EIA-5812

 96

DRG 

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG 

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos / Sommaire

| | | |
|----|--|----|
| 1 | INTENDED USE..... | 3 |
| 2 | PRINCIPLE OF THE ASSAY | 3 |
| 3 | MATERIALS..... | 3 |
| 4 | STABILITY AND STORAGE | 4 |
| 5 | REAGENT PREPARATION | 4 |
| 6 | SAMPLE COLLECTION AND PREPARATION..... | 4 |
| 7 | ASSAY PROCEDURE | 5 |
| 8 | RESULTS..... | 5 |
| 9 | SPECIFIC PERFORMANCE CHARACTERISTICS | 6 |
| 10 | LIMITATIONS OF THE PROCEDURE..... | 7 |
| 11 | PRECAUTIONS AND WARNINGS | 7 |
| | | |
| 1 | VERWENDUNGSZWECK..... | 8 |
| 2 | TESTPRINZIP | 8 |
| 3 | MATERIALIEN | 8 |
| 4 | STABILITÄT UND LAGERUNG | 9 |
| 5 | VORBEREITUNG DER REAGENZIEN..... | 9 |
| 6 | ENTNAHME UND VORBEREITUNG DER PROBEN | 9 |
| 7 | TESTDURCHFÜHRUNG | 10 |
| 8 | BERECHNUNG DER ERGEBNISSE | 10 |
| 9 | TESTMERKMALE | 11 |
| 10 | GRENZEN DES VERFAHRENS..... | 12 |
| 11 | SICHERHEITSMASSNAHMEN UND WARNHINWEISE | 12 |
| | | |
| 1 | USO PREVISTO | 14 |
| 2 | PRINCIPIO DEL TEST..... | 14 |
| 3 | MATERIALI | 14 |
| 4 | MODALITÀ DI CONSERVAZIONE | 15 |
| 5 | PREPARAZIONE DEI REAGENTI..... | 15 |
| 6 | PRELIEVO E PREPARAZIONE DEI CAMPIONI | 15 |
| 7 | PROCEDIMENTO | 16 |
| 8 | RISULTATI..... | 16 |
| 9 | CARATTERISTICHE DEL TEST..... | 17 |
| 10 | LIMITAZIONI | 17 |
| 11 | PRECAUZIONI E AVVERTENZE..... | 18 |

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1 INTENDED USE

The Strongyloides IgG/IgM ELISA is intended for the qualitative determination of antibodies against Strongyloides in human serum or plasma (citrate or heparin).

2 PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microtiterplate reader.

3 MATERIALS

3.1 Reagents supplied

- **Microtiterplate:**
12 breakapart 8-well snap-off strips coated with recombinant Strongyloides antigens in resealable aluminium foil.
- **DIL:**
1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **SOLN STOP:**
1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH BUF 20x:**
1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M); pH 7.2 ± 0.2; for washing the wells; white cap; 0.2% (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Conjugate:**
1 bottle containing 20 mL of peroxidase labelled Protein A/G; coloured blue, ready to use; black cap; ≤ 0.02% (v/v) MIT.
- **SUB TMB:**
1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB) < 0.1%; ready to use; yellow cap.
- **Positive Control:**
1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:**
1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:**
1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 11.1.

3.2 Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)

3.3 Materials and Equipment needed

- ELISA microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing microtiterplate
- Pipettes to deliver volumes between 10 µL and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

4 STABILITY AND STORAGE

Store the kit at 2 °C - 8 °C.

The opened reagents are stable up to the expiry date stated on the label when stored at 2 °C - 8 °C.

5 REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20 °C - 25 °C) and mix them before starting the test run!

5.1 Microtiterplate

The break-apart snap-off strips are coated with recombinant Strongyloides antigens.

Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2 °C - 8 °C.

5.2 WASH BUF 20x

Dilute WASH BUF 20x 1 + 19; e. g. 10 mL WASH BUF 20x + 190 mL distilled water.

The diluted buffer (WASH BUF 1x) is stable for 5 days at room temperature (20 °C - 25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

5.3 SUB TMB

The reagent is ready to use and has to be stored at 2 °C - 8 °C, away from the light.

SUB TMB should be colourless or could have a slight blue tinge. If SUB TMB turns into blue, it may have become contaminated and should be thrown away.

6 SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay.

If the assay is performed within 5 days after sample collection, the samples should be kept at 2 °C - 8 °C; otherwise they should be aliquoted and stored deep-frozen (-70 °C to -20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

6.1 Sample Dilution

Before assaying, all samples should be diluted 1+100 with DIL.

Dispense 10 µL sample and 1 mL DIL into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

7 ASSAY PROCEDURE

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of **WASH|BUF|1x** from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 °C ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 °C ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of **WASH|BUF|1x**. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20 °C - 25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL **SUB|TMB** into all wells.
9. **Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL **SOLN|STOP** into all wells in the same order and at the same rate as for **SUB|TMB**, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of **SOLN|STOP**.

7.1 Measurement

Adjust the ELISA microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

8 RESULTS

8.1 Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < **0.100**
- **Negative Control:** Absorbance value < **0.200** and < **Cut-off**
- **Cut-off Control:** Absorbance value **0.150 – 1.300**
- **Positive Control:** Absorbance value > **Cut-off**

If these criteria are not met, the test is not valid and must be repeated.

8.2 Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86/2 = 0.43
Cut-off = 0.43

8.2.1 Results in Units [DU]

$$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{DRG Units} = \text{DU}]$$

Example:
$$\frac{1.591 \times 10}{0.43} = 37 \text{ DU}$$

8.3 Interpretation of Results

| | | |
|--|-----------|---|
| Cut-off | 10 DU | - |
| Positive | > 11 DU | Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine). |
| Equivocal | 9 – 11 DU | Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. |
| Negative | < 9 DU | The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely. |
| <p>Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.</p> | | |

9 SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

9.1 Precision

| Intra assay | n | Mean (E) | CV (%) |
|-------------|----|-----------|--------|
| # 1 | 24 | 0.499 | 7.27 |
| # 2 | 24 | 1.210 | 4.42 |
| # 3 | 24 | 1.106 | 3.89 |
| Inter assay | n | Mean (DU) | CV (%) |
| # 1 | 12 | 31.77 | 10.72 |
| # 2 | 12 | 24.62 | 9.84 |
| # 3 | 12 | 8.68 | 9.69 |

9.2 Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 94.12% (95% confidence interval: 83.76% - 98.77%).

9.3 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 89.47% (95% confidence interval: 75.2% - 97.06%).

9.4 Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

9.5 Cross Reactivity

Cross reactions with antibodies against other worms and parasites cannot be excluded.

10 LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

11 PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analysers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

11.1 Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

| | |
|-----------|---|
| H317 | May cause an allergic skin reaction. |
| P261 | Avoid breathing spray |
| P280 | Wear protective gloves/ protective clothing. |
| P302+P352 | IF ON SKIN: Wash with plenty of soap and water. |
| P333+P313 | If skin irritation or rash occurs: Get medical advice/ attention. |
| P362+P364 | Take off contaminated clothing and wash it before reuse. |

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Therefore, the following hazard and precautionary statements apply.

Warning

| | |
|----------------|--|
| H315 | Causes skin irritation. |
| H319 | Causes serious eye irritation |
| P280 | Wear protective gloves/ protective clothing. |
| P302+P352 | IF ON SKIN: Wash with plenty of soap and water. |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P337+P313 | If eye irritation persists: Get medical advice/attention. |

Further information can be found in the safety data sheet.

11.2 Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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14 SCHEME OF THE ASSAY


Test Preparation

| |
|--|
| <p>Prepare reagents and samples as described. Establish the distribution and identification plan for all samples and standards/controls. Select the required number of microtiter strips or wells and insert them into the holder.</p> |
|--|

Assay Procedure

| | Substrate Blank (A1) | Negative Control | Cut-off Control | Positive Control | Sample (diluted 1+100) |
|--|----------------------|------------------|-----------------|------------------|------------------------|
| Negative Control | - | 100 µL | - | - | - |
| Cut-off Control | - | - | 100 µL | - | - |
| Positive Control | - | - | - | 100 µL | - |
| Sample (diluted 1+100) | - | - | - | - | 100 µL |
| Cover wells with foil supplied in the kit Incubate for 1 h at 37 °C ± 1 °C Wash each well three times with 300 µL of WASH BUF 1x | | | | | |
| Conjugate | - | 100 µL | 100 µL | 100 µL | 100 µL |
| Incubate for 30 min at room temperature (20 °C - 25 °C) Do not expose to direct sunlight Wash each well three times with 300 µL of WASH BUF 1x | | | | | |
| SUB TMB | 100 µL | 100 µL | 100 µL | 100 µL | 100 µL |
| Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark | | | | | |
| SOLN STOP | 100 µL | 100 µL | 100 µL | 100 µL | 100 µL |
| Photometric measurement at 450 nm (reference wavelength: 620 nm) | | | | | |

SYMBOLS USED

| Symbol | English | Deutsch | Italiano | Español | Français |
|---|---|-------------------------------------|---|--|---|
|  | European Conformity | CE-Konformitäts-kennzeichnung | Conformità europea | Conformidad europea | Conformité normes européennes |
|  | Consult instructions for use * | Gebrauchsanweisung beachten * | Consultare le istruzioni per l'uso | Consulte las instrucciones de uso | Consulter les instructions d'utilisation |
| IVD | <i>In vitro</i> diagnostic medical device * | <i>In-vitro</i> -Diagnostikum * | Dispositivo medico-diagnostico in vitro | Producto sanitario para diagnóstico In vitro | Dispositif médical de diagnostic in vitro |
| REF | Catalogue number * | Artikelnummer * | Numero di Catalogo | Número de catálogo | Référence de catalogue |
| LOT | Batch code * | Chargencode * | Codice del lotto | Código de lote | Numéro de lot |
|  | Contains sufficient for <n> tests * | Ausreichend für <n> Prüfungen * | Contenuto sufficiente per "n" saggi | Contenido suficiente para <n> ensayos | Contenu suffisant pour "n" tests |
|  | Temperature limit * | Temperaturbegrenzung * | Temperatura di conservazione | Temperatura de conservación | Température de conservation |
|  | Use-by date * | Verwendbar bis * | Utilizzare prima del | Establa hasta | Utiliser jusque |
|  | Manufacturer * | Hersteller * | Fabbricante | Fabricante | Fabricant |
|  | Date of manufacture * | Herstellungsdatum * | Data di produzione | Fecha de fabricación | Date de production |
|  | Biological risks * | Biologische Risiken * | Rischi biologici | Riesgos biológicos | Risques biologiques |
|  | Caution * | Achtung * | Attenzione | Precaución | Attention |
| UDI | Unique device Identifier * | eindeutige Produktidentifizierung * | Identificativo unico del dispositivo* | Identificación exclusiva del dispositivo * | Identifiant de dispositif unique* |
| | | | | | |
| RUO | For research use only | Nur für Forschungszwecke | Solo a scopo di ricerca | Sólo para uso en investigación | Seulement dans le cadre de recherches |
|  | Distributed by | Vertreiber | Distributore | Distribuidor | Distributeur |
| <i>Content</i> | Content | Inhalt | Contenuto | Contenido | Contenu |
| <i>Volume/No.</i> | Volume / No. | Volumen/Anzahl | Volume/Quantità | Volumen/Número | Volume/Quantité |
| | | | | | |
| MTP | Microtiterplate | Mikrotiterplatte | Piastre di Microtitolazione | Placa de Microtitulación | Plaque de Microtitrage |
| CONJ | Conjugate | Konjugat | Coniugato | Conjugado | Conjugué |
| CONTROL - | Negative Control | Negativkontrolle | Controllo negativo | Control positivo | Contrôle négatif |
| CONTROL + | Positive Control | Positivkontrolle | Controllo positivo | Control positivo | Contrôle positif |
| CUT OFF | Cut off control | Cut off-Kontrolle | Controllo Cut-off | Control Cut-off | Contrôle Cut-off |
| CAL | Standard or Calibrator | Standard oder Kalibrator | Standard o Calibratore | Estándar o Calibrador | Standard o Etalon |
| DIL | Sample Dilution Buffer | Probenverdünnungspuffer | Tampone di Diluizione del Campione | Tampón de Dilución de Muestra | Tampon de Dilution d'Échantillon |
| DIL A | IgA Sample Dilution Buffer | IgA-Probenverdünnungspuffer | Tampone di Diluizione del Campione IgA | Tampón de Dilución de Muestra IgA | Tampon de Dilution d'Échantillon IgA |
| DIL G | IgG Sample Dilution Buffer | IgG-Probenverdünnungspuffer | Tampone di Diluizione del Campione IgG | Tampón de Dilución de Muestra IgG | Tampon de Dilution d'Échantillon IgG |
| DIL M | IgM Sample Dilution Buffer | IgM-Probenverdünnungspuffer | Tampone di Diluizione del Campione IgM | Tampón de Dilución de Muestra IgM | Tampon de Dilution d'Échantillon IgM |
| SOLN STOP | Stop solution | Stopplösung | Soluzione bloccante | Solución de parada | Solution d'arrêt |