Instructions for Use

Anti Zona Pellucida Ab ELISA

IVD

CE

REF EIA-6138

Σ 96

DRG

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la versión válida de la metodico técnico incluido aquí en el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.

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1	INTENDED USE	2
2	PRINCIPLE OF THE TEST	2
3	WARNINGS AND PRECAUTIONS	2
4	REAGENTS	3
5	SPECIMEN COLLECTION AND PREPARATION	4
6	ASSAY PROCEDURE	5
7	EXPECTED NORMAL VALUES	6
8	QUALITY CONTROL	6
9	PERFORMANCE CHARACTERISTICS	6
10	LIMITATIONS OF USE	7
11	LEGAL ASPECTS	7

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1 INTENDED USE

The DRG Anti Zona Pellucida Ab ELISA is an enzyme immunoassay for the quantitative *in vitro* diagnostic measurement of antibodies directed against human zona pellucida in serum.

1.1 Information

Antibodies directed against zona pellucida antigens may cause infertility in women. The detection of the anti-zona pellucida antibodies may serve as an aid for the diagnosis of immunologically caused disorders of fertility.

2 PRINCIPLE OF THE TEST

The DRG Anti Zona Pellucida Ab ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle.

The microtiter wells are coated with zona pellucida antigen.

During incubation, specific human anti-zona pellucida antibodies in the patient sample bind to the coated surface of the wells.

A washing step removes unbound sample components.

Added enzyme conjugate binds to the immobilized antibody-antigen-complexes.

The conjugate contains polyclonal antibodies directed against epitopes on the Fc part of human immunoglobulins and labelled with horseradish peroxidase (HRP).

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution.

The colorimetric reaction is abruptly stopped by addition of stop solution and optical density (OD) of the resulting yellow product is measured.

The intensity of colour is proportional to the concentration of specific antibodies in the sample.

A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.

19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

4 REAGENTS

4.1 Materials Provided

Symbol	Quantity	Description	Preparation
Microtiterwells	12 x 8 wells (break apart)	Microtiter plate Coated with a mix of zona pellucida antigens	Ready to use
Standard (Standard [0 - n])	4 vials x 0.5 mL	Standards * Concentrations: 6 – 25 – 50 – 100 U/mL	Ready to use
Quality Control	1 vial x 0.5 mL	Controls * <i>For control values and ranges please refer to vial label or Certificate of Analysis.</i> <i>Contains material of animal origin;</i>	Ready to use
Dilution Buffer / Zero Standard	1 x 50 mL	Sample Diluent *	Ready to use
Enzyme Conjugate	1 x 6 mL	Enzyme Conjugate * Anti-human poly Ig conjugated to horseradish peroxidase Colored red. <i>Contains material of animal origin;</i>	Ready to use
Substrate Solution	1 x 14 mL	Substrate Solution Contains 3,3',5,5'-tetramethylbenzidine (TMB). <i>Keep away from direct sun light.</i>	Ready to use
Stop Solution	1 x 14 mL	Stop Solution Contains < 5 % H ₂ SO ₄ . <i>Avoid contact with the stop solution. It may cause skin irritations and burns.</i>	Ready to use
Wash Solution	1 x 30 mL	Wash Solution, 40X concentrate♦	See "Reagent Preparation".
	1 x	Cover foil	
	1 x	Instructions for Use	
	1 x	Certificate of Analysis (CoA)	

* Contain(s) < 0.0015 % CMIT/ MIT (3:1)

♦ Contain(s) 0.0108 % CMIT/ MIT (3:1)

Abbreviations:

CMIT: 5-chloro-2-methyl-4-isothiazolin-3-one

MIT: 2-methylisothiazol-3(2H)-one

BND: 5-Bromo-5-nitro-1,3-dioxane

4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm) (e.g. the DRG Instruments Microtiter Plate Reader)
- Calibrated variable precision micropipettes
- Tubes for sample dilution
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for 4 weeks if stored as described above.

4.4 Reagent Preparation

Wash Solution

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 1 week at room temperature.

Stability after dilution:	at 20 °C to 26 °C	1 week
---------------------------	-------------------	--------

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum can be used in this assay.

Note: Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use haemolytic, icteric, or lipaemic specimens.

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 7 days at 2 °C to 8 °C prior to assaying.

Specimens held for a longer time (up to 12 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

Prior to assaying, dilute each patient specimen **1 / 100** with *Dilution Buffer*.

Example:

Dilution 1:100: 5 µL sample + 495 µL *Dilution Buffer* (mix thoroughly)

Note: The *Quality Control* is ready to use and must not be diluted!

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Optical density is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve. All reagents and the required number of wells should be brought to room temperature (20 °C to 26 °C) before use.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **50 µL** each of **Zero Standard**, **Standard**, **Quality Control** and **diluted samples** with new disposable tips into appropriate wells.
3. Cover with foil and incubate for **60 minutes** at **37 °C**.
4. Rinse the wells **3 times** with **400 µL** diluted *Wash Solution* per well, if a plate washer is used.
- OR -
Briskly shake out the contents of the wells.

Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.

Important note:

The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!

5. Dispense **50 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
6. Cover with foil and incubate for **60 minutes** at **37 °C**.
7. Rinse the wells **5 times** with **400 µL** diluted *Wash Solution* per well, if a plate washer is used.
- OR -
Briskly shake out the contents of the wells.
Rinse the wells **5 times** with **300 µL** diluted *Wash Solution* per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.
8. Add **50 µL** of **Substrate Solution** to each well.
9. Incubate for **30 minutes** at room temperature.
10. Stop the enzymatic reaction by adding **100 µL** of **Stop Solution** to each well.
11. Measure the optical density (OD) of the solution in each well at **450 nm (measurement wavelength)** and at **620 nm or 630 nm (reference wavelength for recommended background subtraction)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve.
6. The standards are already pre-diluted, therefore the 1/100 dilution of the samples must not be taken into account for the final calculation of sample concentrations.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Units (450 nm)
Zero Standard	0.03
Standard 1	0.42
Standard 2	1.08
Standard 3	1.69
Standard 4	2.40

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

Normal values 0 - 10 U/mL

Elevated values > 10 U/mL

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each standard curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Detection Capability

Calculated according to CLSI guideline EP17-A2:2012.

Limit of Blank (LoB)	0.097 U/mL
Limit of Detection (LoD)	4.170 U/mL
Limit of Quantification (LoQ)	15.878 U/mL
Measuring range	4.170 – 100 U/mL

9.2 Repeatability

For the determination of the repeatability (within-run), 6 kits from 6 different batches (produced on different days) were used. One patient sample (OD > 1.0) was applied 96 times per testing procedure.

Mean CV: 6.8% (range from 5.1% – 7.9 %)

9.3 Reproducibility (Between-run Precision)

For the determination of the reproducibility (between-run), one strip each of 12 kits stemming from 6 different batches (produced on different days) were used. One patient sample (OD > 1.0) was applied 72 times per testing procedure.

Mean CV: 7.2% (range from 5.5% – 9.0 %)

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Severely haemolytic or lipaemic sera or sera from patients with liver diseases should not be used.

Results may be adversely affected by certain pathologic conditions, such as poly- and monoclonal gammopathies, autoimmune diseases or by an altered immune status.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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12 REFERENCES / LITERATURE

1. Hasegawa H, Tanaka H and Shibahara H. Infertility and immunocontraception based on zona pellucida.
2. Reprod Med Biol, 2014, 13 :1-9.
3. Sacco A, and Moghissi K. Anti-zona pellucida activity in human sera.
4. Fertility and Sterility, 1979, 31(5), 503-6.
5. Bronson RA, Cooper GW, and Rosenfeld DL. Sperm-specific isoantibodies and autoantibodies inhibit the binding of human sperm to the human zona pellucida.
6. Fertility and Sterility, 1982, 38 (6), 724-9.
7. Mori T et al. A method for specific detection of autoantibodies to the zona pellucida in infertile women.
8. Fertility and Sterility, 1979, 32(1), 67-72.
9. Nishimoto T et al. Autoantibodies to zona pellucida in infertile and aged women.
10. Fertility and Sterility, 1980, 34(6), 552-5.
11. Mohammadi Yeganeh et al. The association of different auto-antibodies against ovarian tissues and gonadotropins and poor ovarian response in intracytoplasmic sperm injection cycles.
12. Human Fertility (Cambridge, England), 2017, 20(2), 126-131.
13. Kamada M et al. Sperm–zona pellucida interaction and immunological infertility.
14. Reproductive Medicine and Biology, 2006; 5, 95 –104.
15. Ivanova M. et al. Zona pellucida autoantibodies in women undergoing ART.
16. Folia Biologica, 1999, 45(2), 59-62.
17. Koyama K et al. Follicular dysfunction induced by autoimmunity to zona pellucida.
18. Reproductive Biology, 2005, 5(3), 269-278.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnóstico in vitro	Producto sanitario para diagnóstico in vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	No. di Cat.	No de catálogo	Référence
	Batch code *	Chargencode *	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservacion	Température de conservation
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Distributor *	Vertriebspartner *	Distributore	Distribuidor	Distributeur
	Date of manufacture *	Herstellungsdatum *	Data di produzione	Fecha de fabricación	Date de production
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution *	Achtung *			
	Unique device Identifier *	eindeutige Produkt-identifizierung *	Identificativo unico del dispositivo*	Identificación exclusiva del dispositivo *	Identifiant de dispositif unique*
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement
<i>Volume/No.</i>	Volume / No.	Volumen / Anzahl	Volume / Quantità	Volumen / Número	Volume / Quantité
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Micropozzetti	Placas multipicillo	Microplaques
<i>Antiserum</i>	Antiserum	Antiserum	Antisiero	Antisuero	Antisérum
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complejo enzimático	Complexe enzymatique
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d'arresto	Solución de parada	Solution d'arrêt

<i>Zero Standard</i>	Zero Standard	Nullstandard	Standard zero	Estándar cero	Zero Standard
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)	1N NaOH	1N NaOH
<i>1 N HCl</i>	1 N HCl	1 N HCl		1 N HCl	1N HCl
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungs-medium	Diluente dei campioni	Solución para dilución de la muestra	Solution pour dilution de l'échantillon
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungs-medium	Diluente del tracciante	Solución para dilución del conjugado	Solution pour dilution du conjugué