## Instructions for Use

# Free Testosterone ELISA

IVD

REF EIA-2924

$\Sigma$  96

**DRG** 

DRG Instruments GmbH, Germany  
Frauenbergstraße 18, 35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

**DRG** 

DRG International, Inc., USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.  
 Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
 Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.  
 Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.  
 Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.  
 Por favor, usar a versão válida das instruções de utilização fornecidas com o kit.**

<b>Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas</b>	
<p>The following changes have been made in comparison to the previous version:            Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen:            Rispetto alla versione precedente, sono state apportate le seguenti modifiche:            Se han introducido los siguientes cambios en comparación con la versión anterior:</p>	
<p>Detailed editorial revision. Changed wording in several chapters.            Change to a new microtiter plate.</p>	
<b>Chapter 3</b>	exchange of wash buffer and Stop Solution information. Returning to previous volumes and concentrations
<b>Chapter 8.2</b>	Adjusting the wash buffer preparation to the above.
<b>Chapter 14.2</b>	correction of reference to chapter 14.5

---

**Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos**

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## 1 INTENDED PURPOSE

Free Testosterone ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of free testosterone in human serum or plasma.

Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis and monitoring of disorders involving the male sex hormones (androgens).

### For *In Vitro* Diagnostic Use

### For Laboratory Professional Use

#### 1.1 CLINICAL SIGNIFICANCE

Testosterone is found in circulation predominantly linked to carrier proteins, the most common of which being sex-hormone binding globulin (SHBG). Testosterone plays a key role in the development of primary and secondary sexual characteristics in males and is involved in the production of female sexual hormones.

Only 1 – 2% of testosterone in circulation is not bound to any protein and is biologically active – this is referred to as ‘free testosterone’ (FT). Bioavailable testosterone refers to the sum of FT and the testosterone bound to serum albumin, since it is bound with low affinity and readily able to dissociate to become available for its biological function.

In males elevated levels of testosterone are associated with several conditions such as, early (precocious) puberty, congenital adrenal hyperplasia (CAH), androgen insensitivity syndrome (AIS), steroid use and testicular or adrenal tumours. Whereas the major causes of suppressed levels include Klinefelter’s syndrome, testicular damage, pituitary disorders etc. In females of all ages, elevated testosterone levels can be associated with a variety of virilising conditions including adrenal tumours and polycystic ovarian syndrome (PCOS).

These clinical conditions are associated with either a lack or excess of testosterone in circulation (hypoandrogenism or hyperandrogenism). Diagnosis of these disorders involve the quantification of total testosterone (TT) in association with other clinical evidence and laboratory data. However, clinical manifestations of androgen disorders are often associated with normal levels of TT. In such cases, additional information may be gained by the assessment of the biologically active, FT level. Several androgen disorders can be caused by alteration of SHBG production which affects the levels of FT available in serum.

Measurement of FT can be considered useful in the diagnosis of several conditions including androgen deficiency in men and androgen excess in women<sup>1</sup>. Assessment of free testosterone levels may prove beneficial<sup>2</sup> and may avoid an incorrect diagnosis of hypogonadism in cases when low concentrations of total testosterone are determined and alterations of SHBG are suspected.

There is an observed and well documented circadian variation of testosterone levels in men with circulating concentrations being higher in the morning and declining throughout the day<sup>3</sup>. Testosterone levels also decline in ageing males (andropause) and is associated with loss of muscle and bone mass, leading to osteoporosis, loss of libido, erectile dysfunction. Depression and impaired cognitive function<sup>4</sup>.

## 2 PRINCIPLE OF THE METHOD

The Free Testosterone ELISA is a competitive enzyme immunoassay (ELISA) where free testosterone (antigen) in the sample competes with the antigenic testosterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti-testosterone coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate ( $H_2O_2$ ) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution ( $H_2SO_4$ ) is added. The colour intensity is inversely proportional to the free testosterone concentration of in the sample.

Free testosterone concentration in the sample is calculated through a calibration curve.

### 3 REAGENTS, MATERIALS, AND INSTRUMENTATION

#### 3.1 Reagents and materials provided with the kit

1. **Standards (Standard 0 – 5)** (6 vials, 1 mL each)
2. **Controls** (Control A and Control B) (2 vials, 1 mL)  
Control concentration is indicated on the Certificate of Analysis
3. **Enzyme Conjugate** (1 vial, 15 mL)  
Testosterone conjugated with Horseradish peroxidase (HRP)
4. **Microtiterwells** (1 breakable microplate)  
Anti-Testosterone antibody adsorbed on microplate
5. **Substrate Solution** (1 vial, 15 mL)  
H<sub>2</sub>O<sub>2</sub>-TMB 0.26 g/L (avoid any skin contact)
6. **Stop Solution** (1 vial, 15 mL)  
Sulphuric acid 0.15 M
7. **Wash Solution**, 10X Concentrated (2 vial, 25 mL);  
0.2M Phosphate Buffer, pH 7.4  
see "Reagent Preparation"

#### 3.2 Materials required but not provided

Distilled water

#### 3.3 Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

### 4 WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in humans or animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV, or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of ProClin™ 300 as preservative. Avoid contact with skin or mucosa.
- The TMB Substrate contains an irritant, which is harmful if inhaled, ingested, or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals, or oxidants. Do not freeze the solution.
- The clinical significance of the determination of free testosterone can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

### 5 PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 °C - 8 °C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 °C - 28 °C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6 REAGENT STORAGE AND STABILITY

Store the kit at 2 °C - 8 °C in the dark.

- The kit is stable at 2 °C - 8 °C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 °C - 8 °C for 6 months.
- The diluted wash solution is stable for 30 days at 2 °C - 8 °C.

**Important note:** open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

## 7 SAMPLE COLLECTION AND STORAGE

The assay should be performed using

**serum** (standard sampling tubes or tubes containing serum separating gel) or **plasma** (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 °C - 8 °C	24 hours
Freeze/thaw cycles	1 cycle

## 8 PROCEDURE

### 8.1 Preparation of Standards and Controls

Before use, mix for 5 minutes with a rotating mixer.

The Standards are ready for use and have the following concentration of Testosterone:

	S0	S1	S2	S3	S4	S5
pg/mL	0	0.2	1.0	4.0	20.0	100.0

### 8.2 Preparation of the Wash Solution

Dilute the content of the two vials "Wash solution" with distilled water to a final volume of 500 mL prior to use (for both vials 25 mL + 25 mL + 450 mL)

For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

### 8.3 Preparation of Samples

The determination of free testosterone can be performed in human serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

The Controls are ready to use.

## 8.4 Procedure

- **Allow all reagents to reach room temperature (22 °C - 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 °C - 8 °C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 °C - 8 °C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (S0 - S5), two for each Control, two for each sample, one for Blank.

Reagent	Standard	Sample / Control	Blank
Standard S0-S5	20 µL		
Sample / Control		20 µL	
Enzyme Conjugate	100 µL	100 µL	
Incubate at 37 °C for 1 hour. Remove the content from each well; wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate Solution	100 µL	100 µL	100 µL
Incubate at room temperature 22 °C – 28 °C for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 9 QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DRG recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories<sup>5</sup>.

## 10 CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls.

A 4-parameter logistic (4PL) curve fit, **including Standard 0 is required.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

### Conversion of units

To convert results to SI units:

$$\text{pmol/L} = \text{pg/mL} \times 3.47$$

To convert results to mass units:

$$\text{pg/mL} = \text{pmol/L} \times 0.289$$

## 11 MEASURING RANGE

The assay measuring range (AMR) is 0.34 – 70.0 pg/mL (1.18 – 242.9 pmol/L).

Any value that reads below 0.34 pg/mL (1.18 pmol/L) should be reported as “< 0.34 pg/mL (< 1.18 pmol/L)”.

Any value that reads above 70.0 pg/mL (242.9 pmol/L) should be reported as “> 70.0 pg/mL (242.9 pmol/L)”.

## 12 METROLOGY AND TRACEABILITY

The Free Testosterone ELISA has been standardised against internal reference standards (serum matrix) which have been value assigned to another commercially available test method.

## 13 EXPECTED VALUES

The following ranges were determined using the Free Testosterone ELISA and are provided for information only.

	No. of subjects	Median pg/mL	Reference Interval (pg/mL)
<b>Males</b>			
<b>21 – 49 years</b>	120	14.13	5.01 – 27.78
<b>&gt; 50 years</b>	120	12.75	4.11 – 21.85
<b>Females</b>			
<b>pre-menopausal</b>	120	0.55	< LOQ – 1.70
<b>post-menopausal</b>	120	0.75	< LOQ – 2.34

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

## 14 PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

### 14.1 Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.10 pg/mL
Limit of Detection (LoD)	0.20 pg/mL
Limit of Quantitation (LoQ)	0.34 pg/mL

## 14.2 Trueness

Trueness has been demonstrated through method comparison of the Free Testosterone ELISA to a commercially available assay using native donor samples – refer to section 14.5

## 14.3 Precision

Precision of the Free Testosterone ELISA was determined by performing a complex precision study.

**Repeatability:** A total of 6 serum samples were assayed in 5 replicates, once per day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (pg/mL)	Within run (Repeatability)	
			SD	CV%
1	75	0.68	0.08	11.3%
2	75	1.62	0.10	5.9%
3	75	5.25	0.19	3.7%
4	75	10.43	0.48	4.6%
5	75	34.99	1.17	3.3%
6	75	68.65	4.62	6.7%

**Reproducibility (between run):** A total of six serum samples were assayed in 5 replicates, once per day for 5 days by 3 operators and with 2 lots.

Sample	n	Mean (pg/mL)	CV%
PS 1	150	0.75	11.1%
PS 2	150	1.74	6.6%
PS 3	150	5.48	3.8%
PS 4	150	10.72	4.1%
PS 5	150	37.41	5.5%
PS 6	150	73.81	6.9%

## 14.4 Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures".

For free testosterone concentration by Free Testosterone ELISA, the measurement procedure shows linearity for the interval from 0.27 to 89.83 pg/mL within the allowable deviation of linearity (ADL) of  $\pm 15\%$ .

## 14.5 Method comparison

The Free Testosterone ELISA was compared against a commercially available quantitative manual ELISA test, following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 95 samples, selected to represent a wide range of free testosterone concentrations, was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
95	0.94 [0.88 to 0.99]	-0.02 [-0.14 to 0.06]	0.93

#### 14.6 Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested	Mean %Cross reactivity
11-keto-testosterone	100 ng/mL	0.1%
11- $\beta$ -hydroxy-testosterone	10 ng/mL	0.2%
17 $\alpha$ OH-Progesterone	500 ng/mL	0.0%
Aldosterone	3000 ng/mL	0.0%
Androstenedione	100 ng/mL	0.0%
Cortisol	1000 ng/mL	0.0%
Cortisone	1000 ng/mL	0.0%
Danazol	1000 ng/mL	0.0%
Dexamethasone	2000 ng/mL	0.0%
DHEA	1000 ng/mL	0.0%
DHEA-S	10000 ng/mL	0.0%
5 $\alpha$ -dihydrotestosterone (DHT)	500 ng/mL	0.0%
Estradiol	1000 ng/mL	0.0%
Estriol	100 ng/mL	0.0%
Estrone	1000 ng/mL	0.0%
Ethisterone	100 ng/mL	0.0%
Norgestrel	100 ng/mL	0.0%
Prednisone	1000 ng/mL	0.0%
Pregnenolone	5000 ng/mL	0.0%
Progesterone	1000 ng/mL	0.0%
Testosterone propionate	1000 ng/mL	0.0%

The following substances do not interfere with a bias of  $> \pm 15\%$  in the Free Testosterone ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	7 g/dL
Triglyceride	500 mg/dL

#### 14.7 Serum-plasma study

The Free Testosterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 21 samples (18 native, 3 spiked) to cover the range were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
SST	0.97 [0.92 to 1.02]	-0.02 [-1.01 to 0.98]	0.99
Lithium Heparin	0.99 [0.97 to 1.02]	-0.01 [-0.40 to 1.02]	1.00
Sodium Heparin	0.93 [0.87 to 0.99]	-0.01 [-1.12 to 1.10]	0.99
EDTA	1.03 [0.97 to 1.08]	0.27 [0.74 to 1.28]	0.99

## 15 LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays<sup>6</sup>. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

## 16 WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

## 17 REPORTING OF SERIOUS INCIDENT

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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**18 BIBLIOGRAPHY**

1. Shea JL, Wong PY, Chen Y. Free testosterone: clinical utility and important analytical aspects of measurement. *Adv Clin Chem*. 2014;63:59-84.
2. Diver MJ. Analytical and physiological factors affecting the interpretation of serum testosterone concentration in men. *Ann Clin Biochem*. 2006 Jan;43(Pt 1):3-12.
3. Brambilla DJ, Matsumoto AM, Araujo AB and McKinlay JB. The Effect of Diurnal Variation on Clinical Measurement of Serum Testosterone and Other Sex Hormone Levels in Men. *J Clin Endocrinol Metab*. 2009 Mar; 94(3): 907–913.
4. Rajfer J. Decreased Testosterone in the Aging Male. *Rev Urol*. 2003;5(suppl 1):S1–S2.
5. Basic QC Practices On-line Course; <http://www.Westgard.com>.
6. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

## 19 SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français	Português
	European Conformity	CE-Konformitätskennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes	Conformidade Europeia
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation	Consultar as instruções de uso
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Diagnostica in vitro	Diagnóstico in vitro	Diagnostic in vitro	Dispositivo médico para diagnóstico in vitro
	Catalogue number *	Katalognummer *	No. di Cat.	No de catálogo	Référence	Número de catálogo
	Batch code *	Chargen-bezeichnung *	Lotto no	Número de lote	No. de lot	Código do lote
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests	Suficiente para <n> determinações
	Temperature limit *	Temperaturgrenzwerte *	Temperatura di conservazione	Temperatura de conservación	Temperature de conservation	Limites de temperatura
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation	Prazo de validade
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant	Fabricante
	Distributor *	Vertriebspartner *	Distributore	Distribuidor	Distributeur	Distribuidor
	Date of manufacture *	Herstellungsdatum *	Data di produzione	Fecha de fabricación	Date de production	Data de fabricação
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques	Riscos biológicos
	Caution *	Achtung *	Attenzione	Precaución	Attention	Cuidado
	Unique device Identifier *	eindeutige Produktidentifizierung *	Identificativo unico del dispositivo*	Identificación exclusiva del dispositivo *	Identifiant de dispositif unique*	Identificador único do dispositivo *
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches	
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement	Conteúdo
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité	Volume / Quantidade