Instructions for Use

RSV IgM ELISA

IVD

CE

REF EIA-3508

Σ 96

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG[®] S[®]
DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the package insert provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Arbeitsanleitung.
Si prega di usare la versione valida dell'inserto del pacco a disposizione con il kit.**

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1	INTRODUCTION	2
2	PRINCIPLE OF THE TEST.....	2
3	WARNINGS AND PRECAUTIONS.....	3
4	KIT COMPONENTS.....	4
5	SPECIMEN COLLECTION AND PREPARATION	5
6	ASSAY PROCEDURE	5
7	RESULTS	8
8	QUALITY CONTROL	8
9	ASSAY CHARACTERISTICS	9
10	LIMITATIONS OF USE	10
11	LEGAL ASPECTS.....	10
12	REFERENCES / LITERATUR.....	11

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1 INTRODUCTION

1.1 Intended Use

The **DRG RSV IgM Enzyme Immunoassay Kit** provides materials for the **qualitative** and **semiquantitative** determination of IgM-class antibodies to Respiratory Syncytial Virus in serum.

This assay is intended for **in vitro diagnostic use only**.

1.2 Summary and Explanation

Respiratory syncytial virus (RSV) is a negative-sense, enveloped RNA virus. The virion is variable in shape and size (average diameter of between 120 and 300 nm), is unstable in the environment (surviving only a few hours on environmental surfaces), and is readily inactivated. It is the causative pathogen for the most common infection of the respiratory tract. Respiratory syncytial virus (RSV) infections usually occur during annual community outbreaks during the late fall, winter, or early spring months. The timing and severity of outbreaks in a community vary from year to year. Most infants become infected with RSV in their first winter season; between 25% and 40% have signs or symptoms of bronchiolitis or pneumonia, and 0.5% to 2% require hospitalization. Most children recover from illness in 8 to 15 days; the majority of children hospitalized for RSV infection is under 6 months of age. Most children will have serological evidence of RSV infection by 2 years of age. RSV also causes repeated infections throughout life, usually associated with moderate-to-severe cold-like symptoms; however, severe lower respiratory tract disease may occur at any age, especially among elderly or among those with compromised cardiac, pulmonary, or immune systems.

2 PRINCIPLE OF THE TEST

The **DRG RSV IgM ELISA Kit** is a solid phase enzyme-linked immunosorbent assay (ELISA)

Patient samples are diluted with *Sample Diluent* and additionally incubated with *IgG-RF-Sorbent*, containing hyper-immune anti-human IgG-class antibody to eliminate competitive inhibition from specific IgG and to remove rheumatoid factors. This pretreatment avoids false negative or false positive results.

Microtiter wells as a solid phase are coated with Respiratory Syncytial Virus antigen.

Pretreated patient specimens and **ready-for-use controls** are pipetted into these wells. During incubation Respiratory Syncytial Virus-specific antibodies of positive specimens and controls are bound to the immobilized antigens.

After a washing step to remove unbound sample and control material horseradish peroxidase conjugated anti-human IgM antibodies are dispensed into the wells. During a second incubation this anti-IgM conjugate binds specifically to IgM antibodies resulting in the formation of enzyme-linked immune complexes.

After a second washing step to remove unbound conjugate the immune complexes formed (in case of positive results) are detected by incubation with TMB substrate and development of a blue color. The blue color turns into yellow by stopping the enzymatic indicator reaction with sulfuric acid.

The intensity of this color is directly proportional to the amount of Respiratory Syncytial Virus-specific IgM antibody in the patient specimen. Absorbance at 450 nm is read using an ELISA microtiter plate reader.

3 WARNINGS AND PRECAUTIONS

- This kit is for in vitro diagnostic use only. For professional use only.
- Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- Avoid contact with Stop Solution containing 0.2 mol/L H₂SO₄. It may cause skin irritation and burns.
- TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
- The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided
- Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
- Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
- Never pipette by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
- Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
- Handling should be in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
- Do not use reagents beyond expiry date as shown on the kit labels.
- All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Chemicals and prepared or used reagents have to be treated as hazardous waste according the national biohazard safety guideline or regulation.
- For information on hazardous substances included in the kit please refer to Safety Data Sheets.
Safety Data Sheets for this product are available upon request directly from DRG.

4 KIT COMPONENTS

4.1 Contents of the Kit

Symbol	Quantity	Description	Preparation
Microtiterwells	12 x 8 wells (break apart)	Microtiter plate Wells coated with Respiratory Syncytial Virus (RSV) Antigen.	Ready to use
Sample Diluent	1 x 100 mL	Sample Diluent * colored yellow, pH 7.2 ± 0.2	Ready to use
IgG-RF-Sorbent	1 x 6,5 mL	IgG-RF-Sorbent * colored yellow; Contains anti-human IgG-class antibody.	Ready to use
Neg. Control	1 x 2 mL	Negative Control * Colored yellow, yellow cap.	Ready to use
Cut-off Control	1 x 2 mL	Cut-off Control * Colored yellow, black cap.	Ready to use
Pos. Control	1 x 2 mL	Positive Control * Colored yellow, red cap.	Ready to use
Enzyme Conjugate	1 x 20 mL	Enzyme Conjugate * Colored red, antibody to human IgM conjugated to horseradish peroxidase.	Ready to use
Substrate Solution	1 x 14 mL	Substrate Solution Contains 3,3',5,5'-tetramethylbenzidine (TMB). <i>Keep away from direct sun light.</i>	Ready to use
Stop Solution	1 x 14 mL	Stop Solution Contains 2 % H ₂ SO ₄ . <i>Avoid contact with the stop solution. It may cause skin irritations and burns.</i>	Ready to use
Wash Solution	1 x 30 mL	Wash Solution, 20X concentrate *	See "Reagent Preparation".
	1 x	Cover foil	
	1 x	Instructions for Use	
	1 x	Certificate of Analysis (CoA)	
		* Contain(s) non-mercury preservative.	

4.1.1 Equipment and material required but not provided

- A microtiter plate calibrated reader (450/620 nm ±10 nm) (e.g. the DRG Instruments Microtiter Plate Reader)
- Calibrated variable precision micropipettes
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Vortex tube mixer
- Deionised or (freshly) distilled water
- Timer
- Absorbent paper

4.2 Storage and Stability of the Kit

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for two months if stored as described above.

4.3 Reagent Preparation

Allow all reagents and required number of strips to reach room temperature prior to use.

Wash Solution

Dilute **Wash Solution 1+19** (e.g. 10 mL + 190 mL) with fresh and germ free redistilled water. This diluted wash solution has a pH value of 7.2 ± 0.2 .

Consumption: ~ 5 mL per determination.

Crystals in the solution disappear by warming up to 37°C in a water bath. Be sure that the crystals are completely dissolved before use.

The diluted Wash Solution is stable for 1 week at 2°C to 8°C .

4.4 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Material Safety Data Sheets.

4.5 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum can be used in this assay.

Do not use haemolytic, icteric or lipaemic specimens.

5.1 Specimen Collection**Serum:**

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

5.2 Specimen Storage

Specimens should be capped and may be stored for up to 5 days at 2°C to 8°C prior to assaying.

Specimens held for a longer time (up to 12 months) should be frozen only once at -20°C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

Prior to assaying each patient specimen is first to be diluted with *Sample Diluent*. For the absorption of rheumatoid factor these prediluted samples then have to be incubated with *IgG-RF-Sorbent*

1. Dilute each patient specimen **1+50** with *Sample Diluent*,
e.g. 10 µL of specimen + 0.5 mL of *Sample Diluent*. **Mix well**.
2. Mix well *IgG-RF-Sorbent* before use.
3. Dilute this prediluted sample **1+1** with *IgG-RF-Sorbent*
e.g. 60 µL prediluted sample + 60 µL *IgG-RF-Sorbent*. **Mix well**
4. **Let stand at room temperature for at least 15 minutes, up to a maximum of 2 hours and mix well again.**
5. Take 100 µL of these pretreated samples for the ELISA.

Please note: Controls are ready for use and must not be diluted!

6 ASSAY PROCEDURE**6.1 General Remarks**

- **It is very important to bring all reagents, samples and controls to room temperature before starting the test run!**
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.

- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- During incubation cover microtiter strips with foil to avoid evaporation.

6.2 Test Procedure

Prior to commencing the assay, dilute *Wash Solution*, **prepare patient samples as described in point 5.3** and establish carefully the **distribution and identification plan** supplied in the kit for all specimens and controls.

1. Select the required number of microtiter strips or wells and insert them into the holder.

Please allocate at least:

1 well (e.g. A1) for the *Neg. Control*,
 2 wells (e.g. B1+C1) for the *Cut-off Control* and
 1 well (e.g. D1) for the *Pos. Control*.

It is left to the user to determine controls and patient samples in duplicate.

2. Dispense

100 µL of *Neg. Control* into well A1
100 µL of *Cut-off Control* into wells B1 and C1
100 µL of *Pos. Control* into well D1 and
100 µL of each pre treated sample with new disposable tips into appropriate wells.

3. Cover wells with foil supplied in the kit. Incubate for **60 minutes at 37 °C**.

4. Briskly shake out the contents of the wells.

Rinse the wells **5 times** with diluted *Wash Solution* (**300 µL per well**). Strike the wells sharply on absorbent paper to remove residual droplets.

Important note:

The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!

5. Dispense **100 µL Enzyme Conjugate** into each well.

6. Incubate for **30 minutes at room temperature (20 °C to 25 °C)**.

Do not expose to direct sun light!

7. Briskly shake out the contents of the wells.

Rinse the wells **5 times** with diluted *Wash Solution* (300 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.

8. Add **100 µL of Substrate Solution** into each well.

9. Incubate for **exactly 15 minutes at room temperature (20 °C to 25 °C) in the dark**.

10. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.

Any blue color developed during the incubation turns into yellow.

Note: Highly positive patient samples can cause dark precipitates of the chromogen!

11. Measure the optical density (OD) of the solution in each well at **450 nm filter (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.

It is recommended that the wells be read **within 30 minutes** after adding the Stop Solution.

6.3 Measurement

Measure the absorbance of all wells **at 450 nm** and record the absorbance values for each control and patient sample in the distribution and identification plan.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable **calculate the mean absorbance values** of all duplicates.

7 RESULTS

7.1 Validation of the Test Run

The test run may be considered valid provided the following criteria are met:

Substrate blank in A1:	Absorbance value lower than 0.100
Neg. Control in A1:	Absorbance value lower than 0.200
Cut-off Control in B1/C1:	Absorbance value between 0.350 – 0.850
Pos. Control in D1:	Absorbance value between 0.650 – 3.000

7.2 Calculation of Results

7.2.1 Semiquantitative method

Where applicable **calculate the mean OD values** of all duplicates.

Calculation of Cut-off [CO]

Calculate the mean OD value of the duplicate determination of the Cut-off Control (e.g. in B1/C1).

$$\text{Example: } \frac{0.440 + 0.460}{2} = 0.450 = \text{CO}$$

Calculation of Semiquantitative Results in DRG Units [DU]

$$\frac{\text{Mean OD}_{\text{sample}} \times 10}{\text{CO}} = \text{DRG Units [DU]}$$

$$\text{Example: } \frac{1.580 \times 10 \text{ DU}}{0.450} = 35 \text{ DU}$$

Interpretation of Semiquantitative Results

Cut-off value: 10 DU	NEGATIVE:	< 9 DU
	GREY ZONE (equivocal):	9 – 11 DU
	POSITIVE:	> 11 DU

Interpretation of Qualitative Results

NEGATIVE Patient sample (mean) OD value more than 10 % below CO
(Mean $\text{OD}_{\text{sample}} < 0.9 \times \text{CO}$)

GREY ZONE Patient sample (mean) OD value from 10 % above to 10 % below CO
repeat test 2 - 4 weeks later - with new patient sample
($0.9 \times \text{CO} \leq \text{Mean } \text{OD}_{\text{sample}} \leq 1.1 \times \text{CO}$)
Results in the second test again in the grey zone \Rightarrow **NEGATIVE**

POSITIVE Patient sample (mean) OD value more than 10 % above CO
(Mean $\text{OD}_{\text{sample}} > 1.1 \times \text{CO}$)

8 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Specificity of Antigen

Sample positive to	n	expected Diagnosis	found Diagnosis	Status
Adenovirus IgM	8	negative	negative	pass
Parvovirus IgM	5	negative	negative	pass
CMV IgM	7	negative	negative	pass
Dengue IgM	8	negative	negative	pass
HSV-1 IgM	6	negative	negative	pass
HSV-2 IgG	5	negative	negative	pass
Rubella	5	negative	negative	pass

9.2 Repeatability, Reproducibility

9.2.1 Repeatability (Within-Run)

Determined by 20 measurements of each sample covering the measuring range of the ELISA. CV must be below 10 % for non-negative samples.

The CV must be below 20 % for negative samples in the OD range of 0.1 to 0.5 and below 40 % for negative samples in the OD range below 0.1.

Sample	n	Diagnosis	Mean OD	CV (%)
1	20	negative	0.409	3.7
2	20	negative	0.338	5.7
3	20	negative	0.407	3.7
4	20	low positive	0.862	2.9
5	20	low positive	0.637	3.5
6	20	low positive	0.651	2.6
7	20	positive	1.483	3.9
8	20	positive	1.291	2.2
9	20	positive	1.525	3.5
10	20	high positive	2.223	2.2
11	20	high positive	2.548	1.8
12	20	high positive	1.999	1.9

9.2.2 Reproducibility (Between-Run)

Determined in 20 independent runs with 2 replicates per run. CV must be < 15 %.

Sample	n	Mean OD	CV (%)
1	40	0.938	5.4
2	40	1.472	4.1
3	40	2.305	3.7

9.2.3 Reproducibility (Between-Lot)

The between-lot variation was determined by 6 measurements of different samples with 3 different kit lots. CV must be < 15 %.

Sample	Between-lot agreement (%)
84 negative samples	100
8 positive samples	100

9.3 Diagnostic Sensitivity and Specificity; Predictive Values and Likelihood Ratios

n = 485		Reference Assay	
		positive	negative
DRG ELISA	positive	12	0
	negative	0	473

Diagnostic Sensitivity:	100 %
Diagnostic Specificity:	100 %
Positive Predictive Value:	100 %
Negative Predictive Value:	100 %
Positive Likelihood Ratio:	-1.01
Negative Likelihood Ratio:	-0.99
Percent agreement:	100 %

Using a cut-off value of 10 DU, the DRG RSV IgM ELISA assay showed a diagnostic sensitivity and diagnostic specificity of 100 % and 100 %, respectively. This indicates that a specificity 100 % can be obtained using a cut-off value of 10 DU. It is important to note that samples in this study included apparently healthy donors as well as donors with potential for assay interferences.

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the instructions for use and in compliance with the laboratory quality assurance guidelines. Any improper handling of samples or modification of this test might influence the results. Bacterial contamination or repeated freeze-thaw cycles of the samples may affect the optical density values. In immunocompromised patients and newborns serological data only have restricted value.

10.1 Interfering Substances

10.1.1 Matrix Interference

No interference (bias < ± 20 %) was found for addition of interferent up to the concentration stated in the table below.

Bilirubin	0.5 mg/mL
Hemoglobin	4.0 mg/mL
Triglyceride	30.0 mg/mL

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12 REFERENCES / LITERATUR

1. Wiegand,R.:Respiratory Syncytial Virus (RSV).In:Brandis,H.,W.Köhler,H.J.Eggers,G.Pulverer (ed.): Med. Mikrobiologie 1994.Gustav Fischer Verlag Stuttgart,Jena,New York (1994) 767-769
2. Hohlberg, C.J., et al.: Risk factors for Respiratory Syncytial Virus associated lower respiratory illness in the year of life. Am. J. Epidemiol. 133 (1991) 1135-1151
3. Hall, C.B., et al Immunity to and frequency of reinfection with Respiratory Syncytial Virus. J. Infect. Dis. 163 (1991) 693-698
4. Stein RT et al. Respiratory Syncytial Virus Hospitalization and Mortality: Systematic Review and Meta-Analysis. Pediatric Pulmonology. 2017; 52:556–569
5. Borchers AT et al. Respiratory syncytial virus--a comprehensive review. Rev Allergy Immunol. 2013, 45(3):331-79.

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SHORT INSTRUCTIONS FOR USE

	All reagents and specimens must be allowed to come to room temperature (18-25°C) before use.
	Dispense 100 µL of Controls into appropriate wells.
	Dispense 100 µL of sample into selected wells. (Please note special sample treatment, point 5.3!)
	Cover wells with foil. Incubate for 60 minutes at 37 °C.
	Briskly shake out the contents of the wells.
	Rinse the wells 5 times with diluted Wash Solution (300 µL per well).
	Strike the wells sharply on absorbent paper to remove residual droplets.
	Dispense 100 µL of Enzyme-Conjugate into each well.
	Incubate for 30 minutes at room temperature.
	Briskly shake out the contents of the wells.
	Rinse the wells 5 times with diluted Wash Solution (300 µL per well).
	Strike the wells sharply on absorbent paper to remove residual droplets.
	Add 100 µL of Substrate Solution to each well.
	Incubate for 15 minutes at room temperature.
	Stop the reaction by adding 100 µL of Stop Solution to each well.
	Determine the absorbance of each well at 450 nm.