Instructions for Use

Toxoplasma gondii IgG ELISA

IVD

CE₀₁₉₇

REF EIA-3519

Σ 96

DRG 

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG 

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.
Utilize apenas a versão válida das Instruções de Utilização fornecidas com o kit.

Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées / Modificações introduzidas	
The following changes have been made in comparison to the previous version: Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen: Rispetto alla versione precedente, sono state apportate le seguenti modifiche: Se han introducido los siguientes cambios en comparación con la versión anterior: Les modifications suivantes ont été apportées par rapport à la version précédente : Foram efetuadas as seguintes alterações em comparação com a versão anterior:	
Detailed editorial revision. Changed wording in several chapters. Ausführliche redaktionelle Überarbeitung. Geänderter Wortlaut in mehreren Kapiteln. Revisione editoriale dettagliata. Modificato il testo in diversi capitoli. Revisión editorial detallada. Se ha cambiado la redacción de algunos capítulos. Révision éditoriale détaillée. Modification de la formulation dans plusieurs chapitres. Revisão editorial detalhada. Texto alterado em vários capítulos.	
7.2 Calculation of Quantitative Results	Changed wording for better understanding. „Values” replaced by „Quantitative Results”

ENGLISH	2
DEUTSCH.....	12
LITERATURE	21
SYMBOLS USED.....	22

1 INTENDED USE

The **DRG Toxoplasma gondii IgG ELISA** is a manual enzyme immunoassay for the **quantitative** and **qualitative** measurement of IgG-class antibodies to *Toxoplasma gondii* in human serum or plasma (EDTA, Li-heparin or citrate plasma).

For *in vitro* diagnostic use. For laboratory professional use.

1.1 Scientific Validity Report

Toxoplasma gondii is a small intracellular parasite, whose life cycle has a sexual and an asexual phase. Sexual development is restricted to the intestinal cells of (probably exclusively) cats; the oocysts formed are excreted and due to their resistant cell walls, they may be infectious under advantageous circumstances for at least 1 year. Animals and man are intermediate hosts for the asexual proliferation of *T. gondii*: the ingested parasites will proliferate explosively within the host cells lysing them eventually. They disseminate throughout the body via circulation and lymphatic system and may infect any cell type.

About half of the *Toxoplasma* infections process clinically unapparent. The other *Toxoplasma* infections show often only unspecific symptoms, after an incubation time from one to three weeks with lightly fever, exhaustion, headaches well as muscle and joint pain. Complications in the form of myocarditis, meningitis or pneumonia occur at 1% of infected children and young adults respectively. After recovery, *Toxoplasma gondii* cells persists in infected tissues by forming cysts which are resistant to the attacks of the immune system. In immunocompetent individuals the latent infection does not reactivate. After suppression of the immune system, activation of latent infections has been observed and can lead to severe complications.

An infection at any stage during pregnancy may result in the transplacental transmission of the parasite to the fetus although the timing of such a transfer does have an effect on the clinical outcome for the fetus.

1. Trimester: 17% → most often aborted, seldom severe damage to the unborn
2. Trimester: 24% → moderate or severe damage to the fetus
3. Trimester: 64% → mild damage or damage appears later in life

When maternal primary infection is detected and subsequently eradicated by chemotherapy, the risk for transmission to the fetus is decreased by about 75%. In AIDS patients, *Toxoplasma* encephalitis is of significant importance as the cause of death.

Interpretation of Results:

Negative results for *Toxoplasma gondii* IgG by ELISA indicate a lack of immunity but do not rule out completely an active infection. Consequently, an additional test for IgM antibodies is recommended.

Seronegative pregnant women should be tested in intervals of max. 12 weeks.

Detection of *Toxoplasma gondii* IgA as additional parameter does have additional, but no exclusive significance. At a detection of *Toxoplasma gondii* IgM, especially with simultaneous high or rising IgG concentration, the possibility of an acute first infection has to be considered. The detection of IgM in comparison with the detection of IgA delivers more certain decision between acute or latent infection. But the obtained result for IgA shows a very good additional decision support, if a suspicion for an acute infection occurs but no high IgM result is found.

Results should be interpreted according to the table below (Robert Koch Institute, Germany)

IgG	IgM	IgG-Avidity	Probable result
positive	negative	-----	Inactive, latent infection
positive	positive	high	Subsiding or latent (inactive) infection
positive	positive	low	Acute infection possible, therefore further clarification methods /clinical monitoring necessary

Infection may be identified by PCR, indirect immunofluorescence (IIF), Serology: detection of antibody production by ELISA.

2 PRINCIPLE OF THE TEST

The **DRG Toxoplasma gondii IgG ELISA** Kit is a solid phase enzyme-linked immunosorbent assay (ELISA)

Microtiter wells as a solid phase are coated with inactivated *Toxoplasma gondii* soluble antigen (strain RH).

Diluted patient specimens and **ready-for-use controls** are pipetted into these wells. During incubation *Toxoplasma gondii*-specific antibodies of positive specimens and controls are bound to the immobilized antigens.

After a washing step to remove unbound sample and control material horseradish peroxidase conjugated anti-human IgG antibodies are dispensed into the wells. During a second incubation this anti-IgG conjugate binds specifically to IgG antibodies resulting in the formation of enzyme-linked immune complexes.

After a second washing step to remove unbound conjugate the immune complexes formed (in case of positive results) are detected by incubation with TMB substrate and development of a blue colour. The blue colour turns into yellow by stopping the enzymatic indicator reaction with sulfuric acid.

The intensity of this colour is directly proportional to the amount of Toxoplasma gondii-specific IgG antibody in the patient specimen. Optical density at 450 nm is read using an ELISA microtiter plate reader.

3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only. For laboratory professional use only.
- Before starting the assay, read the instructions for use completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Do not reuse microtiter wells.
- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact DRG.
- Microbial contamination of reagents or samples may give false results.
- Allow the reagents to reach room temperature (20 °C to 25 °C) before starting the test. Temperature will affect the optical density readings of the assay.
- All indicated volumes must be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into original vials as reagent contamination may occur.

General Precautions

- Follow laboratory quality assurance and laboratory safety guidelines
- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.

Biohazard Information

- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. However, no known test method can offer total assurance that no infectious agent is present.
- The device contains material of animal origin, which is certified apparently free of infectious or contagious diseases and injurious parasites.
- Bovine components originate from countries where BSE (Bovine spongiform encephalopathy) has not been reported.
- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation. Waste must be discarded according to local rules and regulations.

Information to Chemical Hazards and Hazard Classification

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Substrate Solution contains an ingredient in non-declarable concentrations which causes serious eye irritation. In case of possible contact with eyes, rinse immediately carefully and thoroughly with eye wash or water. After contact with skin, wash with plenty of water. Take-off contaminated clothing and wash it before reuse.
- Avoid contact with Stop Solution containing < 2 % H₂SO₄. It may cause skin irritation and burns.
- Chemicals and prepared or used reagents must be treated as hazardous waste according to the national safety guideline or regulation.
- This product does not contain substances which have carcinogenic, mutagenic or toxic for reproduction (CMR) properties.

All reagents of this test kit do NOT contain hazardous substances in concentrations to be declared, a classification and labelling is not required.

For detailed information, please refer to the Safety Data Sheet, which is available upon request directly from DRG.

4 MATERIALS

4.1 Materials Provided with the Kit

Symbol	Quantity	Description	Preparation
Microtiterwells	12 x 8 wells (break apart)	Microtiter plate Wells coated with inactivated Toxoplasma gondii soluble antigen (strain RH).	Ready to use
Sample Diluent	1 x 100 mL	Sample Diluent * colored yellow, pH 7.2 ± 0.2	Ready to use
Standards (Standard 1 – 3)	3 x 1 vial, S1: 2 mL S2: 1 mL S3: 1 mL	Standards * <i>Calibrated against the following reference material: WHO International Standard 3rd IS for Anti-Toxoplasma Serum Ig (NIBSC code: TOXM)</i> Standard 1: 50 IU/mL, colored green, green cap. Standard 2: 100 IU/mL, colored blue, blue cap. Standard 3: 200 IU/mL, colored red, red cap.	Ready to use
Neg. Control	1 x 2 mL	Negative Control * Colored yellow, yellow cap.	Ready to use
Pos. Control	1 x 1 mL	Positive Control * colored purple, black cap.	Ready to use
Enzyme Conjugate	1 x 20 mL	Enzyme Conjugate * Colored red, antibody to human IgG conjugated to horseradish peroxidase.	Ready to use
Substrate Solution	1 x 14 mL	Substrate Solution Contains 3,3',5,5'-tetramethylbenzidine (TMB). <i>Keep away from direct sun light.</i>	Ready to use
Stop Solution	1 x 14 mL	Stop Solution Contains < 2 % H ₂ SO ₄ . <i>Avoid contact with the stop solution. It may cause skin irritations and burns.</i>	Ready to use
Wash Solution	1 x 30 mL	Wash Solution, 20X concentrate *	See "Reagent Preparation".
	1 x	Cover foil	
	1 x	Instructions for Use	
	1 x	Certificate of Analysis (CoA)	
		* Contain(s) non-mercury preservative.	

4.2 Materials Required but Not Provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Incubator for 37 °C
- Manual or automatic equipment for rinsing microtiter plate wells
- Vortex tube mixer
- Absorbent paper
- Freshly distilled water
- Timer
- Graph paper or software for data reduction

4.3 Storage and Stability of the Kit

Unopened kits and reagents as well as **opened reagents** must be stored at 2 °C to 8 °C.

The microplate must always be stored in the resealable aluminum pouch containing a desiccant. Do not open the pouch until it has reached room temperature. The microtiter plate consists of 12 individual strips. Each strip can be divided into 8 individual wells.

Unused wells must be immediately returned to the aluminum pouch with the desiccant and stored again tightly resealed at 2 °C to 8 °C.

Once opened, reagent vials must be closed tightly again.

	Storage Temperature	Stability
Unopened kits and unopened reagents	2 °C to 8 °C	Until the expiration date printed on the label. Do not use reagents beyond this date!
Opened kit	2 °C to 8 °C	8 weeks

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (RT, 20 °C to 25 °C) prior to use.

Wash Solution

Add fresh and germ-free distilled water to the 20X concentrated Wash Solution.

Dilute the complete content of the vial **1 + 19** (30 mL *Wash Solution* + 570 mL distilled water) to a final volume of 600 mL. If crystals have formed in the wash solution concentrate, ensure that they are completely transferred and dissolved into the solution.

This diluted wash solution must have a pH value of 7.2 ± 0.2 .

Stability after dilution:	at 2 °C to 8 °C	1 week
---------------------------	-----------------	--------

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components must not be used for a test run. They have to be stored until a final solution has been found. After this, they must be disposed of according to the official regulations.

5 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

Human serum or plasma (EDTA plasma, lithium heparin plasma or citrate plasma)

Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use hemolytic, icteric, or lipemic samples. For further information refer to chapter "Interfering Substances".

5.1 Sample Collection

Serum: Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma: Whole blood should be collected into centrifuge tubes containing anticoagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

Whole blood should not be frozen before centrifugation.

5.2 Sample Storage

Samples must be stored tightly capped prior to performing the assay. If stored frozen, freeze only once. Thawed samples must be inverted several times prior to testing.

Stability	at 2 °C to 8 °C	5 days
	at -20 °C (in aliquots)	up to 18 months

5.3 Sample Preparation

Before starting the test, dilute each sample **1 + 100** with *Sample Diluent*, e.g. 10 µL sample + 1 mL *Sample Diluent mix well, let stand for 15 minutes and mix well before use.*

Take **100 µL** of the prediluted sample for the ELISA. Use only fresh diluted samples!

For patients with concentrations greater than Standard 3 a second 1:10 dilution of this 1+100 diluted patient sample should be performed;

e.g. 30 µL of first sample dilution + 270 µL Sample Diluent (mix well).

Please note: Standards and Controls are ready to use and must not be diluted!

6 ASSAY PROCEDURE

6.1 Procedural Notes

- All reagents and samples must be allowed to come to room temperature (20 °C to 25 °C) before use.
- All reagents must be mixed without foaming.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- Use new disposal plastic pipette tips for each standard, control, or sample in order to avoid carry-over.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Mix the contents of the microtiter plate wells thoroughly to ensure good test results.
- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Once the test has been started, all steps must be completed without interruption and in the same sequence for each step.
- The enzymatic reaction is linearly proportional to time and temperature.
- Optical density is a function of the incubation time and temperature. Respect the incubations times and temperatures as given in chapter "Test Procedure".
- Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- During the incubation at 37 °C cover microtiter strips with foil to avoid evaporation.
- **Important note to wash procedure:**
Washing is critical. Improperly washed wells will give erroneous results. The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
- **Test performance using fully automated analysis devices:**
Automated test performance using fully automated, open-system analysis devices is possible. However, the combination must be validated by the user.

6.2 Test Procedure

Each run must include a standard curve.

Before starting the assay, dilute *Wash Solution*, **prepare samples as described in section 5.3** and mix well before pipetting. For all standards/controls and samples carefully document the distribution on the plate on the Plate Layout Sheet supplied in the kit.

1. Secure the desired number of microtiter wells in the frame holder.

Please allocate at least:

1 well (e.g. A1) for the *Neg. Control*
4 wells (e.g. B1 to E1) for the *Standard 1 – 3*
1 well (e.g. F1) for the *Pos. Control*.

It is left to the user to determine standards, controls and samples in duplicate.

2. Pipette

100 µL of *Neg. Control* into well A1
100 µL of *Standard 1* into well B1 + C1
100 µL of *Standard 2* into well D1
100 µL of *Standard 3* into well E1
100 µL of *Pos. Control* into well F1 and
100 µL of each diluted sample with new disposable tips into appropriate wells.

3. Cover wells with the foil supplied in the kit. Incubate for **60 minutes at 37 °C**.

4. Wash the wells as follows:

If the wash step is performed manually:

Briskly shake out the contents of the wells.

Rinse the wells **5 times** with **300 µL** diluted *Wash Solution* per well.

If an automated plate washer is used:

Rinse the wells **5 times** with **400 µL** diluted *Wash Solution* per well.

At the end of the washing step, always strike the wells sharply on absorbent paper to remove residual droplets!

5. Pipette **100 µL Enzyme Conjugate** into each well.

6. Incubate for **30 minutes at room temperature (20 °C to 25 °C)**.

7. Wash as described in step 4.

8. Pipette **100 µL of Substrate Solution** into each well.

9. Incubate for exactly **15 minutes at room temperature (20 °C to 25 °C) in the dark**.

10. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.

Any blue color developed during the incubation turns into yellow.

Note: Highly positive samples can cause dark precipitates of the chromogen!

11. Measure the optical density (OD) of the solution in each well at **450 nm (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.

It is recommended that the wells be read **within 30 minutes** after adding the *Stop Solution*.

6.3 Measurement

Measure the optical density (OD) of all wells at **450 nm** and record the absorbance values for each control and patient sample in the Plate Layout Sheet.

Dual wavelength reading using 620 nm as reference wavelength is recommended.

Where applicable **calculate the mean OD values** of all duplicates.

7 RESULTS

The results alone should not be the only reason for any therapeutic consequences. The results must be correlated to other clinical observations and diagnostic tests.

7.1 Validation of the Test Run

The test run may be considered valid provided the following criteria are met:

Neg. Control in A1:	OD value lower than 0.200
Standard 1 (Cut-off) in B1 + C1:	OD value between 0.350 – 0.850
Standard 2 in D1:	OD value between 0.750 – 1.500
Standard 3 in E1:	OD value between 1.000 – 2.000
Pos. Control in F1:	OD value between 0.650 – 3.000

(OD Standard 1 < OD Standard 2 < OD Standard 3).

7.2 Calculation of Quantitative Results

In order to obtain **quantitative results in IU/mL** plot the (mean) OD values of *Neg. Control* and *Standard 1, 2, and 3* on (linear/linear) graph paper in a system of coordinates against their corresponding concentrations (0, 50, 100 and 200 IU/mL) and draw a standard calibration curve (OD values on the vertical y-axis, concentrations on the horizontal x-axis). Read results from this standard curve employing the (mean) OD values of each patient specimen and control.

All validated computer programs can be used, that provide the following function: 4 PL (4 Parameter Logistics).

We use DRG regression program for Windows (4 parameter Rodbart regression).

If other regression software is used, the obtained values have to be validated by the user.

NOTE: Quantitative Results of additionally (1:10, in total 1:1000) diluted patient samples must be multiplied by the appropriate dilution factor in order to obtain correct results! (Dilution: 1:10 = Dilution factor: 10). (See chapter "5.3 Specimen Dilution").

7.3 Interpretation of quantitative Results

The following values should be considered as a guideline:

Cut-off value: 50 IU/mL	NEGATIVE:	< 45	IU/mL
	GREY ZONE(equivocal):	45 – 55	IU/mL
	POSITIVE:	> 55	IU/mL

7.4 Calculation of qualitative Results

7.4.1 Calculation of Cut-off (CO)

Standard 1 is used as Cut-off!

Calculate the mean OD value of the duplicate determination of *Standard 1* (e.g. in B1/C1).

Example: $\frac{0.480 + 0.460}{2} = 0.470 = \text{CO}$

7.5 Interpretation of qualitative Results

NEGATIVE	Mean OD_{patient} < OD_{CO} - 10%
GREY ZONE	OD_{CO} - 10% ≤ Mean OD_{patient} ≤ OD_{CO} + 10% Repeat test 2 - 4 weeks later - with new patient samples Results in the second test again in the grey zone → NEGATIVE
POSITIVE	Mean OD_{patient} > OD_{CO} + 10%

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each standard curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the Quality Control Laboratory are stated in the Certificate of Analyses (CoA) added to the kit. The values and ranges stated on the CoA always refer to the current kit lot and must be used for direct comparison of the results.

If available, it is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above-mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 3.45 – 200 IU/mL.

9.2 Specificity of Antigen

No cross reactivity was found for *Entamoeba histolytica*, *Giardia lamblia*, *Schistosoma mansoni*, *Toxocara canis*, *Strongyloides* and *Echinococcus*.

9.3 Analytical Sensitivity

The analytical sensitivity of the DRG ELISA was calculated by adding 2 standard deviations to the mean OD of 20 replicate analyses of the negative control and was found to be 3.45 IU/mL (OD_{450nm} 0.056).

9.4 Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. (Detected by method comparison with Diamed Eurogen ELISA with three lots of DRG ELISA, 125 samples, therefrom 59 negative samples are assayed.)

It is 98.3% (for all three DRG production lots).

9.5 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. (Detected by method comparison with Diamed Eurogen ELISA with three lots of DRG ELISA, 125 samples, therefrom 66 positive samples are assayed.)

It is 100% (for all three DRG production lots).

9.6 Method Comparison

The DRG ELISA was compared with another *Toxoplasma gondii* IgG ELISA (Diamed Eurogen). 125 serum samples are assayed.

		Other commercially available CE-marked ELISA	
		positive	negative
DRG ELISA	positive	66	1
	negative	0	58

Agreement: 99.20%

9.7 Reproducibility

9.7.1 Intra-assay

The intra-assay (within-run) precision of the DRG Toxoplasma gondii IgG ELISA was determined by 20 measurements of 12 serum samples covering the whole measuring range.

Sample	Mean Conc. (IU/mL)	Intra-Assay CV (%)	n
1	17.93	9.8	20
2	33.80	5.2	20
3	26.25	7.9	20
4	53.29	6.3	20
5	55.21	4.6	20
6	65.34	3.2	20
7	95.23	5.9	20
8	100.15	6.0	20
9	97.82	8.0	20
10	118.89	5.8	20
11	196.32	6.7	20
12	193.35	4.8	20

9.7.2 Inter-assay

The inter-assay variation of the DRG Toxoplasma gondii IgG ELISA was determined with 3 samples with 2 production kits in 10 independent runs with 2 replicates per run.

Sample	Mean Conc. (IU/mL)	Inter-Assay CV (%)	n
1	112.02	6.8	40
2	138.63	5.8	40
3	211.45	5.0	40

9.8 Recovery

Samples have been spiked by adding 3 solutions with known concentrations in a 1:1 ratio.

The % recovery has been calculated by multiplication of the ratio of the measurements and the expected values with 100 (expected value = (endogenous value + added value) / 2; because of a 1:2 dilution of serum with spike material).

	Sample 1	Sample 2	Sample 3
Concentration (IU/mL)	27.00	32.14	37.99
Average Recovery (%)	96.0	97.5	94.2
Range of Recovery (%)	from	87.0	88.9
	to	101.6	104.4

9.9 Linearity

Three samples (serum) containing different amounts of analyte were serially diluted with sample diluent and assayed with the DRG ELISA. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

	Sample 1	Sample 2	Sample 3
Concentration (IU/mL)	200.00	280.00	403.00
Average Recovery (%)	103.6	95.9	103.2
Range of Recovery (%)	from	99.3	88.8
	to	106.8	103.9

11 LIMITATIONS OF THE PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

Bacterial contamination or repeated freeze-thaw cycles of the samples may affect the optical density values.

In immunocompromised patients and newborns serological data only have restricted value.

11.1 Interfering Substances

Hemoglobin (up to 4.0 mg/mL), bilirubin (up to 0.5 mg/mL) and triglyceride (up to 30.0 mg/mL) have no influence on the assay results.

12 LEGAL ASPECTS

12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the laboratory quality assurance guidelines and applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact DRG.

12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12.4 Reporting of Serious Incident

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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13 LITERATURE

1. P. Cassinotti, Human *Toxoplasma gondii* infections and their diagnosis
Alpe Adria Microbiology Journal 1995, 4, 235-246.
2. M. Schleuning, *Toxoplasma gondii* Infektionen
Deutsches Ärzteblatt 93, Heft 43, (Oktober 1996), B2182-B2185.
3. NCCLS. Clinical Use and Interpretation of Serologic Tests for *Toxoplasma gondii*; Approved Guideline.
NCCLS document M36-A [ISBN 1-56238-523-2]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
4. Wilson M, Jones JL, McAuley JM. *Toxoplasma*.
In: Murray PR, Baron EJ, Pfaller MA, Jorgensen JH, Tenover FC, Tenover FC, editors. Manual of Clinical Microbiology. 8th ed. Washington, D.C.: American Society for Microbiology; 2003. p. 1970-1980.
5. Remington JS, McLeod R, Thulliez P, Desmonts G. Toxoplasmosis.
In: Remington JS, Klein JO, editors. Infectious Diseases of the Fetus and Newborn Infant. 5th ed. Philadelphia, PA: The WB Saunders Co.; 2001. p. 205-346.
6. Robert-Koch-Institut (2009) Toxoplasmose, RKI-Ratgeber Infektionskrankheiten – Merkblätter für Ärzte.
7. Ybañez RHD, Ybanez AP and Nishikawa Y. Review on the Current Trends of Toxoplasmosis Serodiagnosis in Humans.
Front Cell Infect Microbiol. 2020; 10: 204.
8. Nayeri T et al. The global seroprevalence of anti-*Toxoplasma gondii* antibodies in women who had spontaneous abortion: A systematic review and meta-analysis
PLoS Negl Trop Dis. Published online 2020 Mar 13.
9. Dasa TT et al. Toxoplasmosis infection among pregnant women in Africa: A systematic review and meta-analysis.
PLoS One. 2021; 16(7), 1-14.
10. Pleyer U, Gross U, Schlüter D, Wilking H, Seeber F: Toxoplasmosis in Germany—epidemiology, diagnosis, risk factors, and treatment.
Dtsch Arztebl Int 2019; 116: 435–44.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français	Português
	European Conformity	CE-Konformitätskennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes	Conformidade Europeia
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter la notice d'utilisation	Consultar as instruções de uso
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico in vitro	Dispositif médical de diagnostic in vitro	Dispositivo médico para diagnóstico in vitro
	Catalogue number *	Artikelnummer *	No. di Cat.	No de catálogo	Référence	Número de catálogo
	Batch code *	Fertigungslosnummer, Charge *	Lotto no	Número de lote	No. de lot	Código do lote
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests	Suficiente para <n> determinações
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation	Limites de temperatura
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation	Prazo de validade
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant	Fabricante
	Date of Manufacture *	Herstellungsdatum *	Data di produzione	Fecha de fabricación	Date de production	Data de fabricação
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques	Riscos biológicos
	Caution *	Achtung *	Attenzione	Precaución	Attention	Cuidado
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches	
	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur	Distribuidor
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement	Conteúdo
<i>Volume/No.</i>	Volume / No.	Volumen / Anzahl	Volume / Quantità	Volumen / Número	Volume / Quantité	Volume / Quantidade
<i>Microtiterwells</i>	Microtiter plate	Mikrotiterplatte	Piastra per microtitolazione	Placa de microtítulo	Microplaque	Placa de microtitulação
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Coniugato enzimatico	Conjugado enzimático	Conjugué enzymatique	Conjugado Enzimático
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution de substrat	Solução de Substrato
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d'arresto	Solución de parada	Solution d'arrêt	Solução de Paragem
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard	Padrão
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle	Controlo
<i>Pos. Control</i>	Positive Control	Positive Kontrolle	Controllo positivo	Control positivo	Positif Contrôle	
<i>Neg. Control</i>	Negative Control	Negative Kontrolle	Controllo negativo	Control negativo	Négatif Contrôle	
<i>Cut-off Control</i>	Cut-off Control	Grenzwert-Kontrolle	Controllo valore limite	Control valor limite	Valeur limite Contrôle	
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage	Solução de Lavagem
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungsmedium	Diluyente per campioni	Diluyente de muestras	Diluant d'échantillon	Diluyente de Amostra
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungsmedium	Diluyente coniugato	Diluyente de conjugados	Diluant de conjugué	Diluyente de Conjugado