Instructions for Use

Toxocara canis IgG ELISA

IVD

CE

REF EIA-3865

Σ 96

DRG Instruments GmbH, Germany
Frauenbergstraße 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

Please use only the valid version of the Instructions for Use provided with the kit.

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1 INTENDED USE

The Toxocara canis IgG ELISA is intended for the qualitative determination of IgG class antibodies against Toxocara canis in human serum or plasma (citrate, heparin).

2 PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

3 MATERIALS

3.1 Reagents supplied

- **Microtiterplate:**
12 break-apart 8-well snap-off strips coated with Toxocara canis antigens; in resealable aluminium foil.
- **DIL:**
1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2 ; coloured yellow; ready to use; white cap; $\leq 0.0015\% \text{ (v/v)}$ CMIT/ MIT (3:1).
- **SOLN|STOP:**
1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH|BUF|20x:**
1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2 , for washing the wells; white cap; 0.2% (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Conjugate:**
1 bottle containing 20 mL of peroxidase Protein A in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **SUB|TMB:**
1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:**
1 vial containing 2 mL control; coloured yellow; ready to use; red cap; $\leq 0.02\% \text{ (v/v)}$ MIT.
- **Cut-off Control:**
1 vial containing 3 mL control; coloured yellow; ready to use; green cap; $\leq 0.02\% \text{ (v/v)}$ MIT.
- **Negative Control:**
1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; $\leq 0.0015\% \text{ (v/v)}$ CMIT/ MIT (3:1).

For hazard and precautionary statements see 11.1

3.2 Materials supplied

- 1 Cover foil
- 1 instructions for use (IFU)

3.3 Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37°C
- Manual or automatic equipment for rinsing Microtiterplate
- Pipettes to deliver volumes between 10 µL and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

4 STABILITY AND STORAGE

Store the kit at 2 °C - 8 °C.

The opened reagents are stable up to the expiry date stated on the label when stored at 2 °C - 8 °C.

5 REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20 °C - 25 °C) and mix them before starting the test run!

5.1 Microtiterplate

The break-apart snap-off strips are coated with Toxocara canis antigens.

Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2 °C - 8 °C.

5.2 **[WASH | BUF | 20x]**

Dilute **[WASH | BUF | 20x]** 1 + 19; e. g. 10 mL **[WASH | BUF | 20x]** + 190 mL distilled water.

The diluted buffer (**[WASH | BUF | 1x]**) is stable for 5 days at room temperature (20 °C - 25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

5.3 **[SUB | TMB]**

The reagent is ready to use and has to be stored at 2 °C - 8 °C, away from the light.

[SUB | TMB] should be colourless or could have a slight blue tinge. If **[SUB | TMB]** turns into blue, it may have become contaminated and should be thrown away.

6 SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay.

If the assay is performed within 5 days after sample collection, the samples should be kept at 2 °C - 8 °C; otherwise they should be aliquoted and stored deep-frozen (-70 °C to -20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

6.1 Sample Dilution

Before assaying, all samples should be diluted **1 + 100** with **[DIL]**.

Dispense 10 µL sample and 1 mL **[DIL]** into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

7 ASSAY PROCEDURE

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems, we recommend increasing the washing steps from three up to five and the volume of [WASH] [BUF] 1x from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 °C ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells.
Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 °C ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of [WASH] [BUF] 1x. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20 °C - 25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL [SUB] TMB into all wells.
9. **Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL [SOLN] STOP into all wells in the same order and at the same rate as for the [SUB] TMB, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the [SOLN] STOP.

7.1 Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

8 RESULTS

8.1 Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < 0.100
- **Negative Control:** Absorbance value < 0.200 and < Cut-off
- **Cut-off Control:** Absorbance value 0.150 - 1.300
- **Positive Control:** Absorbance value > Cut-off

If these criteria are not met, the test is not valid and must be repeated.

8.2 Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86/2 = 0.43
Cut-off = 0.43

8.2.1 Results in Units [DU]

Sample (mean) absorbance value × 10 = [DRG Units = DU]
Cut-off

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ DU}$

8.3 Interpretation of Results

Cut-off	10 DU	-
Positive	> 11 DU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 - 11 DU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks.
Negative	< 9 DU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

9 SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

9.1 Precision

Intra assay	n	Mean (E)	CV (%)
#1	24	0.516	3.61
#2	24	0.904	2.90
#3	24	1.061	2.20
Inter assay	n	Mean (DU)	CV (%)
#1	12	18.81	4.83
#2	12	21.65	7.01
#3	12	2.23	12.17

9.2 Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 98.63% (95% confidence interval: 96.52% - 99.62%).

9.3 Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 96.92% (95% confidence interval: 89.32% - 99.63%).

9.4 Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

9.5 Cross Reactivity

Cross reactivity with antibodies against Ascaris lumbricoides and Schistosoma cannot be excluded.

10 LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

11 PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

11.1 Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning	H317	May cause an allergic skin reaction.
	P261	Avoid breathing spray.
	P280	Wear protective gloves/ protective clothing.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water.
	P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
	P362+P364	Take off contaminated clothing and wash it before reuse.

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Therefore, the following hazard and precautionary statements apply.

Warning	H315	Causes skin irritation.
	H319	Causes serious eye irritation
	P280	Wear protective gloves/ protective clothing.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water.
	P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313

If eye irritation persists: Get medical advice/attention.

Further information can be found in the safety data sheet.

11.2 Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
																		
PAP 21	PAP 21	PAP 22																
<table border="1" style="margin-left: auto; margin-right: auto;"> <tr> <td>SOLN</td><td>STOP</td><td>WASH</td><td>BUF</td><td>20x</td><td>SUB</td><td>TMB</td><td>DIL</td></tr> <tr> <td>CONJ</td><td></td><td>CONTROL</td><td>+</td><td></td><td>CONTROL</td><td>-</td><td>CUT OFF</td></tr> </table>	SOLN	STOP	WASH	BUF	20x	SUB	TMB	DIL	CONJ		CONTROL	+		CONTROL	-	CUT OFF	MTP	
SOLN	STOP	WASH	BUF	20x	SUB	TMB	DIL											
CONJ		CONTROL	+		CONTROL	-	CUT OFF											
		PET / ALU / LDPE 90																
HDPE 2	PP 5																	

12 SCHEME OF THE ASSAY

Test Preparation

Prepare reagents and samples as described.

Establish the distribution and identification plan for all samples and standards/controls.

Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate Blank (A1)	Negative Control	Cut-off Control	Positive Control	Sample (diluted 1+100)
Negative Control	-	100 µL	-	-	-
Cut-off Control	-	-	100 µL	-	-
Positive Control	-	-	-	100 µL	-
Sample (diluted 1+100)	-	-	-	-	100 µL
Cover wells with foil supplied in the kit					
Incubate for 1 h at 37 °C ± 1 °C					
Wash each well three times with 300 µL of WASH BUF 1x					
Conjugate	-	100 µL	100 µL	100 µL	100 µL
Incubate for 30 min at room temperature (20 °C - 25 °C)					
Do not expose to direct sunlight					
Wash each well three times with 300 µL of WASH BUF 1x					
SUB TMB	100 µL	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature (20 °C - 25 °C) in the dark					
SOLN STOP	100 µL	100 µL	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	Numero di Catalogo	Référence de catalogue
	Batch code *	Chargencode *	Codice del lotto	Numéro de lot
	Unique Device Identifier	Eindeutige Produktidentifizierung	identificazione unica del dispositivo	identification unique des dispositifs
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricant
	Caution *	Achtung *		
RUO	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volume/Quantité
MTP	Microtiterplate	Mikrotiterplatte	Piastre di Microtitolazione	Plaque de Microtitrage
CONJ	Conjugate	Konjugat	Coniugato	Conjugué
CONTROL -	Negative Control	Negativkontrolle	Controllo negativo	Contrôle négatif
CONTROL +	Positive Control	Positivkontrolle	Controllo positivo	Contrôle positif
CUT OFF	Cut off control	Cut off-Kontrolle	Controllo cut-off	Contrôle cut-off
CAL	Standard or Calibrator	Standard oder Kalibrator	Standard o Calibratore	Standard o Etalon
DIL	Sample Dilution Buffer	Probenverdünnungspuffer	Tampon diluente	Tampon de dilution d'échantillon
DIL A	IgA Sample Dilution Buffer	IgA-Proben-verdünnungspuffer	Tampon diluente IgA	Tampon de dilution d'échantillon IgA
DIL G	IgG Sample Dilution Buffer	IgG-Proben-verdünnungspuffer	Tampon diluente IgG	Tampon de dilution d'échantillon IgG
DIL M	IgM Sample Dilution Buffer	IgM-Proben-verdünnungspuffer	Tampon diluente IgM	Tampon de dilution d'échantillon IgM
SOLN STOP	Stop solution	Stoplösung	Soluzione bloccante	Solution d'arrêt
SUB TMB	TMB Substrate solution	TMB-Substratlösung	Soluzione substrato TMB	Solution de substrat TMB
WASH BUF 20x	Washing Buffer 20x concentrated	Waschpuffer 20x konzentriert	Tampon de lavage concentré 20 x	Tampon de lavage concentré 20 x
WASH BUF 1x	20-fold dilution of WASH BUF 20x	20-fach Verdünnung von WASH BUF 20x	Diluizione 20 volte del WASH BUF 20x	Dilution 20 fois du WASH BUF 20x