Instructions for Use

Anti-LKM-1 ELISA

IVD

REF EIA-4277

Σ 96

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.
Utilize apenas a versão válida das Instruções de Utilização fornecidas com o kit.

**Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas /
 Modifications apportées / Modificações introduzidas**

The following changes have been made in comparison to the previous version:
 Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen:
 Rispetto alla versione precedente, sono state apportate le seguenti modifiche:
 Se han introducido los siguientes cambios en comparación con la versión anterior:
 Les modifications suivantes ont été apportées par rapport à la version précédente :
 Foram efetuadas as seguintes alterações em comparação com a versão anterior:

Detailed editorial revision. Changed wording in several chapters.
Ausführliche redaktionelle Überarbeitung. Geänderter Wortlaut in mehreren Kapiteln.
Revisione editoriale dettagliata. Modificato il testo in diversi capitoli.
Revisión editorial detallada. Se ha cambiado la redacción de algunos capítulos.
Révision éditoriale détaillée. Modification de la formulation dans plusieurs chapitres.
Revisão editorial detalhada. Texto alterado em vários capítulos.

Table of Contents / Inhaltsverzeichnis / Tabla de Contenidos

1	INTENDED PURPOSE	3
2	DIAGNOSTIC RELEVANCE	3
3	TEST PRINCIPLE	3
4	TEST COMPONENTS	4
5	MATERIALS REQUIRED BUT NOT PROVIDED	4
6	STORAGE AND STABILITY	4
7	GENERAL INFORMATION	4
8	PREPARATION.....	5
9	TEST PERFORMANCE	6
10	TEST EVALUATION	6
11	PERFORMANCE CHARACTERISTICS.....	7
12	WARNINGS AND PRECAUTIONS	8
13	DISPOSAL	8
14	REFERENCES.....	8
1	ZWECKBESTIMMUNG	9
2	DIAGNOSTISCHE BEDEUTUNG	9
3	TESTPRINZIP	9
4	TESTKOMPONENTEN	10
6	LAGERUNG UND HALTBARKEIT	10
7	ALLGEMEINE HINWEISE.....	10
8	VORBEREITUNG.....	11
9	DURCHFÜHRUNG	12
10	AUSWERTUNG	12
11	LEISTUNGSMERKMALE.....	13
12	WARNHINWEISE UND VORSICHTSMAßNAHMEN	14
13	ENTSORGUNG.....	14
14	REFERENZEN	14

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1 INTENDED PURPOSE

The Anti-LKM-1 is a quantitative immunoassay for the determination of IgG antibodies against cytochrome p450 IID6 in human serum.

The Anti-LKM-1 is intended as an aid in the diagnosis of autoimmune hepatitis (AIH) in conjunction with other clinical and laboratory findings.

The immunoassay is designed for manual professional *in vitro* diagnostic use.

2 DIAGNOSTIC RELEVANCE

The group of primary autoimmune liver disease (PAL) comprises AIH, primary biliary cirrhosis (PBC) and primary sclerosing cholangitis (PSC).

The clinical picture of PAL is in most cases not different from other chronic liver diseases. About 15% of all cases with chronic liver diseases show an autoimmune pathogenesis. Therefore, after exclusion of infectious etiology especially by viruses, the determination of different autoantibodies is recommended.

Patients suffering from AIH show a variety of autoantibodies. Due to the appearance of different antibody specificities classification of AIH into different subgroups is discussed. Type I is characterized by the occurrence of antinuclear antibodies (ANA) and antibodies to smooth muscles (ASMA). ASMA recognize antigenic structure formed by polymeric f-Actin. For type II a high prevalence of antibodies to liver and kidney microsomal antigens (LKM) has been described. LKM1 antibodies recognize epitopes of cytochrome P450 IID6, a 50 kDa cytoplasmic protein found in hepatocytes and proximal tubular kidney cells. LC1 antibodies are specific for type II hepatitis, too. The respective antigen is formiminotransferase / cyclodeaminase located in the cytosol of liver cells. Patients with type III autoimmune hepatitis exhibit antibodies to the soluble liver antigen (SLA). This type III is not yet fully accepted as independent subgroup for autoimmune hepatitis by the "International Hepatitis Group".

LKM-1 associated AIH predominantly occurs in girls between 2 and 14 years of age underlining the importance of this parameter for pediatrics.

3 TEST PRINCIPLE

The ELISA (Enzyme Linked Immunosorbent Assay) is an immunoassay for the determination of specific antibodies. The strips of the microtiter plate are coated with test-specific antigens. If antibodies are present in the patient's sample, they bind to the antigens. A secondary antibody conjugated with the enzyme peroxidase detects the generated immune complex. A colorless substrate is converted into the colored product. The signal intensity of the reaction product is proportional to the antibody activity in the sample. After stopping the signal intensity of the reaction product is measured photometrically.

4 TEST COMPONENTS

Component	Description
Microtiter plate A Ag 96 1 piece	12 breakable microtiter strips (ready-to-use), 8 wells per strip, each well coated with recombinant human cytochrome P450 IID6
Calibrator 0 – 4 CAL 5 x 1 mL, white cap	Colored dilutions of human serum (ready-to-use; contains ProClin 950) <i>The antibody activities are indicated on the quality control certificate.</i>
Negative control N CONTROL - 1 x 1 mL, green cap	Colored dilution of human serum (ready-to-use; contains ProClin 950) <i>The antibody activity is indicated on the quality control certificate.</i>
Positive control P CONTROL + 1 x 1 mL, red cap	Colored dilution of human serum (ready-to-use; contains ProClin 950) <i>The antibody activity is indicated on the quality control certificate.</i>
Sample diluent C DIL 5x 1 x 20 mL, black cap	Colored solution (5x; contains ProClin 950)
Wash buffer B BUF WASH 10x 1 x 100 mL, white cap	Concentrated solution (10x; contains ProClin 950)
Conjugate IgG D CONJ 1 x 15 mL, red cap	Colored solution of polyclonal anti-human IgG antibody conjugated to horseradish peroxidase (ready-to-use; contains ProClin 950)
Substrate E SOLN TMB 1 x 15 mL, blue cap	3,3',5,5'-Tetramethylbenzidin (ready-to-use)
Stop solution F H2SO4 0.25M 1 x 15 mL, yellow cap	0.25 M Sulfuric acid (ready-to-use)
Adhesive Foil 2 pieces	-
QC Certificate 1 piece	-
Instructions for Use 1 piece	-

5 MATERIALS REQUIRED BUT NOT PROVIDED

- Common laboratory equipment
- Precision pipettes (5 – 1000 µL), multi-channel pipettes (100 – 1000 µL) and disposable pipette tips
- Graduated cylinders (100 – 1000 mL)
- Sample tubes for the preparation of dilutions
- Vortex mixer or other rotators
- Microtiter plate washer or wash comb
- Microtiter plate reader with optical filters for 450 nm and 620 nm or 690 nm
- Adsorbent paper or paper towel
- Distilled or de-ionized water

6 STORAGE AND STABILITY

Upon receipt, all test components must be stored at 2 °C to 8 °C, preferably in the original kit box. If stored properly in their original containers, all components are stable until their expiry date.

All components are stable for at least 2 months after opening when stored properly at 2 °C to 8 °C.

7 GENERAL INFORMATION

This product is for *in vitro* diagnostic use only. The instructions for use must be carefully read before use. They are valid only for the present product with the given composition and must be strictly followed to ensure reliable test results. Deviations can lead to erroneous test results. Components must not be exchanged by test reagents of different lots or of other manufacturers.

Contamination of reagents must be avoided by use of aseptic techniques when removing aliquots from the vials. After use, reagent vials must be tightly closed with their corresponding caps.

Cross-contamination of samples or reagents can lead to inconsistent test results and must be avoided by use of consistent pipetting techniques.

Exposure of reagents to strong light must be avoided throughout the entire test procedure and storage.

Insufficient washing will result in poor precision and elevated measurement signals. After each washing step any residual fluid has to be removed completely.

8 PREPARATION

8.1 Preparation of Reagents

All components including the microtiter plate must be brought to room temperature (RT: 18 °C to 25 °C) before use for at least 30 minutes. All liquid components must be mixed gently to ensure homogeneity.

8.1.1 Microtiter Plate

The microtiter plate is sealed in an aluminium bag. Unused test strips should always be stored refrigerated and protected from moisture with the desiccant in the properly sealed aluminum bag.

Carefully resealed, the test strips can be used for 8 weeks after opening.

8.1.2 Calibrators

The calibrators are ready-to-use and must not be diluted any further. Calibrators must be used in each test run.

8.1.3 Controls

The positive and the negative controls are ready-to-use and must not be diluted any further. Controls must be used in each test run. Laboratories can also validate their own control samples and use them alternatively.

8.1.4 Sample Diluent

The sample diluent is concentrated and must be diluted 1:5 with distilled water before use (e. g. 20 mL + 80 mL). A sufficient amount of sample diluent solution must be prepared.

The diluted sample diluent solution can be stored at 2 °C to 8 °C up to 30 days.

8.1.5 Wash Buffer

The wash buffer is concentrated and must be diluted 1:10 with distilled water before use (e. g. 100 mL + 900 mL). A sufficient amount of washing solution must be prepared.

The diluted washing solution can be stored at 2 °C to 8 °C up to 30 days.

8.1.6 Conjugate

The conjugate is ready-to-use and stable up to 8 weeks after opening when stored at 2 °C to 8 °C.

8.1.7 Substrate

The substrate is ready-to-use. Exposure of the substrate solution to strong light should be avoided.

8.1.8 Stop Solution

The stop solution is ready-to-use.

8.2 Preparation of Samples

8.2.1 Sample Material

The use of freshly collected **serum** from blood taken by venipuncture is recommended.

The use of icteric, lipemic, hemolytic or bacterially contaminated samples should be avoided. Insoluble substances must be removed from the sample by centrifugation. Samples must not be thermally inactivated.

8.2.2 Sample Dilution

The samples must be **diluted 1:101** (e. g. 10 µL + 1000 µL) with prepared sample diluent solution and mixed thoroughly. Building of foam should be avoided.

8.2.3 Sample Storage

Samples may be kept at 2 °C to 8 °C up to three days.

Long-term storage requires -20 °C. Repeated freezing and thawing should be avoided. For multiple use, samples should be aliquoted and kept at -20 °C.

9 TEST PERFORMANCE

9.1 Pipetting Scheme

The following pipetting scheme is recommended:

	1	2	3	4
A	CAL 0	Sample 2		
B	CAL 1	Sample 3		
C	CAL 2	Sample 4		
D	CAL 3	Sample 5		
E	CAL 4	...		
F	N	...		
G	P	...		
H	Sample 1	...		

9.2 Procedure

The indicated incubation times and temperatures must be adhered to and significant time shifts during pipetting samples and reagents must be avoided. The microtiter plate should be shortly shaken after addition of reagents.

Step	Description
1. Addition of calibrators, controls and diluted samples	Add 100 µL ready-to-use calibrators, controls and diluted samples per well
2. Incubation	Cover the plate and incubate for 60 minutes at RT
3. Wash cycle	Aspirate the solution and wash 3 times with 300 µL washing solution with at least 5 seconds soaking time each; dry by tapping the microtiter plate on a paper towel to remove any residual droplets
4. Addition of conjugate	Add 100 µL ready-to-use conjugate to each well
5. Incubation	Cover the plate and incubate for 30 minutes at RT
6. Wash cycle	Aspirate the solution and wash 3 times with 300 µL washing solution with at least 5 seconds soaking time each; dry by tapping the microtiter plate on a paper towel to remove any residual droplets
7. Addition of substrate	Add 100 µL ready-to-use substrate to each well
8. Incubation	Cover the plate and incubate for 15 minutes in the dark at RT
9. Addition of Stop Solution	Add 100 µL ready-to-use stop solution to each well
10. Analysis	Read optical density (OD) at 450 nm versus 620 or 690 nm within 30 minutes after stopping the reaction

9.3 Automation

Automated processing of the immunoassays must be performed analogous to manual use and validated by the user.

10 TEST EVALUATION

10.1 Metrological Traceability

The immunoassay is calibrated using an internal reference sample. Quantitative results are expressed in U/mL.

10.2 Standard Curve

For generation of a standard curve, the optical signals (optical density, OD) of the calibrators are plotted against their antibody activities and correlated by a 4-parameter logistic (4 PL) fit. Antibody activities of unknown samples can be derived directly from their optical signals by use of the generated standard curve.

10.3 Criteria of Validity

Test runs are only valid if the following criteria of validity are fulfilled:

- OD CAL 0 < CAL 1 < CAL 2 < CAL 3 < CAL 4
- OD CAL 4 > 1.2
- The negative control must be evaluated negative.
- The positive control must be evaluated positive and present an antibody activity within the validity range indicated on the quality control certificate.

If these criteria are not met, the test is not valid and must be repeated.

10.4 Troubleshooting

In case of an invalid test run, the expiry dates and storage conditions, incubation times and temperatures, and precise calibration of all instruments used should be verified. If no reason for an invalid test run could be identified, please contact the supplier or manufacturer of the product.

10.5 Reference Ranges

The reference ranges are indicated below:

	Interpretation
Antibody activity < 15 U/mL	negative
Antibody activity ≥ 15 U/mL	positive

As a result of different seroprevalences in individual regions, each laboratory should verify the reference ranges by own analysis and adapt, if necessary.

10.6 Interpretation of Test Results

A **positive** test result indicates the presence of specific antibodies.

A **negative** result indicates the absence of specific antibodies, but does not exclude the possibility of an autoimmune reaction.

In case of a **borderline** test result, a reliable evaluation is not possible.

10.7 Limitations of the Method

The interpretation of test results must always be considered in combination with the clinical picture of the patient. The diagnosis should not be based on the results of a sole diagnostic method. All clinical and laboratory findings should be evaluated to state a diagnosis. For confirmation, further investigations should be carried out.

11 PERFORMANCE CHARACTERISTICS

11.1 Analytical Performance Characteristics

11.1.1 Precision

The precision of test results was assessed by the determination of the intra- and interassay variation by the analysis of multiple samples with different antibody activities.

	Intra assay Precision		Inter assay Precision	
	U/mL	CV (%)	U/mL	CV (%)
Sample 1	182.5	3.5	190.1	3.5
Sample 2	61.8	3.5	63.8	2.3
Sample 3	19.6	3.1	19.7	1.2

11.1.2 Measurement Range

Reliable accuracy, trueness, precision, linearity and recovery of test results have been observed within the measurement range of the assay from the LoQ to the upper calibrator in comprehensive studies.

Samples with test results above the upper calibrator should be reported as >max.

Samples with test results below the LoQ should be reported as <min.

If test results above the upper calibrator are observed, the samples may be tested at a higher dilution. The resulting antibody activity must be multiplied with the additional dilution factor.

11.2 Diagnostic Performance Characteristics

11.2.1 Diagnostic Sensitivity and Specificity

Sensitivity and specificity of the immunoassay were assessed by the analysis of 135 samples from patients with AIH Type II and unselected blood donors.

	Diagnostic Performance
Sensitivity	> 99 %
Specificity	> 96 %

12 WARNINGS AND PRECAUTIONS

The product is designed exclusively for *in vitro* diagnostic use by qualified, authorized and trained personnel. All test components and human samples should be handled with care as potentially hazardous. Good laboratory practices (GLP) and all relevant regulations should be adhered to.

In case the product is damaged or product information including labelling is wrong or incorrect, please contact the manufacturer or supplier.

This product contains preparations of human and / or animal origin. Any material derived from human body fluids or organs used for the preparation of components were tested and found negative for HBsAg (Hepatitis B-Virus-surface Antigen) and anti-HIV as well as anti-HCV antibodies. However, all components and all patient samples should be handled as potentially hazardous in accordance with national laws and appropriate guidelines on biological safety.

As the product contains potentially hazardous materials, the following precautions should be followed: Do not smoke, eat or drink while handling kit material or samples. Avoid direct contact to kit material or samples by wearing protective gloves laboratory coat and safety glasses. Never pipette material by mouth. Wipe up spills promptly and wash the affected surface thoroughly with a decontaminant. Wash hands thoroughly after use.

Some of the reagents contain ProClin (< 1.0 %) as a preservative, may cause skin sensitization (H317) and must not be swallowed or allowed to come into contact with skin or mucosa (P280, P333+P313).

The information in the safety data sheet on possible hazards, first aid measures, measures in the event of the unintentional release of large quantities, handling and storage, personal protective equipment, information on disposal as well as information on toxicology must be observed.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state in which the user and/or the patient is established.

13 DISPOSAL

For decontamination and disposal the recommendations of the CDC as well as the relevant local and national environmental guidelines and regulations should be adhered to. Samples, potentially contaminated materials and infectious waste must be decontaminated, e.g. by autoclaving for 20 minutes at 121 °C.

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