## Instructions for Use

# Noradrenaline ELISA

IVD

CE

REF EIA-4317

96

DRG Instruments GmbH, Germany  
Frauenbergstraße 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

DRG International, Inc., USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.**  
**Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.**  
**Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.**  
**Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.**  
**Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

<b>Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées</b>	
The following changes have been made in comparison to the previous version: Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen: Rispetto alla versione precedente, sono state apportate le seguenti modifiche: Se han introducido los siguientes cambios en comparación con la versión anterior: Les modifications suivantes ont été apportées par rapport à la version précédente :	
1.1 Intended Use and principle of the test	
2.1 Procedural cautions, guidelines, warning and limitations	
4.1 Contents of the kit	
Throughout the document	

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## 1 INTRODUCTION

### 1.1 Intended use and principle of the test

Enzyme Immunoassay for the quantitative determination of noradrenaline (norepinephrine) in plasma and urine.

Noradrenaline (norepinephrine) is extracted by using a cis-diol-specific affinity gel, acylated and then converted enzymatically.

The subsequent competitive ELISA uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The acylated standards, controls and samples compete with the solid phase bound analytes compete for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate resulting in a colour reaction. The reaction is monitored at a wavelength of 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a reference curve prepared with known standard concentrations. Manual processing of the ELISA is recommended. The use of automatic laboratory equipment is the responsibility of the user. This in-vitro diagnostic is for professional use only.

### 1.2 Clinical application

In humans the catecholamines adrenaline (epinephrine), noradrenaline (norepinephrine) and dopamine are neurotransmitters of the sympathetic nervous system and are involved in many physiological processes. The sympathetic nervous system sets the body to a heightened state of alert, also called as the body's fight-or-flight response.

In the human body the catecholamines and their metabolites indicate the adaptation of the body to acute and chronic stress.

Next to the metanephrine/normetanephrine the catecholamines are important for the diagnosis and the follow-up of tumors of the sympathoadrenal system like the pheochromocytoma. The quantitative determination of catecholamines in urine is preferred for the diagnosis of these tumors, whereas the determination of catecholamines in plasma is medically sensible for the localization of the tumor and for function testing. Values above the cut-off can provide an indication for neuroendocrine tumors.

However, in literature various diseases like hypertension, cardiovascular diseases, schizophrenia and manic depression are described with abnormal low or high levels of catecholamines

Therapeutic consequences should never be based on laboratory results alone even if these results are assessed in accordance quality criteria of the method. Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient it can be used for therapeutic consequences.

## 2 PROCEDURAL CAUTIONS, GUIDELINES, WARNINGS AND LIMITATIONS

### 2.1 Procedural cautions, guidelines and warnings

1. This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and must be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
2. This assay was validated for a certain type of sample as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
3. The principles of Good Laboratory Practice (GLP) must be followed.
4. In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
5. If serious incidents should occur in connection with this product, they should be reported to the manufacturer and the competent national authorities.
6. All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water. Avoid repeated freezing and thawing of reagents and specimens.
7. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided. Microtiter strips which are removed from the frame for usage should be marked accordingly to avoid any mix-up.
8. Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
9. Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
10. Incubation times do influence the results. All wells should be handled in the same order and time intervals.
11. To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
12. A standard curve must be established for each run.
13. The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.

14. Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
15. Avoid contact with Stop Solution containing 0.25 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
16. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Rinse contaminated items before reuse.
17. For information on hazardous substances included in the kit please refer to Safety Data Sheet (SDS). The Safety Data Sheet for this product is made available upon request.
18. Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.
19. The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
20. In case of any severe damage to the test kit or components, the manufacturer has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components must not be used for a test run. They must be stored properly until the manufacturer decides what to do with them. If it is decided that they are no longer suitable for measurements, they must be disposed of in accordance with national regulations.
21. The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence but must be correlated to other diagnostic tests and clinical observations.

## 2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

### 2.2.1 Interfering substances and proper handling of specimens

#### Plasma

Samples containing precipitates or fibrin strands or which are hemolytic or lipemic might cause inaccurate results. Hemolytic samples (up to 4 mg/ml hemoglobin), icteric samples (up to 50 mg/dl bilirubin) and lipemic samples (up to 800 mg/dl triglycerides) have no influence on the assay results.

If the concentrations cannot be estimated and there are doubts as to whether the above limit values for hemolytic, icteric or lipemic samples are complied with, the samples should not be used in the assay.

#### 24-hour urine

Please note the sample collection! If the percentage of the final concentration of acid is too high, this will lead to incorrect results for the urine samples.

### 2.2.2 Drug and food interferences

There are no known substances (drugs) which ingestion interferes with the measurement of noradrenaline level in the sample.

### 2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

## 3 STORAGE AND STABILITY

Store kit and reagents at 2 °C - 8 °C until expiration date.

Do not use kit and components beyond the expiry date indicated on the kit labels.

Once opened the reagents are stable for 2 months when stored at 2 °C - 8 °C.

Once the resealable pouch of the ELISA plate has been opened, care should be taken to close it tightly again including the desiccant.

## 4 MATERIALS

### 4.1 Contents of the kit

	<b>FOILS</b>	<b>Adhesive Foil</b> - ready to use
Content:	Adhesive foils in a resealable pouch	
Number:	1 x 4 foils	
	<b>WASH-CONC 50x</b>	<b>Wash Buffer Concentrate</b> - Concentrated 50x
Content:	Buffer with a non-ionic detergent and physiological pH	
Volume:	1 x 20 mL/vial, purple cap	
	<b>CONJUGATE</b>	<b>Enzyme Conjugate</b> - ready to use
Content:	Goat anti-rabbit immunoglobulins conjugated with peroxidase	
Volume:	1 x 12 mL/vial, red cap	
Description:	Species is goat	
	<b>SUBSTRATE</b>	<b>Substrate</b> - ready to use
Content:	Chromogenic substrate containing 3,3',5,5'-tetramethylbenzidine, substrate buffer and hydrogen peroxide	
Volume:	1 x 12 mL/vial, black cap	
	<b>STOP-SOLN</b>	<b>Stop Solution</b> - ready to use
Content:	0.25 M sulfuric acid	
Volume:	1 x 12 mL/vial, light grey cap	
	<b>■■■ NAD NMN</b>	<b>Noradrenaline Microtiter Strips</b> - ready to use
Content:	1 x 96 well (12x8) antigen precoated microwell plate in a resealable yellow pouch with desiccant	
	<b>NAD-AS</b>	<b>Noradrenaline Antiserum</b> - ready to use
Content:	Rabbit anti-noradrenaline antibody in buffer with proteins and non-mercury preservative, yellow coloured	
Volume:	1 x 6 mL/vial, yellow cap	
Description:	Species of antibody is rabbit, species of protein in buffer is bovine	
	<b>ACYL-REAG</b>	<b>Acylation Reagent</b> - ready to use
Content:	Acylation reagent in DMSO	
Volume:	1 x 3 mL/vial, white cap	
	<b>ADJUST-BUFF</b>	<b>Adjustment Buffer</b> - ready to use
Content:	TRIS buffer	
Volume:	1 x 4 mL/vial, green cap	
	<b>ACYL-BUFF</b>	<b>Acylation Buffer</b> - ready to use
Content:	Buffer with light alkaline pH for the acylation	
Volume:	1 x 20 mL/vial, white cap	

**ASSAY-BUFF****Assay Buffer - ready to use**

Content: 1M hydrochloric acid and a non-mercury preservative  
 Volume: 1 x 6 mL/vial, grey cap  
 Hazard pictograms: 

**GHS05**

Signal word: Danger  
 Hazard statements: H314 Causes severe skin burns and eye damage.  
 Precautionary statements: P280 Wear protective gloves, protective clothing, eye protection.  
 P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.  
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P310 Immediately call a doctor, a POISON CENTER.  
 P501 Dispose of contents/container to an authorised waste collection point.

**COENZYME****Coenzyme - ready to use**

Content: S-adenosyl-L-methionine  
 Volume: 1 x 4 mL/vial, purple cap

**ENZYME****Enzyme - lyophilized**

Content: Catechol-O-methyltransferase  
 Volume: 2 vials, pink cap  
 Description: Catechol-O-methyltransferase from pig liver

**EXTRACT-BUFF****Extraction Buffer - ready to use**

Content: Buffer containing carbonate  
 Volume: 1 x 6 mL/vial, brown cap

**EXTRACT-PLATE 48****Extraction Plate - ready to use**

Content: 2 x 48 well plates coated with boronate affinity gel in a resealable pouch

**HCL****Hydrochloric Acid - ready to use**

Content: 0.025 M Hydrochloric Acid, yellow coloured  
 Volume: 1 x 20 mL/vial, green cap

**Standards and Controls - Ready to use**

Component	Colour/Cap	Concentration ng/mL <b>NAD</b>	Concentration nmol/L <b>NAD</b>	Volume/ Vial
<b>STANDARD A</b>	white	0	0	4 mL
<b>STANDARD B</b>	yellow	5	30	4 mL
<b>STANDARD C</b>	orange	20	118	4 mL
<b>STANDARD D</b>	blue	75	443	4 mL
<b>STANDARD E</b>	grey	250	1478	4 mL
<b>STANDARD F</b>	black	1000	5910	4 mL
<b>CONTROL 1</b>	green	Refer to QC report for expected value and acceptable range!		4 mL
<b>CONTROL 2</b>	red			4 mL

Conversion: Noradrenaline (ng/mL) x 5.91 = Noradrenaline (nmol/L)

Content: Acidic buffer with non-mercury stabilizer, spiked with defined quantity of noradrenaline

#### 4.2 Additional materials required but not provided in the kit

- Water (deionized, distilled, or ultra-pure)
- Absorbent material (paper towel)

#### 4.3 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 10 – 700 µL; 1 mL
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 - 650 nm
- Microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Vortex mixer

### 5 SAMPLE COLLECTION, HANDLING AND STORAGE

#### Plasma

Whole blood should be collected into centrifuge tubes containing EDTA as anti-coagulant and centrifuged according to manufacturer's instructions immediately after collection.

In case of hemolytic, icteric or lipemic samples see 2.2.1.

Storage: up to 6 hours at 2 °C - 8 °C,  
for longer period (up to 6 months) at -20 °C.

Repeated freezing and thawing should be avoided.

#### Urine

Spontaneous urine or 24-hour urine, collected in a bottle containing 10 - 15 mL of 6 M HCl, can be used.

If 24-hour urine is used please record the total volume of the collected urine.

Storage: up to 48 hours at 2 °C - 8 °C,  
up to 24 hours at room temperature,  
for longer periods (up to 6 months) at -20 °C. Repeated freezing and thawing should be avoided.

Avoid exposure to direct sunlight.

### 6 TEST PROCEDURE

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Number the Extraction Plate and microwell plates (Microtiter Strips which are removed from the frame for usage should be marked accordingly to avoid any mix-up). Duplicate determinations are recommended.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the enzyme immunoassay is between 20 °C – 25 °C.

**⚠** *The use of a microtiter plate shaker with the following specifications is mandatory: shaking amplitude 3 mm; approx. 600 rpm. Shaking with differing settings might influence the results.*

#### 6.1 Preparation of reagents

##### Wash Buffer

Dilute the 20 mL **WASH-CONC 50X** with water to a final volume of 1000 mL.

Storage: 2 months at 2 °C - 8 °C

##### Enzyme Solution

Reconstitute the content of the vial labelled **ENZYME** with 1 mL water (deionized, distilled, or ultra-pure) and mix thoroughly. Add 0.3 mL of **COENZYME** followed by 0.7 mL of **ADJUST-BUFF**. The total volume of the Enzyme Solution is 2.0 mL.

**⚠** *The Enzyme Solution has to be prepared freshly prior to the assay (not longer than 10 - 15 minutes in advance). Discard after use!*

## Noradrenaline Microtiter Strips

In rare cases residues of the blocking and stabilizing reagent can be seen in the wells as small, white dots or lines. These residues do not influence the quality of the product.

## Acylation Reagent

The **ACYL-REAG** has a freezing point of 18.5 °C. To ensure that it is liquid when being used, it must be ensured that the Acylation Reagent has reached room temperature and forms a homogeneous, crystal-free solution before being used.

### 6.2 Sample preparation, extraction and acylation

1. Pipette **10 µL** of **standards, controls, urine samples** and **300 µL** of **plasma samples** into the respective wells of the **EXTRACT-PLATE 48**.
2. Add **250 µL** of **water** (deionized, distilled, or ultra-pure) to the wells with **standards, controls** and **urine samples**.
3. Pipette **50 µL** of **ASSAY-BUFF** into all wells.
4. Pipette **50 µL** of **EXTRACT-BUFF** into all wells.
5. Cover plate with **FOILS** and incubate **30 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
6. Remove the foil. Empty plate and blot dry by tapping the inverted plate on absorbent material.
7. Pipette **1 mL** of **Wash Buffer** into all wells. Incubate the plate for **5 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm). Empty plate and blot dry by tapping the inverted plate on absorbent material.
8. Pipette another **1 mL** of **Wash Buffer** into all wells. Incubate the plate for **5 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm). Empty plate and blot dry by tapping the inverted plate on absorbent material.
9. Pipette **150 µL** of **ACYL-BUFF** into all wells.
10. Pipette **25 µL** of **ACYL-REAG** into all wells.
11. Incubate **15 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
12. Empty plate and blot dry by tapping the inverted plate on absorbent material.
13. Pipette **1 mL** of **Wash Buffer** into all wells. Incubate the plate for **10 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm). Empty plate and blot dry by tapping the inverted plate on absorbent material.
14. Pipette **150 µL** of **HCL** into all wells.
15. Cover plate with **FOILS**. Incubate **10 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm). Remove the foil and discard.

**Do not decant the supernatant thereafter!**

The following volumes of the supernatant are needed for the subsequent ELISA:

<b>Noradrenaline</b>	<b>20 µL</b>
----------------------	--------------

### 6.3 Noradrenaline ELISA

1. Pipette **25 µL** of the **Enzyme Solution** (refer to 6.1) into all wells of the **W NAD NMN**.
2. Pipette **20 µL** of the extracted **standards, controls and samples** into the appropriate wells.
3. Incubate for **30 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
4. Pipette **50 µL** of the **NAD-AS** into all wells and cover plate with **FOILS**.
5. Incubate for **2 h at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
6. Remove the foil. Discard or aspirate the content of the wells. Wash the plate **3 x** by adding **300 µL** of **Wash Buffer**, **discarding** the content and **blotting dry each time** by tapping the inverted plate on absorbent material.
7. Pipette **100 µL** of the **CONJUGATE** into all wells.
8. Incubate for **30 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm).
9. Discard or aspirate the content of the wells. Wash the plate **3 x** by adding **300 µL** of **Wash Buffer**, **discarding** the content and **blotting dry each time** by tapping the inverted plate on absorbent material.
10. Pipette **100 µL** of the **SUBSTRATE** into all wells and incubate for **25 ± 5 min at RT** (20 °C - 25 °C) on a **shaker** (approx. 600 rpm). **Avoid exposure to direct sunlight!**
11. Add **100 µL** of the **STOP-SOLN** to all wells and shake the microtiter plate shortly.
12. **Read** the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to **450 nm** (if available a reference wavelength between 620 nm and 650 nm is recommended).

## 7 CALCULATION OF RESULTS

Measuring range		Noradrenaline
	Urine	2.5 - 1000 ng/mL
	Plasma	93 – 33333 pg/mL

The standard curve, which can be used to determine the concentration of the unknown samples, is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis) using a concentration of 0.001 ng/ml for Standard A (this alignment is mandatory because of the logarithmic presentation of the data). Use non-linear regression for curve fitting (e. g. 4-parameter, marquardt).

 *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

### Urine samples and controls

The concentrations of the **urine samples** and the **Controls 1 and 2** can be read directly from the standard curve.

Calculate the 24 h excretion for each urine sample:  $\mu\text{g}/24\text{ h} = \mu\text{g/L} \times \text{L}/24\text{ h}$

### Plasma samples

The read concentrations of the **plasma samples** have to be **divided by 30**.

### Conversion

Noradrenaline [ng/mL]  $\times$  5.91 = Noradrenaline [nmol/L]

#### 7.1 Expected reference value


It is strongly recommended that each laboratory should determine its own normal and abnormal values.

	Noradrenaline
24-hour urine	< 90 $\mu\text{g}/\text{day}$ (535 nmol/day)
Plasma	< 600 pg/mL

#### 7.2 Typical standard curve

*Example, do not use for calculation!*

## 8 CONTROL SAMPLES

It is recommended to use control samples according to national regulations. Use controls at both normal and pathological levels. Commercially obtained control samples should be treated like unknown samples. Control samples should fall within established confidence limits. The confidence limits of the kit controls are indicated on the QC-Report.

## 9 ASSAY CHARACTERISTICS

<b>Analytical Sensitivity</b>		
<b>Noradrenaline</b>		
Limit of Blank (LOB)	Urine [ng/ml]	1.5
	Plasma [pg/ml]	32
Limit of Detection (LOD)	Urine [ng/ml]	1.7
	Plasma [pg/ml]	36
Limit of Quantification (LOQ)	Urine [ng/ml]	2.5
	Plasma [pg/ml]	93

<b>Analytical Specificity (Cross Reactivity)</b>	
<b>Substance</b>	<b>Cross Reactivity [%]</b>
	Noradrenaline
Derivatized Adrenaline	0.08
Derivatized Noradrenaline	100
Derivatized Dopamine	0.03
Metanephrine	< 0.01
Normetanephrine	0.16
3-Methoxytyramine	< 0.01
3-Methoxy-4-hydroxyphenylglycol	< 0.01
Tyramine	< 0.01
Phenylalanine, Caffeinic acid, L-Dopa, Homovanillic acid, Tyrosine, 3-Methoxy-4-hydroxymandelic acid	< 0.01

<b>Precision</b>							
<b>Intra-Assay Urine (n = 60)</b>				<b>Intra-Assay Plasma (n = 60)</b>			
	<b>Sample</b>	<b>Range (ng/mL)</b>	<b>CV (%)</b>		<b>Sample</b>	<b>Range (pg/mL)</b>	<b>CV (%)</b>
Noradrenaline	1	26.1 ± 3.6	13.8	Noradrenaline	1	510 ± 65	12.8
	2	97 ± 12.8	13.4		2	1358 ± 194	14.3
	3	267 ± 35	13.1		3	3363 ± 374	11.1
<b>Inter-Assay Urine (n = 33)</b>				<b>Inter-Assay Plasma (n = 18)</b>			
Noradrenaline	<b>Sample</b>	<b>Range (ng/mL)</b>	<b>CV (%)</b>	Noradrenaline	<b>Sample</b>	<b>Range (pg/mL)</b>	<b>CV (%)</b>
	1	19.5 ± 3.9	20.0		1	445 ± 40.9	9.2
	2	80.6 ± 10.6	13.2		2	1232 ± 134	10.9
	3	226 ± 39.5	17.4		3	3283 ± 302	9.2

<b>Linearity</b>			<b>Serial dilution up to</b>	<b>Range (%)</b>	<b>Mean (%)</b>
	Noradrenaline		Urine	1:512	112
		Plasma	1:512	102 - 125	112

<b>Recovery</b>			<b>Mean (%)</b>	<b>Range (%)</b>	<b>Range</b>
	Noradrenaline		Urine	91 - 113	58.6 – 260 ng/mL
		Plasma	87	75 - 107	51 – 14,251 pg/mL

## 10 REFERENCES/LITERATURE

1. Kim et al. Vitamin C prevents stress-induced damage on the heart caused by the death of cardiomyocytes, through the down-regulation of the excessive production of catecholamine, TNF- $\alpha$ , and ROS production in GULO(-I-) Vit C-Insufficient mice. Free Radical Biology and Medicine, 65:573-583 (2013)
2. Bada et al. Peripheral vasodilatation determines cardiac output in exercising humans: insight from atrial pacing. The Journal of Physiology, 590(8):2051-2060 (2012)
3. Parks et al. Employment and work schedule are related to telomere length in women. Occupational & Environmental Medicine 68(8):582-589 (2011)

*For updated literature or any other information please contact your local supplier.*

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## 10 REFERENZEN/LITERATUR

- (1) Kim et al. Vitamin C prevents stress-induced damage on the heart caused by the death of cardiomyocytes, through the down-regulation of the excessive production of catecholamine, TNF- $\alpha$ , and ROS production in GULO(-I-) Vit C-Insufficient mice. Free Radical Biology and Medicine, 65:573-583 (2013)
- (2) Bada et al. Peripheral vasodilatation determines cardiac output in exercising humans: insight from atrial pacing. The Journal of Physiology, 590(8):2051-2060 (2012)
- (3) Parks et al. Employment and work schedule are related to telomere length in women. Occupational & Environmental Medicine 68(8):582-589 (2011)

**SYMBOLS USED**

<b>Symbol</b>	<b>English</b>	<b>Deutsch</b>	<b>Français</b>	<b>Español</b>	<b>Italiano</b>
	European Conformity	CE-Konformitäts-kennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Para uso Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di Cat.
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeits-datum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore
<i>Content</i>	Content	Inhalt	Contenu	Contenido	Contenuto
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Numéro	Volumen/Número	Volume/Quantità
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Plaques de micro-titration	Placas multipocillo	Micropozzetti
<i>Antiserum</i>	Antiserum	Antiserum	Antisérum	Antisuero	Antisiero
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Conjugué enzymatique	Conjugado enzimático	Tracciante enzimatico
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complexe enzymatique	Complex enzimático	Complesso enzimatico
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Solution substrat	Solución de sustrato	Soluzione di substrato
<i>Stop Solution</i>	Stop Solution	Stoplösung	Solution d'arrêt	Solución de parada	Soluzione d'arresto
<i>Zero Standard</i>	Zero Standard	Nullstandard	Zero Standard	Estándar cero	Standard zero
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Contrôle	Control	Controllo
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampon d'essai	Tampón de ensayo	Tampone del test
<i>Wash Solution</i>	Wash Solution	Waschlösung	Solution de lavage	Solución de lavado	Soluzione di lavaggio
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)
<i>1 N HCl</i>	1 N HCl	1 N HCl	1N HCl	1 N HCl	
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungs-medium	Solution pour dilution de l'échantillon	Solución para dilución de la muestra	Diluente dei campioni
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungs-medium	Solution pour dilution du conjugué	Solución para dilución del conjugado	Diluente del tracciante