## Instructions for Use

# ImAnOx (TAS/TAC) Kit

IVD

REF ENZ-5814

$\Sigma$  96

**DRG**

**DRG Instruments GmbH**, Germany  
Frauenbergstraße 18, 35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by:

**DRG**

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.  
 Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
 Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.  
 Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.  
 Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

<b>Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées</b>	
The following changes have been made in comparison to the previous version: Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen: Rispetto alla versione precedente, sono state apportate le seguenti modifiche: Se han introducido los siguientes cambios en comparación con la versión anterior: Les modifications suivantes ont été apportées par rapport à la version précédente :	
	Changed wording in several chapters.
3 MATERIAL SUPPLIED:	Volume for Reaction buffer A changed from 105 ml to 100 mL.
5 STORAGE AND PREPARATION OF REAGENTS - Calibrators and controls:	The exact concentration is stated in the respective product specifications and no longer on the label.

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## 1 INTENDED USE

The photometric assay is intended for the quantitative determination of the **total antioxidative status/capacity (TAS/TAC)** in EDTA-plasma and serum.

For *in vitro* diagnostic only.

## 2 INTRODUCTION

The human body is constantly under attack from free radicals that occur as part of normal cell metabolism, and by exposure to environmental factors such as UV light, cigarette smoke, environmental pollutants and gamma radiation. The resulting „Reactive Oxygen Species“ (ROS) circulates freely in the body with access to all organs and tissues, which can have serious repercussions throughout the body. The body possesses a number of mechanisms both to control the production of ROS and to cope with free radicals in order to limit or repair the damage to tissues

Overproductions of ROS or insufficient defense mechanisms lead to a dangerous disbalance in the organism. Thereby several pathomechanisms implicated in over 100 human diseases, e.g. cardiovascular disease, cancer, diabetes mellitus, inflammatory disease, aging, etc., were induced.

Determination of the antioxidative capacity becomes of fundamental importance in medical diagnosis and in research. The ImAnOx-Assay is fast, reliable and easy to perform. The total antioxidative capacity is measured.

## 3 MATERIAL SUPPLIED

Content	Kit Components	Quantity
RECSOL	Reconstitution solution	15 mL
PLATE	Microtiter plate	2 x
CAL	Calibrator, lyophilised (see specification data sheet for concentration)	4 x
CTRL1	Control 1; lyophilised	4 x
CTRL2	Control 2; lyophilised	4 x
PER	Peroxide solution	250 µL
REABUF A	Reaction buffer A	100 mL
REABUF B	Reaction buffer B	1.5 mL
ENZ	Enzyme solution	50 µL
STOP	Stop solution	15 mL

## 4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Vortex mixer
- Various pipettes
- Incubation chamber for 37 °C
- Microtiter plate reader

## 5 STORAGE AND PREPARATION OF REAGENTS

### **Calibrators and controls**

The **lyophilised calibrator** (CAL) and the **lyophilised controls 1 and 2** (CTRL1 and CTRL2) are stable at **-20 °C** until the expiry date stated on the label.

Before use, they have to be reconstituted each with **250 µL reconstitution solution** (RECSOL). Allow the vial content to dissolve for 5 min and then mix thoroughly by vortexing.

Aliquots of the **calibrator** (reconstituted CAL) and the **controls 1 and 2** (reconstituted CTRL1 and CTRL2) can be stored at **-20 °C** for **4 weeks**. **Avoid repeated thawing and freezing.**

The concentration of calibrator and controls slightly changes from lot to lot. The exact concentration is stated in the respective product specifications.

### **Storage of the other reagents**

Test reagents are stable until the expiry date (see label of test package) when stored at **2 °C - 8 °C**.

### **Preparation of Reagent 1**

**In order to ensure the functionality of reagent 1, it is important to use light-protected containers (e.g. brown glass) during preparation.**

Reagent 1 has to be prepared directly before use.

Peroxide solution (PER) is diluted in reaction buffer (REABUF).

For example:

5 mL reaction buffer A (REABUF A) + 10 µL peroxide solution (PER) = dilution 1

100 µL of dilution 1 + 4.9 mL reaction buffer A (REABUF A) = reagent 1

**Please note: Reagent 1 is not stable and cannot be stored.**

**Important:** The amounts are sufficient for 50 tests (25 duplicates). For varying sample numbers, the buffer volumes must be adjusted accordingly.

### **Preparation of Reagent 2**

**In order to ensure the functionality of reagent 2, it is important to use light-protected containers (e.g. brown glass) during preparation.**

During the incubation, hydrogen peroxide generates reaction products which absorb at 450 nm. Due to this effect and self-absorption, it is important to measure **with and without addition of enzyme**.

Reagent 2 has to be prepared directly before use.

Reaction buffer A (REABUF A) and reaction buffer B (REABUF B) are mixed with and without addition of enzyme solution (ENZ).

For example:

Reagent 2a: 5 mL reaction buffer A (REABUF A) + 100 µL reaction buffer B (REABUF B) + 5 µL enzyme solution (ENZ)  
= reagent 2a **with** enzyme solution

Reagent 2b: 5 mL reaction buffer A (REABUF A) + 100 µL reaction buffer B (REABUF B)  
= reagent 2b **without** enzyme solution

**Please note: Reagent 2a with and reagent 2b without addition of enzyme solution (ENZ) is not stable and cannot be stored.**

**Important:** To avoid losses, the enzyme solution (ENZ) should be centrifuged prior to use so as not to lose any liquid possibly adhering to the lid. After use, the vial should be immediately and correctly closed to avoid contamination or evaporation (e. g. parafilm).

**Important:** The amounts are sufficient for 50 tests (25 duplicates). For varying sample numbers, the buffer volumes must be adjusted accordingly.

## **6 STORAGE AND PREPARATION OF SAMPLES**

Serum and EDTA-plasma derived from venous fasting blood is suitable for this test system.

Lipemic and hemolytic samples affect the test results. Such samples should not be measured.

### **Sample storage**

The blood sample can be shipped as whole blood at 2 °C - 8 °C within 24 hours.

Serum and EDTA-plasma should be stored at -20 °C until the measurement takes place. They are stable at -20 °C for 4 weeks.

### **Preparation**

Samples with visible amounts of precipitates should be centrifuged (5 min at 10 000 g) prior to measurement, and the resulting supernatant is used in the test.

For testing in duplicates, pipette **4 x 10 µL** of each prepared sample per well.

## 7 ASSAY PROCEDURE

### 7.1 Principle of the test

The determination of the antioxidative capacity is performed by the reaction of antioxidants in the sample with a defined amount of exogenously provided hydrogen peroxide ( $H_2O_2$ ). The antioxidants in the sample eliminate a certain amount of the provided hydrogen peroxide. The residual  $H_2O_2$  is determined photometrically by an enzymatic reaction which involves the conversion of TMB to a colored product.

After addition of a stop solution, the samples are measured at 450 nm in a microtiter plate reader. The quantification is performed by the delivered calibrator.

The difference between applied and measured peroxide concentration in a defined time period is proportional to the reactivity of the antioxidants of the sample (antioxidative capacity). Quantification is performed with the enclosed calibrator.

**Please note:** As the reaction speed of the distinct antioxidants in the sample is different, the measured concentrations of antioxidants are equivalent to the reactivity of the distinct antioxidants and not to their total amount in the sample. Therefore, we use in our test system hydrogen-peroxide equivalents as unit for the antioxidative capacity.

The concentration of total antioxidative status/capacity (TAS/TAC) can be quantified by referring the optical density of the calibrator.

### 7.2 Test procedure

Bring all reagents and samples to room temperature and mix well.

Mark the positions of CAL /SAMPLE /CTRL (calibrator/sample/controls) on a protocol sheet.

#### Example

	with enzyme		without enzyme		with enzyme		without enzyme		with enzyme		without enzyme	
	1	2	3	4	5	6	7	8	9	10	11	12
A	P1	P1	P1	P1	P9	P9	P9	P9	P14	P14	P14	P14
B	P2	P2	P2	P2	Ctrl1	Ctrl1	Ctrl1	Ctrl1	P15	P15	P15	P15
C	P3	P3	P3	P3	P10	P10	P10	P10	P16	P16	P16	P16
D	P4	P4	P4	P4	Cal	Cal	Cal	Cal	P17	P17	P17	P17
E	P5	P5	P5	P5	P11	P11	P11	P11	P18	P18	P18	P18
F	P6	P6	P6	P6	Ctrl2	Ctrl2	Ctrl2	Ctrl2	P19	P19	P19	P19
G	P7	P7	P7	P7	P12	P12	P12	P12	P20	P20	P20	P20
H	P8	P8	P8	P8	P13	P13	P13	P13	P21	P21	P21	P21

Px = patient sample, Cal = calibrator, Ctrl = control

Take as many microtiter strips as needed from kit. Store unused strips covered at 2 °C - 8 °C. Strips are stable until expiry date stated on the label.

We recommend to carry out the tests in duplicate.

1. Add **10 µL of calibrator** (reconstituted CAL), **sample or control 1 or 2** (reconstituted CTRL1 or CTRL2) into the respective wells.  
**Please note:** Always pipet samples, controls and calibrator in 4 wells, 2 wells for enzyme treatment and 2 without enzyme treatment.
2. Add **100 µL** freshly prepared **reagent 1** in each well.
3. Incubate for **10 min** at **37 °C**
4. Add **100 µL** freshly prepared **reagent 2a** (with enzyme), respectively **reagent 2b** (without enzyme) in the corresponding wells.
5. Incubate for **5 min** at **room temperature**.
6. Add **50 µL** stop solution (STOP) and mix well.
7. Determine **absorption immediately** with a microtiter plate reader at **450 nm**.

## 8 RESULTS

### 8.1 Calculation

The difference of the sample values with and without enzyme is inversely proportional to the antioxidative capacity: To get the  $\Delta OD$ , subtract the OD-values of samples without enzyme from the OD-values of samples with enzyme. The antioxidative capacity is calculated according to the following formula:

$$\text{antioxidative capacity of a sample } [\mu\text{mol/L}] = 392 - (392 - \text{calibrator concentration}) \times \frac{\Delta OD_{\text{sample}}}{\Delta OD_{\text{calibrator}}}$$

## 9 LIMITATIONS

The use of heparin plasma results in clouding of the starting solution and thus to incorrect measurement results. Strong haemolytic and lipaemic samples often show pathological antioxidant capacity. We therefore do not recommend the measurement of such samples. Whole blood is not suitable for the measurement.

## 10 QUALITY CONTROL

DRG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### 10.1 Reference range

Based on studies of EDTA-plasma and serum of apparently healthy persons (n = 69), the following reference values have been estimated:

low antioxidative capacity	< 280 $\mu\text{mol/L}$
middle antioxidative capacity	280 - 320 $\mu\text{mol/L}$
high antioxidative capacity	> 320 $\mu\text{mol/L}$
Mean value	305 $\mu\text{mol/L}$

We recommend each laboratory to establish its own reference range. The values mentioned above are only for orientation and can deviate from other published data.

## 11 PERFORMANCE CHARACTERISTICS

### 11.1 Precision and reproducibility

#### Intra-Assay (n = 12)

Sample	ImAnOx [ $\mu\text{mol/L}$ ]	CV [%]
1	213	3.99
2	308	2.04

#### Inter-Assay (n = 12)

Sample	ImAnOx [ $\mu\text{mol/L}$ ]	CV [%]
1	217	2.65
2	285	3.89

### 11.2 Analytical Sensitivity

Limit of detection 130  $\mu\text{mol/L}$

## 12 PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- The stop solution consists of diluted sulphuric acid ( $\text{H}_2\text{SO}_4$ ), a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are harmful to health and the environment. Substrates for enzymatic color reactions can also cause skin and/or respiratory irritation. Any contact with the substances should be avoided. Further safety information can be found in the safety data sheet, which is available on request.
- As a precaution, it is recommended that the human material used is always considered potentially infectious.

## 13 GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which has not been consulted with the producer, may influence the results of the test. DRG can therefore not be held responsible for any damage resulting from incorrect use.
- Do not interchange different lot numbers of any kit component within the same assay.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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