
DRG:HYBRID•XL[®]

Instructions for Use

Aldosterone

IVD

CE

REF

HYE-5338

40 tests

DRG Instruments GmbH, Germany
Frauenbergstraße 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de


Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

ASSAY PROTOCOL BARCODE (APB)

**(Version 2.70 Software or later / ab Softwareversion 2.70 oder höher / Versione Software 2.70 o superiore /
 Versión de Software 2.70 o posterior / Version de logiciel 2.70 et supérieure /
 Wersja oprogramowania 2.70 lub późniejsza)**

436 - ALDO - HYE-5338 - 02.70.01

HYE-5338– v2.70.1


The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the SCAN NEW LOT page.
 Der Barcode muss in dem Menü „NEUES LOT SCANNEN“ eingelesen werden, um das Protokoll in der DRG:HYBRiD-XL-Software zu installieren.

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Please refer to section **3: Routine Procedures: "Installing a new assay product"** of the User Manual v2.60 or later.

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1 INTENDED USE

The DRG:HYBRID-XL Aldosterone is an automated enzyme immunoassay for the quantitative measurement of aldosterone in serum, plasma (EDTA, lithium heparin or citrate plasma) and urine.

For in vitro diagnostic use.

For laboratory professional use only with the DRG:HYBRID-XL Analyzer.

The device is **intended** to be used as an aid to diagnosis of primary and secondary aldosteronism.

The device is **not** intended for the diagnosis of adenomas.

1.1 Scientific Validity

The steroid hormone aldosterone is a potent mineral corticoid that is produced by the zona glomerulosa of the adrenal cortex in the adrenal gland. The synthesis and release are controlled by the renin-angiotensin-aldosterone system (RAAS)¹, as well as by plasma potassium concentration^{2,12}, the pituitary peptide ACTH, and by the blood pressure via pressure sensitive baroreceptors in the vessel walls of nearly all large arteries of the body^{3,12}. Aldosterone binds to mineralocorticoid receptors (MR) and triggers the transcription of hormone-responsive genes. In consequence, aldosterone increases the blood pressure by reabsorption of sodium and water from the distal tubules of the kidney into the blood, secretion of potassium into the urine, and elevation of circulating blood volume. Chronic overproduction and secretion of aldosterone leads to hypertension. Aldosterone activity is reduced in Addison's disease and increased in Conn's syndrome.

Primary hyperaldosteronism, which may be caused by aldosterone-secreting adrenal adenoma/carcinomas or adrenal cortical hyperplasia, is characterized by hypertension accompanied by increased aldosterone levels, hypernatremia, and hypokalemia. Secondary hyperaldosteronism (e.g. in response to renovascular disease, salt depletion, potassium loading, cardiac failure with ascites, pregnancy, Bartter's syndrome) is characterized by increased aldosterone levels and increased plasma renin activity.^{4,5,8,9,10,11,13}

Condition	Serum Aldosterone	Plasma Renin
Primary Aldosteronism	High	Low
Secondary Aldosteronism	High	High

This differentiation is vital in the treatment and management of the disease. The adrenal adenomas respond well to surgery whereas hyperplastic disease of the adrenals is generally better managed medically⁶.

In addition, pharmacological modulation of nuclear hormone receptors is a common strategy for the treatment of cardiovascular disease⁷. Therefore, determining the effects of such treatments on the RAAS is of increasing value in evaluating the safety and efficacy of new therapeutics. In addition, obese subjects often exhibit hyperaldosteronism, with increased salt sensitivity of blood pressure (BP). Systemic RAS, and aldosterone/MR activation plays a key role in the development of hypertension and organ damage in obesity¹⁴⁻¹⁶.

In summary, the precise and accurate measurement of serum aldosterone by enzyme immunoassay can be an important adjunct to a diagnostic laboratory battery for the differential diagnosis of hypertensive disease.

2 PRINCIPLE OF THE TEST

The DRG:HYBRID-XL Aldosterone Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **principle of competitive binding**.

The antibody coated wells (ACW) of the reagent cartridges are coated with a monoclonal [rabbit] antibody directed towards a unique antigenic site of the aldosterone molecule. Endogenous aldosterone of a patient sample competes with an aldosterone-horseradish peroxidase conjugate for binding to the coated antibody. After incubation the unbound conjugate is washed off.

Having added the substrate solution, the intensity of colour developed is inversely proportional to the concentration of aldosterone in the patient sample.

The concentrations of unknown samples are determined from the master curve.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. This kit can only be used in combination with the DRG:HYBRID-XL Analyzer

3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is clear and understood.
4. **Do not remove, exchange, discard, or damage any of the barcode labels provided with each kit and its components. All barcodes build an integral system for the kit lot.**
5. Respect the general safety measures for use of laboratory reagents.
6. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
7. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
9. Wear appropriate disposable gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may cause false results.
10. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
11. Do not use reagents beyond expiry date as shown on the kit labels.
12. Unused reagent cartridges must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant provided.
13. Optimal test results are only obtained when using calibrated pipettes.
14. Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
15. Some reagents contain ProClin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
16. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant amount of water and skin with soap and plenty of water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
17. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.
18. For information on hazardous substances included in the kit please refer to Safety Data Sheets.
For professional users the Safety Data Sheet for this product is available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

4.1.1 Reagent Cartridges

40 pieces containing the following:

- **Antibody Coated Well (ACW)** coated with anti-aldosterone antibody (monoclonal).
- **Enzyme Conjugate**, 160 µL, aldosterone antigen conjugated with horseradish peroxidase; Contains non-mercury preservative.
- **Substrate Solution**, 260 µL Tetramethylbenzidine (TMB).

4.1.2 Re-Calibrator 1 & 2

2 vials, 1.0 mL each, lyophilized;

For re-calibration of the quantitative DRG:HYBRID-XL Aldosterone test.

Concentrations are lot-specific.

Conversion: 1 pg/mL = 2.77 pmol/L

Re-Calibrators are standardised against the following reference material:

Code A-096 (Cerilliant)

See "Reagent Preparation".

Contain non-mercury preservative.

4.1.3 Control 1 & 2

2 vials, 1.0 mL each, lyophilized;

For control values and ranges please refer to the bar code on vial label or to the QC-Datasheet.

See "Reagent Preparation".

Contain non-mercury preservative.

4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:
Resistivity at 25 °C [MΩ·cm]: > 10
Conductivity at 25 °C [µS·cm⁻¹]: < 0.1
Total Organic Carbon/p.p.b [µg/L]: < 50
Colloids [µg/mL]: <0.05
- **[REF]** HYB-5252 DRG:HYBRID-XL Analyzer
- **[REF]** HYI-5392: *System Solution 5L*, 5000 mL;
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:
Resistivity at 25 °C [MΩ·cm]: > 1
Conductivity at 25 °C [µS·cm⁻¹]: < 1
Total Organic Carbon/p.p.b [µg/L]: < 200
Colloids [µg/mL]: <0.1)
- **[REF]** HYI-6234: *Wash Buffer*, 10x concentrate, 100 mL
- **[REF]** HYI-5395: *Needle Cleaning Solution*, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- **[REF]** HYI-5387: *Cuvettes*, (2 x 360 pieces)

For use of the *Secondary Sample Holder* for secondary tubes the following tubes are required:

- **[REF]** HYI-5391: *Sample Tubes (Secondary)*, 2500 pcs.
- **Optional:**
Reagents for determination of **Aldosterone in urine** (**[REF]** EIA-5298-URIN) – **Contents:**
 - 1) **Release Reagent**, 1 vial, 3 mL, ready to use.
Containing 1M HCl.
Avoid contact with *Release Reagent*. It may cause skin irritation.
 - 2) **Neutralization Buffer**, 1 vial, 3 mL, ready to use.
Containing Tris buffer, pH 8.5.
 - 3) **Dilution Buffer**, 2 vials, 25 mL each, ready to use.
Containing PBS.
- **Optional:**
Plastic tubes (e.g. 0.5 – 1.5 mL) for pre-treatment of urine samples

4.3 Storage Conditions

All kit components should be stored at 2 °C to 8 °C to ensure product performance until the defined expiry date.

When stored at 2 °C to 8 °C, **unopened kits** will retain reactivity until expiration date. Do not use reagents beyond this date.

- Cartridges (stored at 2 °C to 8 °C) in the supplied and unopened zip/foil bags will retain reactivity until expiration date.
- Unopened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity until expiration date.

Opened reagents and the reagent cartridges must be stored at 2 °C to 8 °C.

Once the plastic bag has been opened, care should be taken to tightly close it again along with the supplied desiccant bag.

Immediately after end of each run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped and stored at 2 °C to 8 °C.

- Unused cartridges in opened zip/foil bags (stored at 2 °C to 8 °C) will retain reactivity until expiration date, if stored as described above.
- Pierced or open cartridges must be disposed of immediately.
- Opened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity for 8 weeks.

4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 25 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

On-board stability:	24 h
---------------------	------

4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 25 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

Re-Calibrator 1 & 2

Reconstitute the lyophilized content of each vial with 1.0 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Note: *The reconstituted re-calibrators are stable for 2 months at 2 °C to 8 °C.*

For longer storage aliquot and freeze at -20 °C.

Control 1 & 2

Reconstitute the lyophilized content of each vial with 1.0 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

Note: *The reconstituted controls are stable for 2 months at 2 °C to 8 °C.*

For longer storage aliquot and freeze at -20 °C.

Wash Buffer (10x) (HYI-6234; not included in the kit)

For Wash Buffer (1x) dilute 100 mL of Wash Buffer (10x) with 900 mL ultra-pure water to a final volume of 1000 mL.

stability after dilution:	20 °C - 26 °C	4 weeks
---------------------------	---------------	---------

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum, plasma (EDTA-, heparin- or citrate plasma) or urine can be used in this assay.

A minimum of 160 µL of sample is needed for one determination.

This includes 100 µL sample and 60 µL dead volume.

Attention:

- This test was not verified with blood collection tubes of all available manufacturers.
- Sample Collection Systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- Do not use haemolytic, icteric or lipaemic specimens.
- Samples containing precipitates have to be centrifuged prior to the test run.
- Do not use heat inactivated samples.
- Do not use standards or external controls stabilized with azide.

5.1 Serum / Plasma Samples

5.1.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.1.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 4 days at 2 °C to 8 °C prior to performing the assay.

Specimens stored for a longer time (up to 12 months) should be frozen only once at -20 °C prior to the assay. Thawed samples should be inverted several times prior to testing.

5.1.3 Manual Sample Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent* and measured again as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

5.2 Urine Samples

5.2.1 Sample Collection

First clean genital area with mild disinfectant to prevent contamination. Then collect clean-catch midstream urine in an appropriate sterile container.

Since aldosterone secretion follows a circadian rhythm, an urine collection is recommended in a special cooled container over a full 24 hour period (24-hour urine).

Directly after collection, the urine should be centrifuged for 5 – 10 minutes (e.g. at 2000 g) to remove cellular debris.

Use supernatant for analyte quantification.

The **supernatant** may be stored for up to 7 days at 2 °C to 8 °C prior to assaying. Supernatant held for a longer time (up to 7 days) should be frozen at -20 °C. Thawed supernatant should be inverted several times prior to testing.

5.2.2 Protocol for Urine Sample Pre-treatment

1. Secure the desired number of vials (e.g. 0.5 – 1.5 mL plastic tubes; not included in this kit).
2. Dispense **25 µL** of **urine** with new disposable tips into appropriate tubes.
3. Dispense **25 µL Release Reagent** into each tube. Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate overnight at 2 °C to 8 °C.
5. Add **25 µL Neutralization Buffer** to each tube and mix thoroughly.
6. Add **400 µL Dilution Buffer** to each tube and mix thoroughly (This pre-treatment leads to a 1:19 dilution. Therefore the dilution factor 19 has to be taken into account for calculation of the final concentration of the urine sample.)
7. Transfer a minimum of **160 µL of pre-treated and diluted urine sample** into a Sample Tube for one determination.

5.2.3 Storage of Pre-treated Urine Samples

Pre-treated and diluted urine samples should be capped and may be stored for up to 7 days at 2 °C to 8 °C or frozen at -20 °C prior to assaying. Samples should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.2.4 Manual Sample Dilution

If in an initial assay an urine sample is found to contain more than the highest standard, the pre-treated and diluted urine sample can be further diluted with *Dilution Buffer* and measured again as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account too.

Example:

dilution 1:10: 10 µL pre-treated and diluted urine sample
+ 90 µL *Dilution Buffer* (mix thoroughly)
(final dilution factor = 19 x 10 = 190)

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents, such as controls and re-calibrators and specimens, must be allowed to come to room temperature (20 °C to 25 °C) before use. All reagents and samples must be mixed without foaming. Reagent Cartridges can be used directly after storage in the refrigerator.
- Samples, controls and re-calibrators should be measured within 2 hours in order to avoid possible evaporation effects.
- The *Secondary Sample Holder* (HY1-5437) for secondary tubes has the capacity for a maximum of 20 samples including controls and re-calibrators. They all have to be pipetted into the secondary tubes, and the respective barcodes of control/re-calibrator vials and, if available, the sample barcodes have to be read with the external barcode scanner.

6.2 Test Procedure

- **The total assay time for DRG:HYBRID-XL Aldosterone is 90 minutes.**
- To ensure proper operation of the test, the instructions in the user manual for the DRG:HYBRID-XL should strictly be followed.
- All test specific information required for the correct operation is included in the respective barcodes of the reagents. **Take care not to damage these bar codes!**
- **It is recommended to tap the bottom of the Cartridge Segments containing the reagent cartridges once on the bench before placing them on the rotor. This is to avoid foam and adhering of the liquid on the sealing of the reagent cartridge.**
- Place reagent cartridges on the rotor of the unit. The heating to 37 °C incubation temperature is performed automatically in the unit.

6.3 Calibration

Traceability:

This method was standardized against the following reference material: Code A-096 (Cerilliant)

Each DRG:HYBRID-XL reagent contains a barcode with the specific information for recalibration of the reagent lot. The Master Curve is printed as a 2-D barcode on the outer label of the kit package and on the QC-Datasheet and has to be scanned with the external barcode scanner prior to the first use of the respective kit lot.

Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRID-XL's system software.

6.4.1 Final Calculation for Urine Samples

Calculate the 24 hours excretion for each urine sample:

$$\mu\text{g}/24 \text{ h} = \mu\text{g}/\text{L} \times \text{L}/24 \text{ h}$$

Example:

Measured concentration for urine sample = 500 µg/mL
Result after correction with the dilution factor 19 = 9500 µg/mL
9500 µg/mL / 1000 = 9.5 µg/L

Total volume of 24 h-urine = 1.3 L (example)

$$9.5 \mu\text{g}/\text{L} \times 1.3 \text{ L}/24 \text{ h} = 12.35 \mu\text{g}/24 \text{ h}$$

7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

It is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analysing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyser. In addition, it is indicated to perform a Recalibration.

In case of further questions please contact your local distributor or DRG directly.

7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the QC-Datasheet added to the kit. Note that the expected values and acceptance ranges stated in the QC-Datasheet always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if re-calibration is required (if one or both internal controls are out of range)

- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

8 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and pathological values.

In a study conducted with **serum/plasma samples** of apparently healthy **adults**, using the DRG:HYBRID-XL Aldosterone, the following values are observed:

Population	n	Mean (pg/mL)	5 th Percentile (pg/mL)	95 th Percentile (pg/mL)
Supine position	55	52.3	11.5	133.1
Upright position	55	97.8	13.6	261.7
Total	110	75.0	12.3	226.3

These results correspond well to published reference ranges ^{8,9}.

In a study conducted with **serum/plasma samples** of apparently healthy **newborns and children**, using the DRG:HYBRID-XL Aldosterone, the following values are observed:

Age (years)	n	Range (pg/mL)	Mean (pg/mL)	1. Percentile (pg/mL)	Median (pg/mL)	99. Percentile (pg/mL)
Newborns	7	18.8 – 4785	1875.6	21.5	1620	4745.4
1 – 5	9	26.6 – 396.0	185.9	28.2	131.5	389.2
6 – 10	7	31.7 – 178.0	106.6	34.0	96.0	177.1
11 – 15	9	48.3 – 114.5	84.7	49.1	84.4	114.2
16 – 20	7	11.4 – 169.0	73.3	11.5	69.8	167.6

In a study conducted with apparently normal healthy **adults**, using the DRG:HYBRID-XL Aldosterone and the DRG:HYBRID-XL Renin (active) (HYE-5373), the following **Aldosterone-Renin Ratios** were determined in serum/plasma:

Ratio Aldosterone Renin (pg/mL / pg/mL)

n	Mean	Median	1 st Percentile	99 th Percentile
110	6.39	4.07	0.48	54.66

In a study conducted with **urine samples** of apparently normal healthy adults, using the DRG:HYBRID-XL Aldosterone, the following values are observed:

n	Mean (µg/24 h)	Median (µg/24 h)	5 th Percentile (µg/24 h)	95 th Percentile (µg/24 h)
40	11.34	9.40	3.55	23.01

These results correspond well to published reference ranges ⁸.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

9 LIMITATIONS OF USE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

10 PERFORMANCE CHARACTERISTICS

10.1 Assay Dynamic Range

The dynamic range of the assay is defined by the limit of blank and the maximum value of the Master Curve.

Values found below the measuring range are indicated as “< 12.329 pg/mL”.

Values found above the measuring range are indicated as “> 1000 pg/mL”.

The measuring range of the assay is between 12.329 pg/mL – 1000 pg/mL.

The linear measuring range is between 22.376 – 1000 pg/mL.

10.2 Specificity of Antibodies (Cross-Reactivity)

The following substances were tested for cross-reactivity of the assay:

Substance	Concentration Range of Spiked Substance (ng/mL)	Mean Cross-Reactivity (%)
11-Desoxy Cortisol	10 – 1000	0.0
17-OH Progesterone	1.2 – 120	0.0
21-OH Progesterone	35 – 3500	0.0
Androstenedione	0.22 – 22	0.6
Corticosterone	0.50 – 50	0.0
Cortisol	8.0 – 800	0.0
Cortisone	16 – 1600	0.0
DHEA	0.50 – 500	0.0
DHEA-S	30 – 3000	0.0
Estradiol	0.02 – 2	0.0
Estriol	1.5 – 150	0.2
Estrone	0.01 – 1	0.0
Prednisolone	35 – 3500	0.0
Prednisone	35 – 3500	0.0
Pregnenolone	35 – 3500	0.0
Progesterone	42.2 – 4220	0.0
Testosterone	0.010 – 1	0.0

10.3 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2.

The Limit of Blank (LoB) is 12.329 pg/mL.

The Limit of Detection (LoD) is 21.097 pg/mL.

The Limit of Quantification (LoQ) is 35.334 pg/mL.

10.4 Precision Performance

The precision study was designed according to CLSI guideline EP5-A3.

10.4.1 Within-run precision

Calculated from the mean value of 5 replicates on 5 days (serum/plasma) or 4 days (urine) with 1 device.

Sample	n	Mean (pg/mL)	CV (%)
P1	5	112.53	5.7
S2	5	234.67	2.8
P3	5	837.66	3.4
U4	4	180.81	5.9
U5	4	278.83	3.0
U6	4	321.31	5.5

10.4.2 Between-run precision

Calculated by 5 replicates each on 5 days (serum/plasma) or 4 days (urine) with 1 device.

Sample	n	Mean (pg/mL)	CV (%)
P1	25	112.53	5.1
S2	25	234.67	4.5
P3	25	837.66	3.3
U4	20	180.81	6.3
U5	20	278.83	3.5
U6	20	321.31	6.0

10.4.3 Between-device precision

Calculated by 5 replicates each on 5 days (serum/plasma) or 4 days (urine) with 3 devices.

Sample	n	Mean (pg/mL)	CV (%)
P1	75	104.42	11.4
S2	75	232.60	4.4
P3	75	812.00	8.2
U4	60	180.19	6.4
U5	60	285.11	5.5
U6	60	312.36	6.2

10.5 Recovery

Recovery was determined by adding four increasing concentrations of the analyte to four different patient samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was assayed and analyte concentrations of the samples were calculated from the Master Curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

Sample	P1	S2	S3	U4	U5	U6	
Concentration (pg/mL)	67.96	115.35	493.75	304.02	413.66	748.68	
Average Recovery (%)	90.0	91.3	99.4	104.9	108.6	95.2	
Range of Recovery (%)	from	85.3	85.0	95.7	96.6	103.8	93.5
	to	98.2	95.2	102.6	113.3	112.1	97.3

10.6 Linearity

Samples containing different amounts of analyte were serially diluted with *Sample Diluent*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

Serum sample	P1	S2	S3	U4	U5	U6	
Concentration (pg/mL)	182.5 0	201.0 0	698.2 2	651.0 0	876.7 4	1210.0 0	
Average Recovery (%)	104.0	108.6	94.2	93.2	101.3	107.6	
Range of Recovery (%)	from	93.2	103.6	86.9	85.2	98.8	98.4
	to	114.3	113.6	108.3	110.7	106.2	114.2

10.7 Interfering Substances

Hemoglobin (up to 4 mg/mL), bilirubin (up to 0.5 mg/mL) and triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

10.8 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of aldosterone in a sample.


10.9 High-Dose-Hook Effect

"High-Dose Hook Effect" is not detected up to 10 000 pg/mL of aldosterone.

11 METHOD COMPARISON

A comparison of DRG:HYBRiD-XL Aldosterone Test HYE-5338 (y) and DiaSorin LIAISON® Aldosterone (x) using **serum samples** gave the following correlation:

$$\begin{aligned} n &= 35 \\ r &= 0.97 \\ y &= -39.13 + 1.055x \end{aligned}$$

**12 LEGAL ASPECTS**

Only for countries where the declaration of European Conformity (CE mark) is applicable.

12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

For further information please refer to the User Manual of the DRG:HYBRiD-XL, analyser-specific application sheets, product information and package inserts of all necessary components.

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


1. Brown RD, Strott CA, and Liddle GW. Site of stimulation of Aldosterone biosynthesis by angiotensin and potassium. *J Clin Invest.* (1972), 51 (6), 1413–8.
2. Bauer JH, Gauntner WC. Effect of potassium chloride on plasma renin activity and plasma aldosterone during sodium restriction in normal man. *Kidney Int.* (1979), 15 (3): 286–93.
3. Williams GH, Dluhy RG. Aldosterone biosynthesis. Interrelationship of regulatory factors. *Am J Med.* (1972), 53 (5), 595–605.
4. Tiu SC et al. The use of aldosterone-renin ratio as a diagnostic test for primary hyperaldosteronism and its test characteristics under different conditions of blood sampling. *J Clin Endocrinol Metab.* (2005), 90 (1), 72-8.
5. Mulatero P et al. Confirmatory tests in the diagnosis of primary aldosteronism. *Horm Metab Res.* (2010), 42 (6), 406-10.
6. Quillo AR. Primary aldosteronism: results of adrenalectomy for nonsingle adenoma. *J Am Coll Surg.* (2011), 213 (1), 106-12.
7. Grossmann C and Gekle M. New aspects of rapid aldosterone signaling. *Mol Cell Endocrinology* (2009), 308 (1-2), 53-62.
8. Thomas L (editor). *Renin-Angiotensin-Aldosteron-System (RAAS)*. 8th Edition; Labor und Diagnose 2012, 1558.
9. Perschel FH et al. Rapid Screening test for primary hyperaldosteronism: ratio of plasma aldosterone to renin concentration determined by fully automated chemiluminescence immunoassays. *Clin. Chemistry* (2004); 50 (9), 1650-55.
10. Rege J et al. Steroid Biomarkers in Human Adrenal Disease. *J Steroid Biochem Mol Biol.* 2019 190: 273–280.
11. Raff H et al. Physiological Basis for the Etiology, Diagnosis, and Treatment of Adrenal Disorders: Cushing's Syndrome, Adrenal Insufficiency, and Congenital Adrenal Hyperplasia. *Compr Physiol.* 2014; 4(2): 739–769.
12. Tsilosani A, Gao C, Zhang W. Aldosterone-Regulated Sodium Transport and Blood Pressure. *Front Physiol.* 2022; 13: 770375.
13. Yozamp N et al. 2021. Intra-Individual Variability of Aldosterone Concentrations in Primary Aldosteronism: Implications for Case Detection. *Hypertension*, 2021; 77(3), 891-899.
14. Kawarazaki W and Fujita T. The Role of Aldosterone in Obesity-Related Hypertension. *American Journal of Hypertension* 2016; 29(4), 415-423.
15. Eisenhofer G et al. Reference intervals for plasma concentrations of adrenal steroids measured by LC-MS/MS: Impact of gender, age, oral contraceptives, body mass index and blood pressure status. *Clin Chim Acta* 2017; 470, 115-124.
16. Kerstens M et al. Reference Values for Aldosterone–Renin Ratios in Normotensive Individuals and Effect of Changes in Dietary Sodium Consumption. *Clin.Chem.* 2011; 57:11, 1607-1611.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français	Polski
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conforme aux normes européennes	Zgodność z normami europejskimi
	Consult instructions for use*	Gebrauchsanweisung beachten*	Consultare le istruzioni per l'uso	Consulte las Instrucciones	Consulter les instructions d'utilisation	Zapoznać się z instrukcją użytkowania
	<i>In vitro</i> diagnostic medical device*	<i>In-vitro</i> -Diagnostikum*	Per uso Diagnostica in vitro	Para uso Diagnóstico in vitro	Usage Diagnostic in vitro	Wyrób medyczny do diagnostyki in vitro
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches	Tylko do użytku w badaniach
	Catalogue number*	Artikelnummer*	No. di Cat.	Número de catálogo	Référence	Numer katalogowy
	Batch code*	Chargencode*	Lotto no	Número de lote	No. de lot	Numer LOT
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests	Zawartość przeznaczona na <n> testów
	Temperature limit*	Temperaturbegrenzung*	Temperatura di conservazione	Temperatura de conservación	Température de conservation	Przechowywać w temperaturze
	Use-by date*	Verwendbar bis*	Data di scadenza	Fecha de caducidad	Date limite d'utilisation	Data przydatności
	Manufacturer*	Hersteller*	Fabbricante	Fabricante	Fabricant	Producent
	Caution*	Achtung*				UWAGA
<i>Distributed by</i>	Distributor	Vertreiber	Distributore	Distribuidor	Distributeur	Dystrybutor
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu	Zawartość
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Numéro	Objętość / Numer
<i>Reagent Cartridge</i>	Reagent Cartridge	Reagenzien-Cartridge	Cartucce di reagente	Cartucho de reactivo	Cartouche de réactif	Kartridż z odczynnikami
<i>Re-Calibrator</i>	Re-Calibrator	Re-Kalibrator	Re-Calibratore	Re-Calibrador	Re-calibrateurs	Re-kalibrator
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle	Kontrola
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique	Koniugat enzymatyczny
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complex enzimático		Enzymu kompleks
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat	Roztwór Substratu
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungs-medium	Diluyente dei campioni	Solución para dilución de la muestra	Diluant d'échantillon	Rozcieńczalnik Próbk
<i>Reagent 1</i>	Reagent 1	Reagenz 1	Reagente 1	Reactivo 1	Réactif 1	Reagent 1
<i>Reagent 2</i>	Reagent 2	Reagenz 2	Reagente 2	Reactivo 2	Réactif 2	Reagent 2
<i>Wash Buffer</i>	Wash Buffer	Waschpuffer	Tampone di lavaggio	Tampón de lavado	Tampon de lavage	Bufor Płuczący
<i>System Solution</i>	System Solution	Systemlösung		Solución de sistema	Solution système	Roztwór systemowy
<i>Needle Cleaning Solution</i>	Needle Cleaning Solution	Nadel-Reinigungslösung	Soluzione lavaggio ago	Solución de lavado de la aguja	Solution de nettoyage des aiguilles	Roztwór Czyszczący Igły
<i>Denaturation Buffer</i>	Denaturation Buffer	Denaturierungspuffer	Tampone die denaturazione	Tampón de denaturalización	Tampon de dénaturation	
<i>Neutralization Buffer</i>	Neutralization Buffer	Neutralisierungspuffer	Tampone di neutralizzazione	Tampón de neutralización	Tampon de neutralisation	
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai	Bufor testowy
<i>Secondary Sample Holder</i>	Secondary Sample Holder	Sekundärprobenhalter	Sostegno secondario dei campioni	Soporte de muestras para tubos secundarios	Support d'échantillons secondaires	Statyw na próbki drugorzędowe
<i>Sample Tubes (Secondary)</i>	Sample Tubes (Secondary)	Sekundärprobenröhrchen	Tubetti di campioni secondari	Tubos de muestra secundarios	Tubes échantillon secondaires	Probówki drugorzędowe
<i>Dilution Cartridge</i>	Dilution Cartridge				Cartouches de dilution	Kartridż do rozcieńczeń
<i>Vial Adapter</i>	Vial Adapter	Fläschchen-Adapter	l'adattatore dei tubetti	Adaptador de tubos	l'Adaptateur de flacons	Nakładki na butelki dla kontroli i kalibratorów