
# DRG:HYBRID•XL®

## Instructions for Use

# LH

**IVD**

**REF**

**HYE-5362**

**40 tests**

**DRG Instruments GmbH, Germany**  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)


**Distributed by:**

**DRG International, Inc., USA**  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**ASSAY PROTOCOL BARCODE (APB)**

(Version 2.61 Software or later / ab Softwareversion 2.61 oder höher / Versione Software 2.61 o superiore /  
 Versión de Software 2.61 o posterior / Version de logiciel 2.61 et supérieure /  
 Wersja oprogramowania 2.61 lub późniejsza)

HYE-5362 – v2.61.1

The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the SCAN NEW LOT page.  
 Der Barcode muss in dem Menü „NEUES LOT SCANNEN“ eingelesen werden, um das Protokoll in der DRG:HYBRiD-XL-Software zu installieren.

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Please refer to section **3: Routine Procedures: "Installing a new assay product"** of the User Manual v2.60 or later.

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1	INTRODUCTION.....	4
2	PRINCIPLE OF THE TEST.....	4
3	WARNINGS AND PRECAUTIONS.....	4
4	REAGENTS.....	5
5	SPECIMEN COLLECTION AND PREPARATION.....	5
6	ASSAY PROCEDURE.....	6
7	QUALITY CONTROL.....	6
8	EXPECTED NORMAL VALUES.....	6
9	LIMITATIONS OF USE.....	6
10	PERFORMANCE CHARACTERISTICS.....	7
11	METHOD COMPARISON.....	7
12	LEGAL ASPECTS.....	8

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## 1 INTRODUCTION

### 1.1 Intended Use

The **DRG:HYBRID-XL LH** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of luteinizing hormone (LH) in serum and plasma (EDTA-, heparin- or citrate plasma). Only for use with the DRG:HYBRID-XL Analyzer.

### 1.2 Summary and Explanation

Luteinizing hormone (LH) is produced in both men and women from the anterior pituitary gland in response to luteinizing hormone-releasing hormone (LH-RH or Gn-RH), which is released by the hypothalamus (1-3). LH, also called interstitial cell-stimulating hormone (ICSH) in men, is a glycoprotein with a molecular weight of approximately 30,000 daltons (4). It is composed of two non-covalently associated dissimilar amino acid chains, alpha and beta (5). The alpha chain is similar to that found in human thyroid-stimulating hormone (TSH), follicle stimulating hormone (FSH), and human chorionic gonadotropin (hCG), while the beta subunits differ, which accounts for their immunological differentiation (6-8).

The basal secretion of LH in men is episodic and has the primary function of stimulating the interstitial cells (Leydig cells) to produce androgenic hormone. The primary and most well-known androgen is testosterone. The variation in LH concentrations in women is subject to the complex ovulatory cycle of healthy menstruating women, and depends upon a sequence of hormonal events along the gonado-hypothalamic-pituitary axis. The decrease in progesterone and estradiol levels from the preceding ovulation initiates each menstrual cycle (9,10). As a result of the decrease in hormone levels, the hypothalamus increases the secretion of gonadotropin-releasing factors (GnRF), which in turn stimulates the pituitary to increase FSH production and secretion (4). Increasing FSH levels stimulate several follicles during the follicular phase, one of these will mature to contain the egg. As the follicle develops, estradiol secretion rises constantly, peaking on day 12 or 13 of a normal cycle. LH is released as a result of high estradiol levels because of direct stimulation of the pituitary and increasing GnRF and FSH levels. These events constitute the pre-ovulatory phase (11). Ovulation occurs approximately 12 to 18 hours after the LH reaches a maximum level. After the egg is released, corpus luteum is formed which secretes progesterone and estrogen - two feedback regulators of LH (3,10).

The luteal phase rapidly follows this ovulatory phase, and is characterized by high progesterone levels, a second estradiol increase, and low LH and FSH levels (12). Low LH and FSH levels are the result of the negative feedback effects of estradiol and progesterone on the hypothalamic-pituitary axis.

After conception, the developing embryo produces hCG, which causes the corpus luteum to continue producing progesterone and estradiol. The corpus luteum regresses if pregnancy does not occur, and the corresponding drop in progesterone and estradiol levels results in menstruation. The hypothalamus initiates the menstrual cycle again as a result of these low hormone levels (12).

Patients suffering from hypogonadism show increased concentrations of serum LH. A decrease in steroid hormone production in females is a result of immature ovaries, primary ovarian failure, polycystic ovary disease, or menopause; in these cases, LH secretion is not regulated (10,13). A similar loss of regulatory hormones occurs in males when the testes develop abnormally or in anorchia. High concentrations of LH may also be found in primary testicular failure and Klinefelter syndrome, although LH levels will not necessarily be elevated if the secretion of androgens continues. Increased concentrations of LH are also present during renal failure, cirrhosis, hyperthyroidism, and severe starvation (10,14).

A lack of secretion by the anterior pituitary may cause lower LH levels. As may be expected, low levels may result in infertility in both males and females. Low levels of LH may also be due to the decreased secretion of GnRH by the hypothalamus, although the same effect may be seen by a failure of the anterior pituitary to respond to GnRH stimulation. Low LH values may therefore indicate some dysfunction of the pituitary or hypothalamus, but the actual source of the problem must be confirmed by other tests (10).

In the differential diagnosis of hypothalamic, pituitary, or gonadal dysfunction, assays of LH concentration are routinely performed in conjunction with FSH assays since their roles are closely interrelated. Furthermore, the hormone levels are used to determine menopause, pinpoint ovulation, and monitor endocrine therapy.

## 2 PRINCIPLE OF THE TEST

The DRG:HYBRID-XL LH Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **sandwich principle**.

The antibody coated wells (ACW) of the reagent cartridges are coated with a monoclonal (mouse) antibody directed towards a unique antigenic site of the LH molecule. An aliquot of patient sample containing endogenous LH is incubated in the coated well with enzyme conjugate, which is an anti-LH antibody conjugated with biotin (*Enzyme Conjugate*).

After incubation the unbound conjugate is washed off.

Bound biotinylated antibodies are detected by consequent incubation of enzyme complex (*Enzyme Complex*), which is Streptavidin-HRP.

After washing off unbound complex, TMB substrate is added.

The intensity of colour developed is proportional to the concentration of LH in the patient sample.

## 3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only.  
For professional use only.
2. This kit can only be used in combination with the DRG:HYBRID-XL Analyzer
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is clear and understood.
4. **Do not remove, exchange, discard or damage any of the barcode labels provided with each kit and its components. All barcodes build an integral system for the kit lot.**
5. Respect the general safety measures for use of laboratory reagents.
6. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
7. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
9. Wear appropriate disposable gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may cause false results.
10. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
11. Do not use reagents beyond expiry date as shown on the kit labels.
12. Unused reagent cartridges must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant provided.
13. Optimal test results are only obtained when using calibrated pipettes.
14. Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
15. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
16. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant amount of water and skin with soap and plenty of water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
17. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.
18. For information on hazardous substances included in the kit please refer to Safety Data Sheets.  
For professional users the Safety Data Sheet for this product is available upon request directly from DRG.

**4 REAGENTS**

**4.1 Reagents provided**

**4.1.1 Reagent Cartridges**

40 pieces containing the following:

- **Antibody Coated Well (ACW)**  
coated with anti-LH antibody (monoclonal).
- **Enzyme Conjugate**, 170 µL,  
Anti-LH antibody conjugated with biotin;  
Contains non-mercury preservative.
- **Enzyme Complex**, 170 µL,  
Streptavidin-HRP  
Contains non-mercury preservative.
- **Substrate Solution**, 270 µL  
Tetramethylbenzidine (TMB).

**4.1.2 Re-Calibrator 1 & 2**

2 vials, 1 mL each, lyophilized;

For re-calibration of the quantitative DRG:HYBRiD-XL LH test.

Concentrations are lot-specific.

Conversion: 1 mIU/mL = 0.176 ng/mL

Re-Calibrators are standardised against the following reference material: WHO International Standard LH human (NIBSC code 81/535) See "Reagent Preparation".

Contain non-mercury preservative.

**4.1.3 Control 1 & 2**

2 vials, 1 mL each, lyophilized;

For control values and ranges please refer to the bar code on vial label or to the QC-Datasheet.

See "Reagent Preparation".

Contain non-mercury preservative.

**4.2 Materials required but not provided**

- General needed laboratory equipment
- Ultra-pure water  
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:  
Resistivity at 25 °C [MΩ·cm]: > 10  
Conductivity at 25 °C [µS·cm<sup>-1</sup>]: < 0.1  
Total Organic Carbon/p.p.b.[µg/L] : < 50  
Colloids [µg/mL]: <0.05
- [REF] HYB-5252 DRG:HYBRiD-XL Analyzer
- [REF] HYI-5392: System Solution 5L, 5000 mL;  
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:  
Resistivity at 25 °C [MΩ·cm]: > 1  
Conductivity at 25 °C [µS·cm<sup>-1</sup>]: < 1  
Total Organic Carbon/p.p.b.[µg/L] : < 200  
Colloids [µg/mL]: <0.1)
- [REF] HYI-6234: Wash Buffer, 10x concentrate, 100 mL
- [REF] HYI-5395: Needle Cleaning Solution, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- [REF] HYI-5387: Cuvettes, (2 x 360 pieces)

For use of the Secondary Sample Holder for secondary tubes the following tubes are required:

- [REF] HYI-5391: Sample Tubes (Secondary), 2500 pcs.

**4.3 Storage Conditions**

All kit components should be stored at 2 °C to 8 °C to ensure product performance until the defined expiry date.

When stored at 2 °C to 8 °C, **unopened kits** will retain reactivity until expiration date. Do not use reagents beyond this date.

- Cartridges (stored at 2 °C to 8 °C) in the supplied and unopened zip/foil bags will retain reactivity until expiration date.
- Unopened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity until expiration date.

**Opened reagents and the reagent cartridges** must be stored at 2 °C to 8 °C.

Once the plastic bag has been opened, care should be taken to tightly close it again along with the supplied desiccant bag.

**Immediately after end of each** run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped and stored at 2 °C to 8 °C.

- Unused cartridges in opened zip/foil bags (stored at 2 °C to 8 °C) will retain reactivity until expiration date, if stored as described above.
- Pierced or open cartridges must be disposed of immediately.
- Opened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity for 8 weeks. (check in Technical Documentation)

**4.3.1 On-board Stability**

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 25 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

On-board stability:	8 h
---------------------	-----

**4.4 Reagent Preparation**

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 25 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

**Re-Calibrator 1 & 2**

Reconstitute the lyophilized content of each vial with 1 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

**Note:** The reconstituted re-calibrators are stable for 2 months at 2 °C to 8 °C.

For longer storage aliquot and freeze at -20 °C.

**Control 1 & 2**

Reconstitute the lyophilized content of each vial with 1 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

**Note:** The reconstituted controls are stable for 2 months at 2 °C to 8 °C.

For longer storage aliquot and freeze at -20 °C.

**Wash Buffer (10x) (HYI-6234; not included in the kit)**

For Wash Buffer (1x) dilute 100 mL of Wash Buffer (10x) with 900 mL ultra-pure water to a final volume of 1000 mL.

stability after dilution:	20 °C - 26 °C	4 weeks
---------------------------	---------------	---------

**4.5 Disposal of the Kit**

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet.

**4.6 Damaged Test Kits**

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

**5 SPECIMEN COLLECTION AND PREPARATION**

Serum or plasma (EDTA-, heparin- or citrate plasma) can be used in this assay.

**A minimum of 85 µL of sample is needed for one determination.** This includes 25 µL sample and 60 µL dead volume.

**Attention:**

- This test was not verified with blood collection tubes of all available manufacturers.
- Sample Collection Systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- Do not use haemolytic, icteric or lipaemic specimens.
- Samples containing precipitates have to be centrifuged prior to the test run.
- Do not use heat inactivated samples.
- Do not use standards or external controls stabilized with azide.

**5.1 Specimen Collection**

**Serum:**

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

**Plasma:**

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

**5.2 Specimen Storage and Preparation**

Specimens should be capped and may be stored for up to 7 days at 2 °C to 8 °C prior to performing the assay. Specimens stored for a longer time (up to 2 months) should be frozen only twice at -20 °C prior to the assay. Thawed samples should be inverted several times prior to testing.

**5.3 Specimen Dilution**

**5.3.1 Manual Sample Dilution**

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent\** and reassayed as described in Assay Procedure. For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 10 µL sample + 90 µL *Sample Diluent* (mix thoroughly)
- b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL *Sample Diluent* (mix thoroughly).

\* *Sample Diluent for manual dilution is not included in this kit, but can be ordered on request (REF HYE-5362-DIL, 20 mL).*

**6 ASSAY PROCEDURE**

**6.1 General Remarks**

- All reagents, such as controls and re-calibrators and specimens, must be allowed to come to room temperature (20 °C to 25 °C) before use. All reagents and samples must be mixed without foaming. Reagent Cartridges can be used directly after storage in the refrigerator.
- Samples, controls and re-calibrators should be measured within 2 hours in order to avoid possible evaporation effects.
- The *Secondary Sample Holder* (HY1-5437) for secondary tubes has the capacity for a maximum of 20 samples including controls and re-calibrators. They all have to be pipetted into the secondary tubes, and the respective barcodes of control/re-calibrator vials and, if available, the sample barcodes have to be read with the external barcode scanner.

**6.2 Test Procedure**

- **The total assay time for DRG:HYBRiD-XL LH is 90 minutes.**
- To ensure proper operation of the test, the instructions in the user manual for the DRG:HYBRiD-XL should strictly be followed.
- All test specific information required for the correct operation is included in the respective barcodes of the reagents. **Take care not to damage these bar codes!**
- **It is recommended to tap the bottom of the Cartridge Segments containing the reagent cartridges once on the bench before placing them on the rotor. This is to avoid foam and adhering of the liquid on the sealing of the reagent cartridge.**
- Place reagent cartridges on the rotor of the unit. The heating to 37 °C incubation temperature is performed automatically in the unit.

**6.3 Calibration**

**Traceability:**

This method was standardized against the WHO International Standard LH human (NIBSC code 81/535)

Each DRG:HYBRiD-XL reagent contains a barcode with the specific information for recalibration of the reagent lot. The Master Curve is printed as a 2-D barcode on the outer label of the kit package and on the QC-Datasheet and has to be scanned with the external barcode scanner prior to the first use of the respective kit lot.

**Recalibration is recommended:**

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.

- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

**6.4 Calculation of Results**

The analyte concentrations are calculated automatically by the DRG:HYBRiD-XL's system software.

**7 QUALITY CONTROL**

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

It is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results. Apply appropriate statistical methods for analysing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyser. In addition, it is indicated to perform a Recalibration. In case of further questions please contact your local distributor or DRG directly.

**7.1 Internal Controls**

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the QC-Datasheet added to the kit. Note that the expected values and acceptance ranges stated in the QC-Datasheet always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if re-calibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

**7.2 External Controls**

Use controls at both normal and pathological levels. The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits. Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

**8 EXPECTED NORMAL VALUES**

It is strongly recommended that each laboratory should determine its own normal and pathological values.

In a study conducted with apparently healthy subjects, using the DRG:HYBRiD-XL LH the following data were observed:

Population	Females				Males
	cycle phase			post-menopausal	
	follicular	ovulation	luteal		
n	19	24	19	17	57
Range (mIU/mL)	3.11 - 14.10	3.38 - 51.40	<1.61 - 9.77	2.24 - 55.85	2.01 - 12.30
Mean (mIU/mL)	8.24	17.95	4.90	27.79	7.03
2.5 <sup>th</sup> - 97.5 <sup>th</sup> Percentile (mIU/mL)	3.23 - 13.34	3.39 - 44.16	<1.61 - 9.62	3.32 - 53.97	0.82 - 11.64
Median (mIU/mL)	8.83	16.65	4.88	26.60	7.01

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

**9 LIMITATIONS OF USE**

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

**10 PERFORMANCE CHARACTERISTICS**

**10.1 Assay Dynamic Range**

The dynamic range of the assay is defined by the limit of detection and the maximum value of the Master Curve.

Values found below the measuring range are indicated as “< 1.61 mIU/mL”.

Values found above the measuring range are indicated as “> 200 mIU/mL”.

The measuring range of the assay is between 1.61 mIU/mL – 200 mIU/mL.

**10.2 Specificity of Antibodies (Cross-Reactivity)**

The following substances were tested for cross-reactivity of the assay:

Sample	Expected conc. ng/mL	Cross-reactivity %
FSH	166.7	2.8
Prolactin	200.0	0.0
HCG	109.2	0.0
Beta-HCG	13.3	0.0
TSH	5217.4	0.0

**10.3 Sensitivity**

The sensitivity study was designed according to CLSI guideline EP17-A2.

The Limit of Blank (LoB) is 0.9 mIU/mL.

The Limit of Detection (LoD) is 1.61 mIU/mL.

The Limit of Quantification (LoQ) is 3.03 mIU/mL.

**10.4 Precision Performance**

The precision study was designed according to CLSI guideline EP5-A2.

**10.4.1 Within-Run (Intra Assay)**

Samples were analysed 2 times daily in duplicates for 20 days on 3 different instruments. CV (%) is given as mean of 40 duplicate CVs.

The within-run variability is shown below:

Sample	n	Mean (mIU/mL)	CV (%)
1	40	6.0	6.4
2	40	8.0	3.9
3	40	14.4	3.3
4	40	74.9	5.3

**10.4.2 Total Precision (Inter Assay)**

Samples were analysed 2 times daily in duplicates for 20 days (n=80) on 3 different instruments. The total precision is shown below:

Sample	n	Mean (mIU/mL)	CV (%)
1	80	6.0	13.1
2	80	8.0	7.9
3	80	14.4	6.0
4	80	74.9	10.2

**10.4.3 Inter-Lot**

The between-lots variation was determined by 6 measurements of 4 samples with 3 different kit lots.

Sample	n	Mean (mIU/mL)	CV (%)
1	18	6.2	11.0
2	18	8.9	6.4
3	18	15.9	4.9
4	18	86.7	12.2

**10.5 Recovery**

Recovery was determined by adding four increasing concentrations of the analyte to four different patient samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was assayed and analyte concentrations of the samples were calculated from the Master Curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

Sample	1	2	3	4
Concentration (mIU/mL)	5.1	7.7	14.4	42.8
Average Recovery (%)	93.4	99.4	99.0	104.6
Range of Recovery (%)	from	87.1	92.2	93.2
	to	97.7	109.8	109.6

**10.6 Linearity**

Four samples containing different amounts of analyte were serially diluted with *Sample Diluent*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

Sample	1	2	3	4
Concentration (mIU/mL)	90.7	136.5	108.0	138.0
Average Recovery (%)	104.1	96.8	108.4	102.0
Range of Recovery (%)	from	102.4	95.8	100.9
	to	105.8	98.6	114.2

**10.7 Interfering Substances**

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results. Until today no substances (drugs) are known to us, which have an influence to the measurement of LH in a sample.

**10.8 High-Dose-Hook Effect**

Hook effect was not observed in this test up to a concentration of 4000 mIU/mL of LH.


**11 METHOD COMPARISON**

A comparison of DRG:HYBRID-XL LH Test HYE-5362 (y) and Roche cobas® 8000 (+ e 602 module) (x) using clinical samples gave the following correlation:

$$n = 52$$

$$r = 0.973$$

$$y = 0.934x + 2.094$$

---

## 12 LEGAL ASPECTS

**Only for countries where the declaration of European Conformity (CE mark) is applicable.**

### 12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

### 12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

**For further information please refer to the User Manual of the DRG:HYBRiD-XL, analyser-specific application sheets, product information and package inserts of all necessary components.**

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