
# DRG:HYBRID•XL®

Instructions for Use

## Renin (active)

IVD

CE

REF

HYE-5373

40 tests

**DRG Instruments GmbH**, Germany  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)


Distributed by:

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**ASSAY PROTOCOL BARCODE (APB)**

(Version 2.61 Software or later / ab Softwareversion 2.61 oder höher / Versione Software 2.61 o superiore /  
 Versión de Software 2.61 o posterior / Version de logiciel 2.61 et supérieure /  
 Wersja oprogramowania 2.61 lub późniejsza)

HYE-5373 – v2.61.1


The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the SCAN NEW LOT page.  
 Der Barcode muss in dem Menü „NEUES LOT SCANNEN“ eingelesen werden, um das Protokoll in der DRG:HYBRiD-XL-Software zu installieren.

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Please refer to section 3: **Routine Procedures: "Installing a new assay product"** of the User Manual v2.60 or later.

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|    |  |   |
|----|--|---|
| 1  | INTRODUCTION.....                        | 4 |
| 2  | PRINCIPLE OF THE TEST.....               | 4 |
| 3  | WARNINGS AND PRECAUTIONS.....            | 4 |
| 4  | REAGENTS.....                            | 5 |
| 5  | SPECIMEN COLLECTION AND PREPARATION..... | 5 |
| 6  | ASSAY PROCEDURE.....                     | 6 |
| 7  | QUALITY CONTROL.....                     | 6 |
| 8  | EXPECTED NORMAL VALUES.....              | 7 |
| 9  | LIMITATIONS OF USE.....                  | 7 |
| 10 | PERFORMANCE CHARACTERISTICS.....         | 7 |
| 11 | METHOD COMPARISON.....                   | 8 |
| 12 | LEGAL ASPECTS.....                       | 8 |

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## 1 INTRODUCTION

### 1.1 Intended Use

The **DRG:HYBRiD-XL Renin (active)** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of renin in serum or EDTA plasma.

Only for use with the DRG:HYBRiD-XL Analyzer.

### 1.2 Summary and Explanation

Renin is an enzyme (Mw of 37 kDa) that belongs to the aspartic acid protease family. Renin is a member of Renin-Angiotensin-Aldosterone System (RAAS) that controls blood pressure, renal blood flux, glomerular filtration, and sodium/potassium homeostasis.

Renin is produced constitutively as prorenin, an inactive precursor with 386 amino acids, in the juxtaglomerular cells of the kidney (1). In response to low intra-renal blood pressure, reduced sodium reabsorption, hypokalemia or activity of the sympathetic nervous system, active renin can be released either from a depot in the kidney or generated from prorenin by cleavage of 46 amino acids at the N-terminus of prorenin (2,3). Prorenin secretion into the blood is continuous, in contrast to the tightly controlled release of renin, and blood concentration of prorenin is approx. 100-fold higher than active renin (4,5). After release and activation, soluble renin mediates cleavage of the  $\alpha_2$ -globulin angiotensinogen into the precursor peptide angiotensin I, which ultimately is processed by angiotensin converting enzyme (ACE) to the octapeptide angiotensin II. All actions of angiotensin II are mediated by the G protein-coupled angiotensin type 1 (AT1) and angiotensin type 2 (AT2) receptors (6). Direct physiological effects of Angiotensin II include vasoconstriction, increase of tubular reabsorption of sodium and chloride, water retention, and release of the hormones aldosterone from adrenal cortex, antidiuretic hormone (ADH, Vasopressin) from posterior pituitary, and adrenocorticotropic hormone (ACTH, Corticotropin) from anterior pituitary. Release of these hormones further supports sodium retention and secretion of potassium/H<sup>+</sup> in the kidney, and increases thirst sensation and the desire for salt through the subfornical organ of the brain (7,8). In a negative feedback loop, renin secretion is reduced by high concentration of angiotensin II (9), and release of aldosterone is lowered by potassium depletion (10). Beside the action of soluble renin, binding of renin and prorenin to the membrane-bound renin receptor ATP6AP2 in brain, heart, placenta, liver, kidney and pancreas enhances efficiency of angiotensinogen cleavage and induces angiotensin-independent intracellular effects by activating mitogen activated kinases ERK1 and ERK2 (11).

Plasma renin is a good index for the activity of the RAAS. In case of dysfunction of the RAAS, the Renin assay will allow clinical implications for diagnosis, treatment, and follow up. Active renin should be measured in:

- Diagnosis of hypertension (high blood pressure: if diastolic blood pressure is > 90 mm Hg and systolic blood pressure is > 140 mm Hg; guideline of the European Society of Cardiology and the European Society of Hypertension)
- Differential diagnosis of hyperaldosteronism (primary hyperaldosteronism, secondary hyperaldosteronism with or without hypertension, pseudo-hyperaldosteronism)
- Diagnosis of isolated deficit in mineral corticoids
- Differential diagnosis of hypokalemia (secondary hyperaldosteronism or primary hypermineralcorticism)
- Detection of Renin producing tumors in the kidney
- Monitoring of glucocorticoid therapy
- Diagnosis of insufficient response to antihypertensive treatment

## 2 PRINCIPLE OF THE TEST

The DRG:HYBRiD-XL Renin (active) Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **sandwich principle**.

The antibody coated wells (ACW) of the reagent cartridges are coated with a monoclonal [mouse] antibody directed towards a unique antigenic site of the renin molecule. An aliquot of patient sample containing endogenous renin is incubated in the coated well with enzyme conjugate, which is an anti-renin antibody conjugated with horseradish peroxidase. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase conjugate is proportional to the concentration of renin in the sample.

Having added the substrate solution, the intensity of colour developed is proportional to the concentration of renin in the patient sample.

## 3 WARNINGS AND PRECAUTIONS

1. This kit is for *in vitro* diagnostic use only.  
For professional use only.
2. This kit can only be used in combination with the DRG:HYBRiD-XL Analyzer
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is clear and understood.
4. **Do not remove, exchange, discard or damage any of the barcode labels provided with each kit and its components. All barcodes build an integral system for the kit lot.**
5. Respect the general safety measures for use of laboratory reagents.
6. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
7. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
9. Wear appropriate disposable gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may cause false results.
10. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
11. Do not use reagents beyond expiry date as shown on the kit labels.
12. Unused reagent cartridges must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant provided.
13. Optimal test results are only obtained when using calibrated pipettes.
14. Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
15. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
16. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant amount of water and skin with soap and plenty of water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
17. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.
18. For information on hazardous substances included in the kit please refer to Safety Data Sheets.  
For professional users the Safety Data Sheet for this product is available upon request directly from DRG.

## 4 REAGENTS

### 4.1 Reagents provided

#### 4.1.1 Reagent Cartridges

40 pieces containing the following:

- **Antibody Coated Well (ACW)**  
coated with anti-renin antibody (monoclonal).
- **Enzyme Conjugate**, 160 µL,  
Anti- renin antibody conjugated with horseradish peroxidase;  
Contains non-mercury preservative.
- **Assay Buffer**, 210 µL  
Contains non-mercury preservative.
- **Substrate Solution**, 260 µL  
Tetramethylbenzidine (TMB).

#### 4.1.2 Re-Calibrator 1 & 2

2 vials, 1 mL each, lyophilized;

For re-calibration of the quantitative DRG:HYBRID-XL Renin (active) test.

Concentrations are lot-specific.

Conversion: 1 pg/mL = 1.44 µIU/mL

See "Reagent Preparation".

Contain non-mercury preservative.

#### 4.1.3 Control 1 & 2

2 vials, 1 mL each, lyophilized;

For control values and ranges please refer to the bar code on vial label or to the QC-Datasheet.

See "Reagent Preparation".

Contain non-mercury preservative.

### 4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water  
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:  
Resistivity at 25 °C [MΩ·cm]: > 10  
Conductivity at 25 °C [µS·cm<sup>-1</sup>]: < 0.1  
Total Organic Carbon/p.p.b.[µg/L] : < 50  
Colloids [µg/mL]: <0.05
- **REF** HYB-5252 DRG:HYBRID-XL Analyzer
- **REF** HYI-5392: *System Solution 5L*, 5000 mL;  
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:  
Resistivity at 25 °C [MΩ·cm]: > 1  
Conductivity at 25 °C [µS·cm<sup>-1</sup>]: < 1  
Total Organic Carbon/p.p.b.[µg/L] : < 200  
Colloids [µg/mL]: <0.1)
- **REF** HYI-6234: *Wash Buffer*, 10x concentrate, 100 mL
- **REF** HYI-5395: *Needle Cleaning Solution*, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- **REF** HYI-5387: *Cuvettes*, (2 x 360 pieces)

For use of the *Secondary Sample Holder* for secondary tubes the following tubes are required:

- **REF** HYI-5391: *Sample Tubes (Secondary)*, 2500 pcs.

### 4.3 Storage Conditions

All kit components should be stored at 2 °C to 8 °C to ensure product performance until the defined expiry date.

When stored at 2 °C to 8 °C, **unopened kits** will retain reactivity until expiration date. Do not use reagents beyond this date.

- Cartridges (stored at 2 °C to 8 °C) in the supplied and unopened zip/foil bags will retain reactivity until expiration date.
- Unopened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity until expiration date.

**Opened reagents and the reagent cartridges** must be stored at 2 °C to 8 °C.

Once the plastic bag has been opened, care should be taken to tightly close it again along with the supplied desiccant bag.

**Immediately after end of each run** the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped and stored at 2 °C to 8 °C.

- Unused cartridges in opened zip/foil bags (stored at 2 °C to 8 °C) will retain reactivity until expiration date, if stored as described above.
- Pierced or open cartridges must be disposed of immediately.
- Opened Re-Calibrators and Controls (stored at 2 °C to 8 °C) will retain reactivity for 8 weeks. (check in Technical Documentation)

#### 4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 25 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

|                     |     |
|---------------------|-----|
| On-board stability: | 8 h |
|---------------------|-----|

### 4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 25 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

#### Re-Calibrator 1 & 2

Reconstitute the lyophilized content of each vial with 1 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

**Note:** *The reconstituted re-calibrators are stable for 2 months at 2 °C to 8 °C.*

*For longer storage the reconstituted re-calibrators should be aliquoted and stored at -20 °C.*

#### Control 1 & 2

Reconstitute the lyophilized content of each vial with 1 mL ultra-pure water and wait for at least 10 minutes at room temperature. Mix several times before use.

**Note:** *The reconstituted controls are stable for 2 months at 2 °C to 8 °C. For longer storage the reconstituted controls should be aliquoted and stored at -20 °C.*

#### Wash Buffer (10x) (HYI-6234; not included in the kit)

For Wash Buffer (1x) dilute 100 mL of Wash Buffer (10x) with 900 mL ultra-pure water to a final volume of 1000 mL.

|                           |               |         |
|---------------------------|---------------|---------|
| stability after dilution: | 20 °C - 26 °C | 4 weeks |
|---------------------------|---------------|---------|

### 4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet.

### 4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

## 5 SPECIMEN COLLECTION AND PREPARATION

Serum or EDTA plasma can be used in this assay.

**A minimum of 110 µL of sample is needed for one determination.**

This includes 50 µL sample and 60 µL dead volume.

#### Attention:

- This test was not verified with blood collection tubes of all available manufacturers.
- Sample Collection Systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- Do not use haemolytic, icteric or lipaemic specimens.
- Samples containing precipitates have to be centrifuged prior to the test run.
- Do not use heat inactivated samples.
- Do not use standards or external controls stabilized with azide.

Conditions under which samples are collected must be carefully controlled, since a number of physiological factors can influence the renin secretion. These include:

- Posture: the patient must have been lying down for more than 1 hour or upright for more than 1 hour
- Daily Renin oscillations: sampling is to be done between 7 AM and 10 AM if possible.
- Diet: sodium content in the diet must be known and eventually verified by the measurement of natriuria over a period of 24 hours
- Medication: the level of active renin can be affected by antihypertensive medication (e.g. diuretics, ACE inhibitors, beta adrenergic blocking agents, vasodilators, renin inhibitors)
- Pregnancy: the level of inactive and active renin increases during pregnancy
- Menstrual cycle: the level of active renin increases on the second phase of the cycle (sampling is to be done if possible during the first phase)
- Age: active renin level decreases with age

#### NOTE:

Samples from tumor patients may contain elevated levels of Renin.

### 5.1 Specimen Collection

#### Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

#### Plasma:

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

### 5.2 Specimen Storage and Preparation

Specimens should be capped and *stored at room temperature and NOT stored at 2 °C to 8 °C prior to processing*, since cryo-activation of prorenin may occur in the temperature range of 2 °C to 8 °C, giving false positive active renin values (12,13).

*If samples cannot be tested within 4 hours of primary collection, store frozen at -20 °C or below.*

*It is recommended to rapidly freeze and thaw processed samples avoiding the temperature range of 2 °C to 8 °C.*

Thawed samples should be inverted several times prior to testing.

### 5.3 Specimen Dilution

#### 5.3.1 Manual Sample Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent\** and reassayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

#### Example:

- a) dilution 1:2: 75 µL sample + 75 µL *Sample Diluent* (mix thoroughly)
- b) dilution 1:5: 30 µL sample + 120 µL *Sample Diluent* (mix thoroughly).

\* *Sample Diluent for manual dilution is not included in this kit, but can be ordered on request (REF HYE-5373-DIL, 20 mL).*

## 6 ASSAY PROCEDURE

### 6.1 General Remarks

- All reagents, such as controls and re-calibrators and specimens, must be allowed to come to room temperature (20 °C to 25 °C) before use. All reagents and samples must be mixed without foaming. Reagent Cartridges can be used directly after storage in the refrigerator.
- Samples, controls and re-calibrators should be measured within 2 hours in order to avoid possible evaporation effects.
- The *Secondary Sample Holder* (HY1-5437) for secondary tubes has the capacity for a maximum of 20 samples including controls and re-calibrators. They all have to be pipetted into the secondary tubes, and the respective barcodes of control/re-calibrator vials and, if available, the sample barcodes have to be read with the external barcode scanner.

### 6.2 Test Procedure

- **The total assay time for DRG:HYBRID-XL Renin (active) is 150 minutes.**
- To ensure proper operation of the test, the instructions in the user manual for the DRG:HYBRID-XL should strictly be followed.
- All test specific information required for the correct operation is included in the respective barcodes of the reagents. **Take care not to damage these bar codes!**
- **It is recommended to tap the bottom of the Cartridge Segments containing the reagent cartridges once on the bench before placing them on the rotor. This is to avoid foam and adhering of the liquid on the sealing of the reagent cartridge.**
- Place reagent cartridges on the rotor of the unit. The heating to 37 °C incubation temperature is performed automatically in the unit.

### 6.3 Calibration

#### Traceability:

This method was standardized against the DRG Renin (active) ELISA (EIA-5125).

Each DRG:HYBRID-XL reagent contains a barcode with the specific information for recalibration of the reagent lot. The Master Curve is printed as a 2-D barcode on the outer label of the kit package and on the QC-Datasheet and has to be scanned with the external barcode scanner prior to the first use of the respective kit lot.

#### Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

### 6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRID-XL's system software.

## 7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

It is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analysing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyser. In addition, it is indicated to perform a Recalibration.

In case of further questions please contact your local distributor or DRG directly.

### 7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the QC-Datasheet added to the kit. Note that the expected values and acceptance ranges stated in the QC-Datasheet always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if re-calibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

### 7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

## 8 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and pathological values.

In a study conducted with **EDTA plasma samples** of apparently healthy **adults**, using the DRG:HYBRiD-XL Renin (active) the following values are observed:

| Population       | n  | Mean (pg/mL) | Median (pg/mL) | 5 <sup>th</sup> Percentile (pg/mL) | 95 <sup>th</sup> Percentile (pg/mL) |
|------------------|----|--------------|----------------|------------------------------------|-------------------------------------|
| Supine position  | 55 | 14.14        | 10.20          | 3.58                               | 36.62                               |
| Upright position | 56 | 15.26        | 11.80          | 4.00                               | 45.98                               |

These values are also valid for serum.

In a study conducted with **serum/plasma samples** of apparently healthy **newborns and children**, using the DRG:HYBRiD-XL Renin (active), the following values are observed:

| Age (years) | n  | Range (pg/mL) | Mean (pg/mL) | 1 <sup>st</sup> Percentile (pg/mL) | Median (pg/mL) | 99 <sup>th</sup> Percentile (pg/mL) |
|-------------|----|---------------|--------------|------------------------------------|----------------|-------------------------------------|
| Newborns    | 9  | 16.40 - 6400  | 1281         | 17.07                              | 84.80          | 6224                                |
| 1 - 5       | 12 | 2.41 - 92.90  | 30.07        | 2.41                               | 19.00          | 91.67                               |
| 6 - 18      | 19 | 2.41 - 168    | 38.23        | 2.45                               | 35.50          | 153.06                              |

In a study conducted with apparently normal healthy **adults**, using the DRG:HYBRiD-XL Aldosterone (HYE 5338) and the DRG:HYBRiD-XL Renin (active) (HYE-5373), the following **Aldosterone-Renin Ratios** were determined in plasma:

### Ratio Aldosterone-Renin (pg/mL / pg/mL)

| n   | Mean | Median | 1 <sup>st</sup> Percentile | 99 <sup>th</sup> Percentile |
|-----|------|--------|----------------------------|-----------------------------|
| 110 | 6.39 | 4.07   | 0.48                       | 54.66                       |

These values are also valid for serum.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

## 9 LIMITATIONS OF USE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

## 10 PERFORMANCE CHARACTERISTICS

### 10.1 Assay Dynamic Range

The dynamic range of the assay is defined by the limit of blank and the maximum value of the Master Curve.

Values found below the measuring range are indicated as "< 2.41 pg/mL". Values found above the measuring range are indicated as "> 128 pg/mL". The measuring range of the assay is between 2.41 pg/mL – 128 pg/mL.

### 10.2 Specificity of Antibodies (Cross-Reactivity)

Mean cross reactivity with prorenin was 0.14% (100 – 8000 pg/mL). No influence on the assay was found for the human proteins gamma globulin, serum albumin, pepsin, cathepsin D, and cathepsin B.

### 10.3 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2.

The Limit of Blank (LoB) is 2.41 pg/mL.

The Limit of Detection (LoD) is 4.14 pg/mL.

The Limit of Quantification (LoQ) is 6.18 pg/mL.

### 10.4 Precision Performance

The precision study was designed according to CLSI guideline EP5-A2.

### 10.4.1 Within-Run (Intra Assay)

Samples were analysed 2 times daily in duplicates for 20 days on 3 different instruments. CV (%) is given as mean of 40 duplicate CVs. The within-run variability is shown below:

| Sample | n  | Mean (pg/mL) | CV (%) |
|--------|----|--------------|--------|
| 1      | 40 | 10.16        | 6.22   |
| 2      | 40 | 18.97        | 3.74   |
| 3      | 40 | 63.30        | 3.07   |
| 4      | 40 | 87.19        | 3.02   |

### 10.4.2 Total Precision (Inter Assay)

Samples were analysed 2 times daily in duplicates for 20 days (n=80) on 3 different instruments. The total precision is shown below:

| Sample | n  | Mean (pg/mL) | CV (%) |
|--------|----|--------------|--------|
| 1      | 80 | 10.16        | 14.53  |
| 2      | 80 | 18.97        | 11.86  |
| 3      | 80 | 63.30        | 12.00  |
| 4      | 80 | 87.19        | 11.59  |

### 10.4.3 Inter-Lot

The between-lots variation was determined by 6 measurements of 4 samples with 3 different kit lots.

| Sample | n  | Mean (pg/mL) | CV (%) |
|--------|----|--------------|--------|
| 1      | 18 | 6.93         | 12.75  |
| 2      | 18 | 17.17        | 12.50  |
| 3      | 18 | 67.03        | 5.01   |
| 4      | 18 | 93.00        | 9.78   |

## 10.5 Recovery

Recovery was determined by adding four increasing concentrations of the analyte to four different patient samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was assayed and analyte concentrations of the samples were calculated from the Master Curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

| Sample                | 1     | 2     | 3     | 4     |
|-----------------------|-------|-------|-------|-------|
| Concentration (pg/mL) | 32.50 | 46.95 | 62.25 | 98.70 |
| Average Recovery (%)  | 98.4  | 110.3 | 107.5 | 100.9 |
| Range of Recovery (%) | from  | 93.1  | 106.0 | 98.9  |
|                       | to    | 105.8 | 112.3 | 112.2 |

## 10.6 Linearity

Four samples containing different amounts of analyte were serially diluted with *Sample Diluent*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

| Sample                | 1     | 2     | 3     | 4      |
|-----------------------|-------|-------|-------|--------|
| Concentration (pg/mL) | 41.05 | 79.35 | 84.95 | 129.85 |
| Average Recovery (%)  | 105.5 | 103.7 | 104.6 | 112.7  |
| Range of Recovery (%) | from  | 98.4  | 99.2  | 95.4   |
|                       | to    | 112.5 | 110.4 | 113.0  |

## 10.7 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

The measured renin concentration will increase by approx. 40% after addition of 10 - 100 µg/mL of the renin inhibitor Aliskiren.

No interference was found for:

Captopril (0.5 - 5 µg/mL),  
Furosemid (6 - 60 µg/mL),  
Enalapril (0.05 - 0.5 µg/mL),  
Nicardipin (5 - 50 µg/mL),  
Cathepsin B (0.01 - 0.1 U/mL), and  
Cathepsin D (0.05 - 0.5 U/mL)


## 10.8 High-Dose-Hook Effect

Hook effect was not observed in this test up to a concentration of 2560 pg/mL of renin.

## 11 METHOD COMPARISON

A comparison of DRG:HYBRiD-XL Renin (active) Test HYE-5373 (y) and DRG ELISA EIA-5125 (x) using clinical samples gave the following correlation:

$$\begin{aligned} n &= 151 \\ r &= 0.984 \\ y &= 1.053x - 1.992 \end{aligned}$$

## 12 LEGAL ASPECTS

**Only for countries where the declaration of European Conformity (CE mark) is applicable**

### 12.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws.

This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

### 12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

**For further information please refer to the User Manual of the DRG:HYBRiD-XL, analyser-specific application sheets, product information and package inserts of all necessary components.**

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
**13 REFERENCES / LITERATURE / LITERATUR / BIBLIOGRAFIA / LITERATURA**

1. Imai T, Miyazaki H, Hirose S, et al. Cloning and sequence analysis of cDNA for human renin precursor. *Proc. Natl. Acad. Sci.* (1983) 80, 7405–7409.
2. Reudelhuber TL, Ramla D, Chiu L, et al. Proteolytic processing of human prorenin in renal and non-renal tissues. *Kidney Int.* (1994) 46, 1522–1524.
3. Neves FA, Duncan KG, Baxter JD. Cathepsin B is a prorenin processing enzyme. *Hypertension* (1996) 27, 514–517.
4. Müller DN, Luft FC. Direct renin inhibition with aliskiren in hypertension and target organ damage. *Clin J Am Soc Nephrol.* (2006) 1, 221-8.
5. Toffelmire EB, Slater K, Corvol P, et al. Response of plasma prorenin and active renin to chronic and acute alterations of renin secretion in normal humans. Studies using a direct immunoradiometric assay. *J Clin Invest.* (1989) 83, 679–687.
6. Carey RM, Padia SH. Angiotensin AT2 receptors: control of renal sodium excretion and blood pressure. *Trends Endocrinol Metab.* (2008) 19, 84-7.
7. Koeppen BM, Stanton BA. *Renal Physiology* (4th ed.). Philadelphia, PA. Mosby Physiology Monograph Series, 2007.
8. Navar LG, Inscho EW, Majid DSA, et al. Paracrine regulation of the renal microcirculation. *Physiol. Rev.* (1996) 76, 425–536.
9. Müller MW, Todorov V, Krämer BK, Kurtz A. Angiotensin II inhibits renin gene transcription via the protein kinase C pathway. *Pflugers Arch.* (2002) 444, 499-505.
10. Spät A, Hunyady L. Control of aldosterone secretion: a model for convergence in cellular signaling pathways. *Physiol Rev.* (2004) 84, 489-539.
11. Nguyen G., Delarue F., Burcklé C., et al. Pivotal role of the renin/prorenin receptor in angiotensin II production and cellular responses to renin. *J Clin Invest.* (2002) 109, 1417–1427.
12. Pitarresi TM., Rubattu S, Henrikson R, Sealey JE. Reversible cryoactivation of recombinant human prorenin. *J Biol. Chem.* (1992) 267, 11753-9.
13. Nicari MJ. Specimen processing and renin activity in plasma. *Clin. Chem.* (1992) 38, 598.

## SYMBOLS USED

| Symbol   | English                                    | Deutsch                         | Italiano                            | Español                                    | Français                                 | Polski  |
|--|--|---------------------------------|-------------------------------------|--|--|---|
|  | European Conformity                        | CE-Konformitätskennzeichnung    | Conformità europea                  | Conformidad europea                        | Conforme aux normes européennes          | Zgodność z normami europejskimi                 |
|  | Consult instructions for use*              | Gebrauchsanweisung beachten*    | Consultare le istruzioni per l'uso  | Consulte las Instrucciones                 | Consulter les instructions d'utilisation | Zapoznać się z instrukcją użytkowania           |
|  | <i>In vitro</i> diagnostic medical device* | <i>In-vitro</i> -Diagnostikum*  | Per uso Diagnostica in vitro        | Para uso Diagnóstico in vitro              | Usage Diagnostic in vitro                | Wyrób medyczny do diagnostyki in vitro          |
|  | For research use only                      | Nur für Forschungszwecke        | Solo a scopo di ricerca             | Sólo para uso en investigación             | Seulement dans le cadre de recherches    | Tylko do użytku w badaniach                     |
|  | Catalogue number*                          | Artikelnummer*                  | No. di Cat.                         | Número de catálogo                         | Référence                                | Numer katalogowy                                |
|  | Batch code*                                | Chargencode*                    | Lotto no                            | Número de lote                             | No. de lot                               | Numer LOT                                       |
|  | Contains sufficient for <n> tests *        | Ausreichend für <n> Prüfungen * | Contenuto sufficiente per "n" saggi | Contenido suficiente para <n> ensayos      | Contenu suffisant pour "n" tests         | Zawartość przeznaczona na <n> testów            |
|  | Temperature limit*                         | Temperaturbegrenzung*           | Temperatura di conservazione        | Temperatura de conservación                | Température de conservation              | Przechowywać w temperaturze                     |
|  | Use-by date*                               | Verwendbar bis*                 | Data di scadenza                    | Fecha de caducidad                         | Date limite d'utilisation                | Data przydatności                               |
|  | Manufacturer*                              | Hersteller*                     | Fabbricante                         | Fabricante                                 | Fabricant                                | Producent                                       |
|  | Caution*                                   | Achtung*                        |                                     |  |  | UWAGA   |
| <i>Distributed by</i>  | Distributor                                | Vertreiber                      | Distributore                        | Distribuidor                               | Distributeur                             | Dystrybutor                                     |
| <i>Content</i>   | Content                                    | Inhalt                          | Contenuto                           | Contenido                                  | Contenu                                  | Zawartość                                       |
| <i>Volume/No.</i>  | Volume / No.                               | Volumen/Anzahl                  | Volume/Quantità                     | Volumen/Número                             | Volume/Numéro                            | Objętość / Numer                                |
| <i>Reagent Cartridge</i>   | Reagent Cartridge                          | Reagenzien-Cartridge            | Cartucce di reagente                | Cartucho de reactivo                       | Cartouche de réactif                     | Kartridż z odczynnikami                         |
| <i>Re-Calibrator</i>   | Re-Calibrator                              | Re-Kalibrator                   | Re-Calibratore                      | Re-Calibrador                              | Re-calibrateurs                          | Re-kalibrator                                   |
| <i>Control</i>   | Control                                    | Kontrolle                       | Controllo                           | Control                                    | Contrôle                                 | Kontrola  |
| <i>Enzyme Conjugate</i>  | Enzyme Conjugate                           | Enzymkonjugat                   | Tracciante enzimatico               | Conjugado enzimático                       | Conjugué enzymatique                     | Koniugat enzymatyczny                           |
| <i>Enzyme Complex</i>  | Enzyme Complex                             | Enzymkomplex                    | Complesso enzimatico                | Complex enzimático                         |  | Enzymu kompleks                                 |
| <i>Substrate Solution</i>  | Substrate Solution                         | Substratlösung                  | Soluzione di substrato              | Solución de sustrato                       | Solution substrat                        | Roztwór Substratu                               |
| <i>Sample Diluent</i>  | Sample Diluent                             | Probenverdünnungsmedium         | Diluente dei campioni               | Solución para dilución de la muestra       | Diluant d'échantillon                    | Rozcieńczalnik Próbkki                          |
| <i>Reagent 1</i>   | Reagent 1                                  | Reagenz 1                       | Reagente 1                          | Reactivo 1                                 | Réactif 1                                | Reagent 1                                       |
| <i>Reagent 2</i>   | Reagent 2                                  | Reagenz 2                       | Reagente 2                          | Reactivo 2                                 | Réactif 2                                | Reagent 2                                       |
| <i>Wash Buffer</i>   | Wash Buffer                                | Waschpuffer                     | Tampone di lavaggio                 | Tampón de lavado                           | Tampon de lavage                         | Bufor Płuczący                                  |
| <i>System Solution</i>   | System Solution                            | Systemlösung                    |                                     | Solución de sistema                        | Solution système                         | Roztwór systemowy                               |
| <i>Needle Cleaning Solution</i>  | Needle Cleaning Solution                   | Nadel-Reinigungslösung          | Soluzione lavaggio ago              | Solución de lavado de la aguja             | Solution de nettoyage des aiguilles      | Roztwór Czyszczący Igły                         |
| <i>Denaturation Buffer</i>   | Denaturation Buffer                        | Denaturierungspuffer            | Tampone die denaturazione           | Tampón de denaturalización                 | Tampon de dénaturation                   |   |
| <i>Neutralization Buffer</i>   | Neutralization Buffer                      | Neutralisierungspuffer          | Tampone di neutralizzazione         | Tampón de neutralización                   | Tampon de neutralisation                 |   |
| <i>Assay Buffer</i>  | Assay Buffer                               | Assaypuffer                     | Tampone del test                    | Tampón de ensayo                           | Tampon d'essai                           | Bufor testowy                                   |
|  |  |                                 |                                     |  |  |   |
|  |  |                                 |                                     |  |  |   |
| <i>Secondary Sample Holder</i>   | Secondary Sample Holder                    | Sekundärprobenhalter            | Sostegno secondario dei campioni    | Soporte de muestras para tubos secundarios | Support d'échantillons secondaires       | Stawy na próbki drugorzędowe                    |
| <i>Sample Tubes (Secondary)</i>  | Sample Tubes (Secondary)                   | Sekundärprobenröhrchen          | Tubetti di campioni secondari       | Tubos de muestra secundarios               | Tubes échantillon secondaires            | Probówki drugorzędowe                           |
| <i>Dilution Cartridge</i>  | Dilution Cartridge                         |                                 |                                     |  | Cartouches de dilution                   | Kartridż do rozcieńczeń                         |
| <i>Vial Adapter</i>  | Vial Adapter                               | Fläschchen-Adapter              | l'adattatore dei tubetti            | Adaptador de tubos                         | l'Adaptateur de flacons                  | Nakładki na butelki dla kontroli i kalibratorów |