DRG:HYBRID•XL®

Instructions for Use

Anti-DGP IgA

IVD

REF

HYE-5819

40 tests

CE

DRG 

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de


Distributed by:

DRG 

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

ASSAY PROTOCOL BARCODE (APB)
ASSAYPROTOKOLL-BARCODE / CODICE A BARRE DEL PROTOCOLLO DEL DOSAGGIO / CÓDIGO DE BARRAS DEL PROTOCOLO DEL ENSAYO / CODE-BARRES DU PROTOCOLE DE TEST

(Version 2.61 Software or later / ab Softwareversion 2.61 oder höher / Versione Software 2.61 o superiore / Versión de Software 2.61 o posterior / Version de logiciel 2.61 et supérieure / Wersja oprogramowania 2.61 lub późniejsza)

HYE-5819 – v2.61.1


The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the menu "SCAN NEW LOT".
 Der Barcode muss in dem Software-Menü „NEUES LOT SCANNEN“ eingelesen werden, um das Protokoll in der DRG:HYBRiD-XL-Software zu installieren.

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Please refer to section 6.2.6 "Assay – Scanning a new kit" of the User Manual v2.61 or later.

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4.4	Reagent Preparation:	For wash buffer stability after dilution changed to 1 week at 20 °C to 26 °C (before: 2 weeks).
8	REFERENCE VALUES:	Added: "Since no international reference material is available, the immunoassay is calibrated using external reference samples. The results are given in arbitrary units (U/mL)."
10.2	Diagnostic Specificity and Sensitivity:	Added.
11	METHOD COMPARISON:	Changed to Passing–Bablok regression.
12.4	Reporting of Serious Incident:	Added.

ENGLISH.....	4
DEUTSCH	9
ITALIANO	14
ESPAÑOL.....	19
FRANÇAIS.....	24
REFERENCES / LITERATURE	29
SYMBOLS USED	30

1 INTENDED USE

The DRG:HYBRiD-XL Anti-DGP IgA is an **automated** enzyme immunoassay or the **quantitative** measurement of IgA class antibodies against deamidated gliadin protein (anti-DGP IgA) in human serum or plasma (EDTA, lithium heparin or citrate plasma).

For *in vitro* diagnostic use.

For laboratory professional use, only with the DRG:HYBRiD-XL Analyzer.

1.1 Scientific Validity Report

Celiac disease or gluten sensitive enteropathy is caused by a hypersensitive reaction of genetically predisposed individuals in response to gliadin (1). Gliadin is a protein present in many kinds of cereal grain, e.g. wheat, oats barley and rye. This disease affects the upper small intestine leading to a more or less complete atrophy of the villi of the mucous membrane and a hyperplasia of the crypts of the upper intestine, causing malabsorption problems such as chronic vitamin deficiency (2). Accordingly, patients suffering from celiac disease must maintain a gluten free diet for the rest of their life.

Gliadins are proteins containing high amounts of the amino acids prolin and glutamine. Undigested gliadin peptides are deamidated by the enzyme tissue transglutaminase (tTG) in the intestinal epithelial cells resulting in deamidated gliadin peptides (3). These peptides trigger an immunological response in which antibodies against deamidated gliadin peptide, as well as tTG and gliadin are produced, resulting in a hypersensitive reaction that leads to celiac disease (4).

The serological determination of IgA antibodies directed at deamidated gliadin peptides ensures a highly specific and sensitive diagnosis of celiac disease specially when diagnosing very young children, avoiding the invasive procedures of biopsies (5). It has been well documented that gliadin antibodies level falls very quickly after a gluten free diet has begun and rises immediately after restoring gluten to the diet. Thus, the serological test represents a reliable method to monitor patients, and in particular children and teenagers, for their adherence to the gluten-free diet (6).

2 PRINCIPLE OF THE TEST

The DRG:HYBRiD-XL Anti-DGP IgA Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **indirect immunoassay principle**.

The wells (ACW) of the reagent cartridges are coated with deamidated gliadin peptide (DGP).

An aliquot of automatically pre-diluted patient sample containing endogenous Anti-DGP IgA antibodies is incubated in the coated well. After incubation, unbound sample components are washed off. Thereafter, enzyme conjugate, which is an anti-human IgA antibody conjugated with horseradish peroxidase is incubated in the coated well. After incubation the unbound conjugate is washed off.

Having added the substrate solution, the intensity of color developed is proportional to the concentration of anti-DGP IgA in the patient sample.

The concentrations of unknown samples are determined from the master curve.

3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only.
For laboratory professional use only.
- This kit can only be used in combination with the DRG:HYBRiD-XL Analyzer.
- Before starting the assay, read the instructions for use completely and carefully. Use the valid version of the instructions for use provided with the kit. Be sure that everything is understood.
- Do not remove, exchange, discard, or damage any of the barcode labels provided with each kit and its components.
All barcodes build an integral system for the kit lot.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact DRG.

- Optimal test results are only obtained when using calibrated pipettes.
- Microbial contamination of reagents or samples may give false results.

General precautions

- Follow laboratory quality assurance and laboratory safety guidelines.
- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.

Biohazard information

- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. However, no known test method can offer total assurance that no infectious agent is present.
- The device contains material of animal origin, which is certified apparently free of infectious or contagious diseases and injurious parasites.
- Bovine components originate from countries where BSE (Bovine spongiform encephalopathy) has not been reported.
- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation.
Waste must be discarded according to local rules and regulations.

Information to chemical hazards and hazard classification

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Substrate Solution contains an ingredient in non-declarable concentrations which causes serious eye irritation. In case of possible contact with eyes, rinse immediately carefully and thoroughly with eye wash or water. After contact with skin, wash with plenty of water. Take-off contaminated clothing and wash it before reuse
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national safety guidelines or regulations.
- This product does not contain substances which have carcinogenic, mutagenic, or toxic for reproduction (CMR) properties.

All reagents of this test kit do NOT contain hazardous substances in concentrations to be declared, a classification and labelling is not required.

A safety data sheet is not required for this product.

For more information, please refer to the Safety Information, which is available upon request directly from DRG.

4 MATERIALS

4.1 Materials provided with the kit

Symbol	Quantity/ Preparation	Description
Reagent Cartridges	2 x 20 pieces in a resealable plastic bag with desiccant Ready to use	Reagent Cartridge containing the following: – Antibody Coated Well (ACW) coated with DGP. – Enzyme Conjugate * , 170 µL, Anti-human IgA antibody conjugated with horseradish peroxidase – Substrate Solution , 270 µL Contains 3,3',5,5'-tetramethylbenzidine. – Sample Diluent * , 300 µL
Additional Cartridges	2x 20 pieces Ready to use	Additional Cartridges Red; empty cartridges for automatically predilution of samples.
Re-Calibrator R1 & R2	2 x 1.0 mL Ready to use	Recalibrators * For recalibration of the quantitative DRG:HYBRID-XL Anti-DGP IgA test. <u>Concentrations are lot-specific.</u>
Control C1 & C2	2 x 1.0 mL Ready to use	Controls * For control values and ranges please refer to the vial label or to the CoA.
	1 x	Instructions for Use
	1 x	Certificate of Analysis (CoA)
* Contain non-mercury preservative.		

4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:
Resistivity at 25 °C [MΩ·cm]: > 10
Conductivity at 25 °C [µS·cm⁻¹]: < 0.1
Total Organic Carbon/p.p.b.[µg/L] : < 50
Colloids [µg/mL]: <0.05
- REF HYB-5252 DRG:HYBRID-XL Analyzer
- REF HYI-5392: *System Solution 5L*, 5000 mL;
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:
Resistivity at 25 °C [MΩ·cm]: > 1
Conductivity at 25 °C [µS·cm⁻¹]: < 1
Total Organic Carbon/p.p.b.[µg/L] : < 200
Colloids [µg/mL]: <0.1)
- REF HYI-6234: *Wash Buffer*, 10x concentrate, 100 mL
- REF HYI-5395: *Needle Cleaning Solution*, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- REF HYI-5387: *Cuvettes*, (2 x 360 pieces)

For use of the *Secondary Sample Holder* for secondary tubes the following tubes are required:

- REF HYI-5391: *Sample Tubes (Secondary)*, 2500 pieces

4.3 Storage and Stability of the Kit

All kit components must be handled and stored as described below.

Unopened and unused components or kits

	Storage Temperature	Stability
– Unopened kits – Unopened reagents – Cartridges in the supplied <u>unopened</u> plastic bag – Unused reagent cartridges	2 °C to 8 °C	Until the expiration date printed on the label. Do not use reagents beyond this date!

Unused cartridges must be stored in the blister pack in the sealed foil pouch including the desiccant. Once the foil bag has been opened, care must be taken to close it tightly again.

Opened components

	Storage Temperature	Stability
Opened – Re-Calibrators – Controls	2 °C to 8 °C	8 weeks

Immediately after end of each run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped, and stored at 2 °C to 8 °C.

Pierced or open cartridges must be disposed of immediately.

4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 26 °C). Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

On-board stability:	24 h
---------------------	------

4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 26 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

Wash Buffer (10x) (HYI-6234; not included in the kit)

For Wash Buffer (1x) dilute 100 mL of Wash Buffer (10x) with 900 mL ultra-pure water to a final volume of 1000 mL.

stability after dilution:	20 °C - 26 °C	4 weeks
---------------------------	---------------	---------

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

Human serum or plasma (EDTA plasma, lithium heparin plasma or citrate plasma)

A minimum volume of 68 µL of sample is needed for one determination. This includes 8 µL sample and 60 µL dead volume.

A minimum of 160 µL of Control and Re-Calibrator (ready-to-use) is needed for one determination. This includes each 100 µL Control and Re-Calibrator and 60 µL dead volume.

Attention:

- This test was not verified with blood collection tubes of all available manufacturers.
- Sample collection systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- In general, it should be avoided to use hemolytic, icteric, or lipemic samples. For further information refer to chapter "Interfering Substances".
- Samples containing precipitates have to be centrifuged prior to the test run.
- Do not use heat inactivated samples.
- Samples or external controls containing sodium azide should not be used in the assay.

5.1 Sample Collection

Serum: Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Samples of patients receiving anticoagulant therapy may require increased clotting time.

Plasma: Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.2 Sample Storage

Samples must be stored tightly capped prior to performing the assay. If stored frozen, freeze only once. Thawed samples must be inverted several times prior to testing.

Stability	at 2 °C to 8 °C	7 days
	at -20 °C (in aliquots)	up to 12 months

5.3 Sample Preparation

Samples can be assayed without additional preparation.

5.3.1 Fully Automated Predilution of Samples

Each sample will be **prediluted fully automated 1:100** using the additional cartridge (red) included in the kit (1 cartridge per 1:100 dilution). Sample Diluent for automated predilution is included in the Reagent Cartridge.

6 ASSAY PROCEDURE

6.1 Procedural Notes

- All reagents, such as controls and re-calibrators, and samples, must be allowed to come to room temperature (20 °C to 26 °C) before use. All reagents must be mixed without foaming. Reagent Cartridges can be used directly from the cold storage (refrigerator).
- Avoid bubbles/foam in the calibrator and control vials, especially when using them directly in the *Primary Sample Holder* (HY1-5436). Bubbles/foam can lead to wrong liquid level detection and hence false results.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- To avoid possible evaporation effects of the recalibrators, controls and samples placed on the instrument, the chapter "On-Board Stability" must be considered.
- For primary samples tubes the *Primary Sample Holder* with a capacity for 12 samples can be used. For low volume samples the *Secondary Sample Holder* (HY1-5437) can be used. It has the capacity for a maximum of 20 samples for secondary tubes. The samples must be transferred manually into the secondary tubes (*Sample Tubes (Secondary)*).

6.2 Test Procedure

Total assay time:	120 min
--------------------------	----------------

Detailed information about operating the DRG:HYBRiD-XL analyzer must be taken from the system's User Manual.

Test specific information required for automatic processing is contained in distinct barcodes of the kit.

- **Assay Protocol Barcode (APB):**

2D Data Matrix code containing the specific assay protocol/method. The APB is found on page 2 of this IFU.

- **Lot barcode:**

2D Data Matrix code containing lot specific information including the master curve (factory calibration). The lot barcode is found on the label of the kit package and CoA.

- **Re-calibrator and control vial barcodes:**

2D Data Matrix code to identify the re-calibrators and controls. This code is found on the vial labels and CoA.

- **Cartridge barcode:**

Reagent cartridges are sealed with a 1D barcode for cartridge identification. Take care not to damage this barcode when inserting the cartridges.

- If it is the first use of this test on the analyzer the **Assay Protocol Barcode (APB)** on page 2 of this instructions for use must be loaded with the external barcode scanner and stored. The APB includes all information about the test procedure itself (specific protocol data, including volumes, incubation times, units etc.). Once installed it will remain until a new APB will become available by the manufacturer.
- If the lot is used the first time or a new lot of this test is available, the lot barcode on the outer label of the kit package or on the CoA must be loaded with the external barcode scanner and stored.
- Program the number of required samples and continue with the test assignment.
- Insert the number of required cartridges and cuvettes with care in the segment. **Make sure cartridges and cuvettes are inserted correctly.**
- Load the segments by placing them on the reagent rotor. Make sure the inserted segment is seated correctly (not tilted!).
- Continue to start the run.

6.3 Recalibration

Detailed information how to perform the recalibration is found in the system's User Manual.

Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

Always both recalibrators (R1 & R2) are required for recalibration.

Always two internal controls are required to validate the recalibration.

Recalibrators and controls are run together in the same run. In addition, samples can be included.

In the "Load-Samples" menu scan the respective barcodes for re-calibrators and controls found on the vials or CoA.

Test assignment for re-calibrators and controls is done automatically.

Continue with loading the cartridges and cuvettes required.

Start the run.

6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRiD-XL's system software.

For the result calculation a four-parameter logistic (4PL) curve fit (4PL Rodbard) is used.

If in an initial assay, a sample is found to contain concentrations above the measuring range, the samples can be diluted manually with *Sample Diluent* and measured again.

For the calculation of the concentrations this dilution factor must be considered.

7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results.

If available, it is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyzer. In addition, it is indicated to perform a Recalibration.

After checking the above-mentioned items without finding any error contact your distributor or DRG directly.

7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the Certificate of Analysis added to the kit. Note that the expected values and acceptance ranges stated in the Certificate of Analysis always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if recalibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

8 REFERENCE VALUES

Since no international reference material is available, the immunoassay is calibrated using external reference samples. The results are given in arbitrary units (U/mL).

It is strongly recommended that each laboratory should determine its own reference values.

In a study conducted with apparently healthy individuals and additional with patients with coeliac disease, and a positive anti-DGP IgA result (by EIA-5909) using the DRG:HYBRiD-XL Anti-DGP IgA the following data were observed:

Population	Apparently Healthy (U/mL)	Coeliac Disease (U/mL)
n	77	58
Mean (U/mL)	4.89	66.94
Median (U/mL)	4.52	35.93
2.5 th - 97.5 th Percentile (U/mL)	1.45 – 10.51	10.29 – 291.08
Range (min. - max.) (U/mL)	1.20 – 15.20	8.70 – 328.00

Values above or below the reference range should be considered as suspicious and require additional testing.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated with other clinical observations and diagnostic tests.

8.1 Cut-off Value

The Cut-off Value of the DRG:HYBRiD-XL Anti-DGP IgA was determined by measuring the values of apparently healthy individuals and patients with coeliac disease and a positive anti-DGP IgA result (obtained by the reference method DRG MGP IgA ELISA (Modified Gliadin Peptide) (EIA-5909))

The Cut-off value was defined as 10.2 U/mL.

Using a Cut-off value of 10.2 U/mL, the specificity was calculated to 97.4% and sensitivity to 96.6%.

8.1.1 Interpretation of results

Normal (negative)	< 8.5 U/mL
Cut-off	10.2 U/mL
Grey zone	8.5 U/mL – 15.2 U/mL
Positive	> 15.2 U/mL

9 LIMITATIONS OF THE PROCEDURE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and in compliance with the laboratory quality assurance guidelines.

Any improper handling of samples or modification of this test might influence the results.

10 PERFORMANCE CHARACTERISTICS

10.1 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2.

Limit of Blank (LoB)	0.570 U/mL
Limit of Detection (LoD)	0.906 U/mL
Limit of Quantification (LoQ)	1.683 U/mL

Measuring range	0.906 U/mL – 100 U/mL
Linear range	0.903 U/mL – 100 U/mL

	Indicated as:
Values below the measuring range	"< 0.906 U/mL"
Values above the measuring range	"> 100 U/mL"

10.2 Diagnostic Specificity and Sensitivity

Diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

Diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

Diagnostic Sensitivity	96.6 %
Diagnostic Specificity	97.4 %
Positive Predictive Value	96.6 %
Negative Predictive Value	97.4 %
Likelihood Ratio	37.17

10.3 Precision Performance

The precision study was designed on the basis of to CLSI guideline EP5-A2.

10.3.1 Within-run precision

The Intra-assay precision was determined with 4 patient samples covering the complete measuring range in 40 runs with 2 replicates per run. CV was calculated as mean CV of 40 runs.

The within-run variability is shown below:

Sample	n	Mean (U/mL)	CV (%)
1	40	3.24	4.7
2	40	11.45	4.1
3	40	19.00	4.0
4	40	28.09	3.7

10.3.2 Between-run precision

The Between-run precision was determined with 4 patient samples covering the measuring range in 40 independent runs on 20 days in duplicate determinations. CV was calculated from 80 determinations.

The total precision is shown below:

Sample	n	Mean (U/mL)	CV (%)
1	80	3.24	13.9
2	80	11.45	7.6
3	80	19.01	8.2
4	80	28.09	7.5

10.3.3 Between-lot precision

The between-lot variation was determined by 6 measurements of 4 samples with 3 different kit lots.

Sample	n	Mean (U/mL)	CV (%)
1	18	2.95	14.1
2	18	11.69	2.9
3	18	14.74	7.2
4	18	26.85	7.7

10.4 Recovery

Recovery was determined by adding increasing amounts of the analyte to different patient samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was assayed, and analyte concentrations of the samples were calculated from the master curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

Sample	1	2	3	4	
Concentration (U/mL)	8.75	10.70	13.55	34.50	
Average Recovery (%)	105.8	100.4	97.5	107.1	
Range of Recovery (%)	from	101.4	96.7	97.0	98.5
	to	113.6	106.6	98.0	112.4

10.5 Linearity

Samples containing different amounts of analyte were serially diluted with *Sample Diluent*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

Sample	1	2	3	4	
Concentration (U/mL)	18.05	18.95	14.10	9.46	
Average Recovery (%)	115.4	116.7	113.5	108.3	
Range of Recovery (%)	from	108.3	113.1	106.1	100.6
	to	119.3	119.9	117.9	115.1

10.6 Interfering Substances

Hemoglobin (up to 4 mg/mL), bilirubin (up to 0.5 mg/mL) and triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

10.7 High-Dose Hook Effect

"High-Dose Hook Effect" is not detected up to 2000 U/mL of anti-DGP-IgA.

11 METHOD COMPARISON

A comparison of DRG:HYBRiD-XL Anti-DGP IgA Test HYE-5819 (y) and the DRG MGP IgA ELISA, EIA-5909 (x) using clinical samples gave the following correlation:

$$\begin{aligned} n &= 111 \\ r &= 0.975 \\ y &= 0.991x + 1.066 \end{aligned}$$

12 LEGAL ASPECTS**12.1 Reliability of Results**

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the laboratory quality assurance guidelines and applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact DRG.

12.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

12.4 Reporting of Serious Incident

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

For further information please refer to the User Manual of the DRG:HYBRiD-XL, analyzer-specific application sheets, product information and package inserts of all necessary components.

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