
DRG:HYBRID•XL[®]

Instructions for Use

Linezolid

IVD

REF

HYE-6163

40 tests

CE

DRG 

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de


Distributed by:

DRG 

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

ASSAY PROTOCOL BARCODE (APB)
ASSAYPROTOKOLL-BARCODE / CODICE A BARRE DEL PROTOCOLLO DEL DOSAGGIO / CÓDIGO DE BARRAS DEL PROTOCOLO DEL ENSAYO / CODE-BARRES DU PROTOCOLE DE TEST

(Software Version 2.70 or later / ab Softwareversion 2.70 oder höher / Versione software 2.70 o successiva / Versión 2.70 del software o posterior / Version de logiciel 2.70 ou ultérieure / Wersja oprogramowania 2.70 lub późniejsza)

459 - LINE - HYE-6163 - 02.70.02

HYE-6163 – v2.70.2


The barcode must be used to install the assay protocol into the DRG:HYBRiD-XL software via the menu "SCAN NEW LOT".
 Der Barcode muss in dem Software-Menü „NEUES LOT SCANNEN“ eingelesen werden, um das Protokoll in der DRG:HYBRiD-XL-Software zu installieren.

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Please refer to section 6.2.6 "Assay – Scanning a new kit" of the User Manual v2.61 or later.

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ASSAY PROTOCOL BARCODE (APB) changed from v2.70.1 to v2.70.2.

4.4 Reagent Preparation:	Stability for diluted wash buffer changed to 1 week at 20 °C to 26 °C (old: 2 weeks).
9.1 Sensitivity:	Measuring range is updated to 52.0 ng/mL – 3500 ng/mL (old; 13.340 ng/mL – 3000 ng/mL).
9.8 High-Dose Hook Effect:	updated

ENGLISH.....	4
DEUTSCH.....	9
ITALIANO.....	14
ESPAÑOL.....	18
FRANÇAIS.....	23
REFERENCES / LITERATURE	28
SYMBOLS USED	29

1 INTENDED USE

The DRG:HYBRiD-XL Linezolid is an **automated** enzyme immunoassay or the **quantitative** measurement of Linezolid in human serum or plasma (EDTA, lithium heparin or citrate plasma).

For *in vitro* diagnostic use.

For laboratory professional use, only with the DRG:HYBRiD-XL Analyzer.

The device is **intended to be used** for therapeutic drug monitoring (TDM) of patients under Linezolid treatment. TDM will ensure tight control of Linezolid levels during treatment.

The DRG:HYBRiD-XL Linezolid can be used as an aid for immediate dose adjustment to avoid toxic effects by Linezolid overdosing and ineffective bacterial elimination by Linezolid underdosing.

The device is **not intended** for estimating anti-bacterial efficacy of Linezolid treatment.

1.1 Scientific Validity Report

Linezolid is an antibiotic from the oxazolidinone substance class. This group is part of the beta-lactam antibiotics.

Linezolid is prescribed for pneumonia but also for skin and pelvic infections [1,2,11,12]. Linezolid is also effective against methicillin-resistant *Staphylococcus aureus* (MRSA) [3].

Linezolid inhibits bacterial protein synthesis by a unique mechanism of action. The antibiotic targets a binding site on the bacterial ribosome (23S of the 50S subunit, preventing the formation of a functional 70S initiation complex, which is an essential part of the translation process.[4]

Linezolid is rapidly and extensively absorbed after oral administration, such that peak plasma concentrations are reached within two hours.[5]

The key element for successful antibiotic therapy is to adjust the dosage of the antibiotic patient-specifically in such a way that its concentration in the blood remains in a microbial effective concentration range over the entire duration of treatment (Therapeutic Drug Monitoring, TDM) [6,7]. Especially intensive care patients are at risk of insufficient dose of the drug due to restricted bodily functions [8]. The underdosing of antibiotics as well as the general environmental pollution promote the increase of antibiotic resistance [9,10].

In summary, the precise and accurate measurement of Linezolid by enzyme immunoassay can allow the determination of antibiotic concentration in the blood for immediate dose adjustment.

2 PRINCIPLE OF THE TEST

The DRG:HYBRiD-XL Linezolid Kit is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **principle of competitive binding**.

The wells (ACW) of the reagent cartridges are coated with Linezolid coupled to a carrier protein.

During the first incubation, Linezolid in the added sample competes with the bound Linezolid on the coated well for binding to the horseradish peroxidase conjugated anti-Linezolid-antibody.

After incubation unbound substances are washed off.

Having added the substrate solution, the intensity of color developed is inversely proportional to the concentration of Linezolid in the sample.

The concentrations of unknown samples are determined from the master curve.

3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only.
For laboratory professional use only.
- This kit can only be used in combination with the DRG:HYBRiD-XL Analyzer.
- Before starting the assay, read the instructions for use completely and carefully. Use the valid version of the instructions for use provided with the kit. Be sure that everything is understood.
- **Do not remove, exchange, discard, or damage any of the barcode labels provided with each kit and its components.**
All barcodes build an integral system for the kit lot.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange reagent cartridges of different kits even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the wells in the reagent cartridges may differ slightly.
- Do not use reagents beyond expiry date as shown on the kit labels.

- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact DRG.
- Optimal test results are only obtained when using calibrated pipettes.
- Microbial contamination of reagents or samples may give false results.

General precautions

- Follow laboratory quality assurance and laboratory safety guidelines.
- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.

Biohazard information

- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. However, no known test method can offer total assurance that no infectious agent is present.
- The device contains material of animal origin, which is certified apparently free of infectious or contagious diseases and injurious parasites.
- Bovine components originate from countries where BSE (Bovine spongiform encephalopathy) has not been reported.
- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation. Waste must be discarded according to local rules and regulations.

Information to chemical hazards and hazard classification

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Substrate Solution contains an ingredient in non-declarable concentrations which causes serious eye irritation. In case of possible contact with eyes, rinse immediately carefully and thoroughly with eye wash or water. After contact with skin, wash with plenty of water. Take-off contaminated clothing and wash it before reuse
- Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national safety guidelines or regulations.
- This product does not contain substances which have carcinogenic, mutagenic, or toxic for reproduction (CMR) properties.

All reagents of this test kit do NOT contain hazardous substances in concentrations to be declared, a classification and labelling is not required.

A safety data sheet is not required for this product.

For more information, please refer to the Safety Information, which is available upon request directly from DRG.

4 MATERIALS

4.1 Materials provided with the kit

Symbol	Quantity/ Preparation	Description
Reagent Cartridges	2 x 20 pieces in a resealable plastic bag with desiccant Ready to use	Reagent Cartridge containing the following: – Coated Well (ACW) coated with Linezolid coupled to a carrier protein. – Sample Diluent *, 230 µL Enzyme Conjugate *, 170 µL, Anti-Linezolid antibody conjugated with horseradish peroxidase. – Substrate Solution , 270 µL Contains 3,3',5,5'-tetramethylbenzidine.
Dilution Cartridges	2 x 20 pieces Ready to use	Dilution Cartridge Yellow, for fully automated sample dilution in this test.
Re-Calibrator R1 & R2	2 x 1 mL Ready to use	Recalibrators * For recalibration of the quantitative DRG:HYBRiD-XL Linezolid test. Concentrations are lot-specific. <i>Standardized against the following reference material: Eureka Kit, code LC79016</i>
Control C1 & C2	2 x 1 mL Ready to use	Controls * For control values and ranges please refer to the vial label or to the CoA.
	1 x	Instructions for Use
	1 x	Certificate of Analysis (CoA)
* Contain non-mercury preservative.		

4.2 Materials required but not provided

- General needed laboratory equipment
- Ultra-pure water
DRG recommends to use Clinical Laboratory Reagent Water (CLRW) according to CLSI guideline 3C-A4 with the following specifications:
Resistivity at 25 °C [MΩ·cm]: > 10
Conductivity at 25 °C [µS·cm⁻¹]: < 0.1
Total Organic Carbon/p.p.b.[µg/L]: < 50
Colloids [µg/mL]: <0.05
- **REF** HYB-5252 DRG:HYBRiD-XL Analyzer
- **REF** HYI-5392: *System Solution 5L*, 5000 mL;
(Instrument Feed Water according to CLSI guideline 3C-A4 with the following specification can also be used:
Resistivity at 25 °C [MΩ·cm]: > 1
Conductivity at 25 °C [µS·cm⁻¹]: < 1
Total Organic Carbon/p.p.b.[µg/L]: < 200
Colloids [µg/mL]: <0.1)
- **REF** HYI-6234: *Wash Buffer*, 10x concentrate, 100 mL
- **REF** HYI-5395: *Needle Cleaning Solution*, 30 mL. Cleaning solution for the pipetting needle (daily and weekly maintenance, see also user manual)
- **REF** HYI-5387: *Cuvettes*, (2 x 360 pieces)

For use of the *Secondary Sample Holder* for secondary tubes the following tubes are required:

- **REF** HYI-5391: *Sample Tubes (Secondary)*, 2500 pieces

4.3 Storage and Stability of the Kit

All kit components must be handled and stored as described below.

Unopened and unused components or kits

	Storage Temperature	Stability
– Unopened kits – Unopened reagents – Cartridges in the supplied <u>unopened</u> plastic bag – Unused reagent cartridges	2 °C to 8 °C	Until the expiration date printed on the label. Do not use reagents beyond this date!

Unused cartridges must be stored in the blister pack in the sealed foil pouch including the desiccant. Once the foil bag has been opened, care must be taken to close it tightly again.

Opened components

	Storage Temperature	Stability
Opened – Re-Calibrators – Controls	2 °C to 8 °C	8 weeks

Immediately after end of each run the Re-Calibrator and Control vials have to be removed from the instrument, tightly capped, and stored at 2 °C to 8 °C.

Pierced or open cartridges must be disposed of immediately.

4.3.1 On-board Stability

For Re-Calibrators and Controls the on-board stability has been evaluated under controlled laboratory conditions at room temperature (20 °C to 26 °C).

Due to the differences in laboratory environmental conditions and reagent volumes, the on-board stability may deviate from the declared value.

On-board stability:	32 h
---------------------	------

4.4 Reagent Preparation

Bring all reagents, such as controls and re-calibrators, to room temperature (20 °C to 26 °C) prior to use. Reagent Cartridges can be used directly after storage in the refrigerator.

Wash Buffer (10x) (HYI-6234; not included in the kit)

For Wash Buffer (1x) dilute 100 mL of Wash Buffer (10x) with 900 mL ultra-pure water to a final volume of 1000 mL.

stability after dilution:	20 °C - 26 °C	4 weeks
---------------------------	---------------	---------

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

Human serum or plasma (EDTA plasma, lithium heparin plasma or citrate plasma)

For **patients samples** a minimum volume of **76 µL** is needed for one determination. This includes **16 µL** sample and 60 µL dead volume.

For **Re-Calibrators and Controls**, a minimum volume of **80 µL** is needed for one determination. This includes **20 µL** sample and 60 µL dead volume.

Attention:

- This test was not verified with blood collection tubes of all available manufacturers.

- Sample collection systems of some manufacturers may contain different materials which in isolated cases could affect the test results.
- If primary tubes for sample collection are used, please follow the instructions of the manufacturer.
- In general, it should be avoided to use hemolytic, icteric, or lipemic samples. For further information refer to chapter "Interfering Substances".
- Samples containing precipitates have to be centrifuged prior to the test run.
- Do not use heat inactivated samples.
- Samples or external controls containing sodium azide should not be used in the assay.

5.1 Sample Collection

Serum: Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Samples of patients receiving anticoagulant therapy may require increased clotting time.

Plasma: Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.2 Sample Storage

Samples must be stored tightly capped prior to performing the assay. If stored frozen, freeze only once. Thawed samples must be inverted several times prior to testing.

Stability	at 2 °C to 8 °C	7 days
	at -20 °C (in aliquots)	up to 12 months

5.3 Sample Preparation

Samples can be assayed without additional preparation.

5.3.1 Fully Automated Sample Dilution

A fully automatic sample dilution (dilution 1:10) is included in this test procedure.

The required additional dilution cartridges (1 cartridge per sample) are included in this test kit.

Sample diluent for automated sample dilution is included in the Reagent Cartridge.

6 ASSAY PROCEDURE

6.1 Procedural Notes

- All reagents, such as controls and re-calibrators, and samples, must be allowed to come to room temperature (20 °C to 26 °C) before use. All reagents must be mixed without foaming. Reagent Cartridges can be used directly from the cold storage (refrigerator).
- Avoid bubbles/foam in the calibrator and control vials, especially when using them directly in the *Primary Sample Holder* (HY1-5436). Bubbles/foam can lead to wrong liquid level detection and hence false results.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- To avoid possible evaporation effects of the recalibrators, controls and samples placed on the instrument, the chapter "On-Board Stability" must be considered.
- For primary samples tubes the *Primary Sample Holder* with a capacity for 12 samples can be used. For low volume samples the *Secondary Sample Holder* (HY1-5437) can be used. It has the capacity for a maximum of 20 samples for secondary tubes. The samples must be transferred manually into the secondary tubes (*Sample Tubes (Secondary)*).

6.2 Test Procedure

Total assay time: 60 min

Detailed information about operating the DRG:HYBRiD-XL analyzer must be taken from the system's User Manual.

Test specific information required for automatic processing is contained in distinct barcodes of the kit:

- **Assay Protocol Barcode (APB):**
2D Data Matrix code containing the specific assay protocol/method. The APB is found on page 2 of this IFU.

- **Lot barcode:**
2D Data Matrix code containing lot specific information including the master curve (factory calibration). The lot barcode is found on the label of the kit package and CoA.
- **Re-calibrator and control vial barcodes:**
2D Data Matrix code to identify the re-calibrators and controls. This code is found on the vial labels and CoA.
- **Cartridge barcode:**
Reagent cartridges are sealed with a 1D barcode for cartridge identification. Take care not to damage this barcode when inserting the cartridges.
- If it is the first use of this test on the analyzer the **Assay Protocol Barcode (APB)** on page 2 of this instructions for use must be loaded with the external barcode scanner and stored. The APB includes all information about the test procedure itself (specific protocol data, including volumes, incubation times, units etc.). Once installed it will remain until a new APB will become available by the manufacturer.
- If the **lot** is used the first time or a new lot of this test is available, the lot barcode on the outer label of the kit package or on the CoA must be loaded with the external barcode scanner and stored.
- Program the number of required samples and continue with the test assignment.
- Insert the number of required cartridges and cuvettes with care in the segment. **Make sure cartridges and cuvettes are inserted correctly.**
- Load the segments by placing them on the reagent rotor. Make sure the inserted segment is seated correctly (not tilted!).
- Continue to start the run.

6.3 Recalibration

Detailed information how to perform the recalibration is found in the system's User Manual.

Recalibration is recommended:

- if a new kit lot is used. Each new lot should be verified by running the kit internal re-calibrators and controls before routine use.
- if one or both assay controls are found outside the specified range.
- after 4 weeks of use of the same reagent kit on the unit.

Always both recalibrators (R1 & R2) are required for recalibration.

Always two internal controls are required to validate the recalibration.

Recalibrators and controls are run together in the same run. In addition, samples can be included.

In the "Load-Samples" menu scan the respective barcodes for re-calibrators and controls found on the vials or CoA.

Test assignment for re-calibrators and controls is done automatically.

Continue with loading the cartridges and cuvettes required.

Start the run.

6.4 Calculation of Results

The analyte concentrations are calculated automatically by the DRG:HYBRiD-XL's system software.

For the result calculation a four-parameter logistic (4PL) curve fit (4PL Rodbard) is used.

7 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results.

If available, it is also recommended to participate in national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid. In this case, please check the following: expiration dates and storage conditions of reagents, operational reliability of the analyzer. In addition, it is indicated to perform a Recalibration.

After checking the above-mentioned items without finding any error contact your distributor or DRG directly.

7.1 Internal Controls

For Quality Control it is necessary to use the two internal controls provided with each kit.

Acceptance ranges for both internal controls (*Control 1 & 2*) were established by the manufacturer and are summarized in the Certificate of Analysis added to the kit. Note that the expected values and acceptance ranges stated in the Certificate of Analysis always refer to the current kit lot.

Internal controls should be run in single determination:

- on a routine basis (e.g. once per 24 h)
- if recalibration is required (if one or both internal controls are out of range)
- if a new kit lot is used (in order to avoid any negative impact on the kit performance by improper transport or to detect improper storage during transport).

7.2 External Controls

Use controls at both normal and pathological levels.

The control intervals and control ranges for external controls should be adapted to the individual requirements of each laboratory. All results must be within the defined limits.

Each laboratory should establish corrective measures to be taken if values of external controls are not found in the acceptance range.

8 LIMITATIONS OF THE PROCEDURE

Reliable and reproducible results will be obtained, when the assay procedure is performed with a complete understanding of the instructions for use and in compliance with the laboratory quality assurance guidelines.

Any improper handling of samples or modification of this test might influence the results.

9 PERFORMANCE CHARACTERISTICS

9.1 Sensitivity

The sensitivity study was designed according to CLSI guideline EP17-A2.

Limit of Blank (LoB)	8.369 ng/mL
Limit of Detection (LoD)	13.340 ng/mL
Limit of Quantification (LoQ)	16.980 ng/mL

Linear measuring range	52.0 ng/mL – 3500 ng/mL
Final linear measuring range, considering the automated pre-dilution of 1:10	0.520 µg/mL – 35 µg/mL

	Indicated as:
Values below the measuring range	"< 52.0 ng/mL"
Values above the measuring range	"> 3500 ng/mL"

9.2 Specificity of Antibodies (Cross-Reactivity)

The following substances were tested for cross-reactivity of the assay:

Substance	Conc. Range of Added Substance	Mean Cross-Reactivity
	ng/mL	%
Amoxicillin	100 - 10000	0.0
Cefepime HCl/Di-HCl	100 - 10000	0.0
Cefixime Trihydrat	10 - 1000	0.0
Cefotaxine	10 - 1000	0.0
Cefpodoxine	10 - 1000	0.0
Ceftazidime	100 - 10000	0.0
Ceftriaxone	10 - 1000	0.1
Flucoxillin	10 - 1000	0.0
Levofloxacin	100 - 10000	0.0
Meropenem	100 - 10000	0.0
Moxifloxacin	100 - 10000	0.0
Penicillin Na-Salt	10 - 1000	0.0
Piperacillin	10 - 1000	0.0
Tazobactam	1000 - 100000	0.0
Tedizolid	10 - 1000	0.0

Substance	Conc. Range of Added Substance	Mean Cross-Reactivity
	ng/mL	%
Ciprofloxacin	100 - 10000	0.0
Daptomycin	100 - 10000	0.0

9.3 Precision Performance

The precision study was designed on the basis of to CLSI guideline EP5-A3.

9.3.1 Within-run precision

Calculated from the mean value of 5 replicates on 5 days with 1 device.

Sample	n	Mean (ng/mL)	CV (%)
1	5	681.2	6.3
2	5	8924.8	4.8
3	5	8221.2	2.4
4	5	10152.0	4.3

9.3.2 Between-run precision

Calculated by 5 replicates each on 5 days with 1 device.

Sample	n	Mean (ng/mL)	CV (%)
1	25	705.2	9.0
2	25	6259.3	7.4
3	25	8435.2	8.2
4	25	10168.1	9.0

9.3.3 Between-device precision

Calculated by 5 replicates each on 5 days with 3 devices.

Sample	n	Mean (ng/mL)	CV (%)
1	75	705.2	10.4
2	75	6259.3	8.5
3	75	8435.2	10.2
4	75	10168.1	10.0

9.3.4 Between-lot precision

Determined with 3 different kit lots.

Sample	n	Mean (ng/mL)	CV (%)
1	18	348.6	9.8
2	18	790.0	10.2
3	18	929.4	7.1
4	18	2298.7	4.9

9.4 Recovery

Recovery was determined by adding increasing amounts of the analyte to different patient samples containing different amounts of endogenous analyte. Each sample (non-spiked and spiked) was assayed, and analyte concentrations of the samples were calculated from the master curve. The percentage recoveries were determined by comparing expected and measured values of the samples.

Sample	1	2	3	4	
Concentration (ng/mL)	892.0	1205.0	1280.0	2180.0	
Average Recovery (%)	99.5	94.0	93.0	90.7	
Range of Recovery (%)	from	88.0	92.8	88.3	87.2
	to	108.3	96.4	102.8	94.3

9.5 Linearity

Samples containing different amounts of analyte were serially diluted with *Sample Diluent*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

Sample	1	2	3	4	
Concentration (ng/mL)	951.0	3300.0	6905.0	8445.0	
Average Recovery (%)	92.9	90.09	96.7	97.5	
Range of Recovery (%)	from	88.3	86.5	86.2	89.9
	to	99.4	95.1	107.3	103.3

9.6 Interfering Substances

Hemoglobin (up to 4 mg/mL), bilirubin (up to 0.5 mg/mL) and triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

9.7 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence on the measurement of Linezolid in a sample.


9.8 High-Dose Hook Effect

"High-Dose Hook Effect" is not detected up to 50 000 ng/mL of Linezolid.

10 METHOD COMPARISON

A comparison of DRG:HYBRiD-XL Linezolid Test (HYE-6163) (y) and LC-MS measurements (x) using clinical samples gave the following correlation:

$$\begin{aligned} n &= 40 \\ r &= 0.98 \\ y &= -9.412 + 0.9798x \end{aligned}$$

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the laboratory quality assurance guidelines and applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different kit lots could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

11.4 Reporting of Serious Incident

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

For further information please refer to the User Manual of the DRG:HYBRiD-XL, analyzer-specific application sheets, product information and package inserts of all necessary components.

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11.4 Notification des incidents graves

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
1. Yue_J, et al.: Linezolid versus vancomycin for skin and soft tissue infections. Cochrane Database of Systematic Reviews 2016, Issue 1. Art. No.: CD008056.
2. Wunderkink et al. Linezolid in Methicillin-Resistant Staphylococcus aureus Nosocomial Pneumonia: A Randomized, Controlled Study. *Clinical Infectious Diseases*. 2012 ;54(5): 621–9.
3. Li et al. Linezolid vs. vancomycin in treatment of methicillin-resistant staphylococcus aureus infections: a meta-analysis. *European Review for Medical and Pharmacological Sciences*. 2017; 21: 3974-3979.
4. Almeida et al.: Ribosome-targeting antibiotics impair T cell effector function and ameliorate autoimmunity by blocking mitochondrial protein synthesis. *Immunity*. 2021, 54: 68–83.
5. Hashemian et al. Linezolid: a review of its properties, function, and use in critical care. *Drug Design, Development and Therapy*. 2018;12: 1759–1767.
6. Richards and Brink. Therapeutic drug monitoring: linezolid too? *Critical Care*. 2014; 18: 525.
7. Esteve-Pitarch E et al. OTAC: Optimization of Antibiotic Therapy in Critically ill patients. Using beta-lactam antibiotics by continuous infusion. *Farm Hosp*. 2019.; 43(5):151-157.
8. Mitchell M Levy MM et al. Outcomes of the Surviving Sepsis Campaign in intensive care units in the USA and Europe: a prospective cohort study. *The Lancet Infectious Diseases*. 2021; 21, 12: 919-924.
9. Noll, I. et al. Daten zur Antibiotikaresistenzlage in Deutschland. *Bundesgesundheitsbl*. 2012; 55: 1370–1376.
10. Friedrichs A, et al. Antibiotic Stewardship: Antibiotic stewardship: Structure and implementation. *Med Klin Intensivmed Notfmed*. 2021; 116(1): 81-92.
11. Hashemian SMR et al. Linezolid: a review of its properties, function, and use in critical care. *Drug Des Devel Ther*. 2018; 12: 1759–1767.
12. Li Y, Xu W. Efficacy and safety of linezolid compared with other treatments for skin and soft tissue infections: a meta-analysis. *Bioscience Report*. 2018; 38(1): 1-9.

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