Instructions for Use

Salivary Estradiol HS ELISA

IVD

REF SLV-6145

Σ 96

DRG

DRG Instruments GmbH, Germany
Frauenbergstraße 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

DRG

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos / Table des matières

1	INTENDED USE	2	1	FINALIDAD PREVISTA.....	22
2	PRINCIPLE OF THE TEST	2	2	FUNDAMENTO DEL ENSAYO	22
3	WARNINGS AND PRECAUTIONS	3	3	PRECAUCIONES	22
4	REAGENTS	4	4	COMPONENTES DEL KIT.....	23
5	SPECIMEN COLLECTION AND PREPARATION.....	5	5	MUESTRAS	24
6	ASSAY PROCEDURE	5	6	PROCEDIMIENTO DE ENSAYO	24
7	EXPECTED NORMAL VALUES	7	7	VALORES ESPERADOS	26
8	QUALITY CONTROL	7	8	CONTROL DE CALIDAD	26
9	PERFORMANCE CHARACTERISTICS.....	7	9	CARACTERÍSTICAS DEL ENSAYO.....	26
10	LIMITATIONS OF USE	9	10	LIMITACIONES DE USO	27
11	LEGAL ASPECTS.....	9	11	ASPECTOS LEGALES	27
1	ZWECKBESTIMMUNG	10	1	UTILISATION PRÉVUE	28
2	TESTPRINZIP	10	2	PRINCIPE DU TEST	28
3	VORSICHTSMAßNAHMEN	10	3	AVERTISSEMENTS ET PRECAUTIONS	28
4	BESTANDTEILE DES KITS	11	4	RÉACTIFS	29
5	PROBENVORBEREITUNG	12	5	PRÉLÈVEMENT ET PRÉPARATION DES SPÉCIMENS	30
6	TESTDURCHFÜHRUNG	12	6	RÉALISATION DU TEST	30
7	ERWARTETE WERTE	14	7	VALEURS NORMALES ATTENDUES	32
8	QUALITÄTSKONTROLLE	14	8	CONTROLE DE QUALITE	32
9	ASSAY-CHARAKTERISTIKA	14	9	CARACTERISTIQUES DES PERFORMANCES	32
10	GRENZEN DES TESTS	15	10	LIMITES D'UTILISATION.....	33
11	RECHTLICHE GRUNDLAGEN.....	15	11	ASPECTS JURIDIQUES.....	33
1	DESTINAZIONE D'USO	16	12	REFERENCES / LITERATURE	34
2	PRINCIPIO DEL TEST	16			
3	PRECAUZIONI	16			
4	COMPONENTI DEL KIT	17			
5	CAMPIONI	18			
6	ATTUAZIONE DEL TEST	18			
7	VALORI NORMALI	20			
8	CONTROLLO QUALITÀ	20			
9	CARATTERISTICHE DEL TEST.....	20			
10	LIMITAZIONE DEL TEST	21			
11	ASPETTI LEGALI	21			
				SYMBOLS USED	35

1 INTENDED USE

The **DRG Salivary Estradiol HS ELISA** is an enzyme immunoassay for the quantitative *in vitro* diagnostic measurement of active free estradiol in saliva.

Results may be used to assess fertility problems in women; to diagnose menopause, and to monitor hormone replacement therapy.

1.1 Summary and Explanation

Estradiol (1,3,5(10)-estratriene-3,17 β -diol; 17 β -estradiol; E2) is a C18 steroid hormone with a molecular weight of 272.4 Dalton. It is the most potent natural estrogen, produced mainly by granulosa cells of the female ovary and the placenta by the aromatization of androstenedione to estrone, followed by conversion of estrone to estradiol by 17 β -HSD. Estradiol is also synthesized in other tissues including testicles, adrenal gland, fat tissue, liver, breast, and brain (1-5).

In plasma, estradiol is largely bound to SHBG and albumin. Only a fraction of 2.21 % is free and biologically active, the percentage remaining constant throughout the menstrual cycle (6-8). Estradiol acts primarily as an agonist of the estrogen receptor (ER) subtypes ER α and ER β , nuclear steroid hormone receptors which trigger the appropriate response at the nuclear level in the target sites. These sites include follicles, uterus, breast, vagina, urethra, hypothalamus, pituitary and to a lesser extent liver and skin.

In the female, estradiol acts as a growth hormone for tissue of the reproductive organs. During the menstrual cycle, estradiol secretion follows a cyclic, biphasic pattern with the highest concentration found immediately prior to ovulation. This estradiol peak stimulates the hypothalamic-pituitary axis to secrete the gonadotropins follicle stimulating hormone (FSH) and luteinizing hormone (LH), which are essential for follicular maturation and ovulation. In the luteal phase, estradiol, in conjunction with progesterone, prepares the endometrium for implantation (9-15). During pregnancy, estradiol concentration increases due to placental production and high levels are sustained throughout pregnancy (16.)

Salivary steroid levels can reflect the circulating level of free steroid rather than total levels in serum, which is explained by the absence of circulating estradiol binding proteins. Therefore, saliva provides an excellent specimen for monitoring estradiol levels across the menstrual cycle or during hormone replacement therapy (17-19).

2 PRINCIPLE OF THE TEST

The DRG Salivary Estradiol HS ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **principle of competitive binding**.

The microtiter wells are coated with a polyclonal (rabbit) antibody directed towards antigenic sites of the estradiol molecule.

During incubation, the estradiol in the added sample competes with the added enzyme conjugate, which is estradiol conjugated to horseradish peroxidase, for binding to the coated antibody.

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of colour is inversely proportional to the concentration of the analyte in the sample.

A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with Stop Solution containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with anti-estradiol antibody (polyclonal).
2. **Standard (Standard 0 - 5)**, 6 vials, 1 mL each, ready to use;
Concentrations: 0 – 1 – 5 – 10 – 50 – 100 pg/mL
The standards are calibrated against the following reference material: Cerilliant E-060
Contain non-mercury preservative.
3. **Control Low & High**, 2 vials, 1 mL each, ready to use;
For control values and ranges please refer to vial label or Certificate of Analysis.
Contain non-mercury preservative.
4. **Enzyme Conjugate**, 1 vial, 14 mL, ready to use;
estradiol conjugated with horseradish peroxidase;
Contains non-mercury preservative.
5. **Substrate Solution**, 1 vial, 14 mL, ready to use;
Tetramethylbenzidine (TMB).
6. **Stop Solution**, 1 vial, 14 mL, ready to use;
Contains 0.5 M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.
7. **Wash Solution**, 1 vial, 30 mL (40X concentrated);
See "Reagent Preparation".

4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm) (e.g. the DRG Instruments Microtiter Plate Reader)
- Calibrated variable precision micropipettes
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

Wash Solution

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 1 week at room temperature.

4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Saliva can be used in this assay.

The sampling must be done in the morning **before** food intake, drinking or smoking. It is recommended to rinse mouth thoroughly with cold water 5 minutes prior to sampling.

If sampling must be done during the day, eating, drinking, smoking, chewing gums or brushing teeth should be avoided at least 2 hours before sampling.

Do not collect samples when oral diseases, inflammation or lesions exist (blood contamination).

Note: Samples containing sodium azide should not be used in the assay.

Please refer to chapters “Interfering Substances” and “Drug Interferences”.

5.1 Specimen Collection

Saliva samples should be collected using SALI-TUBES 100 (REF SLV-4158, available from DRG).

Other saliva sampling devices have not been tested and should be validated under the responsibility of the user.

5.2 Specimen Storage and Preparation

Fresh saliva samples

Immediately after arrival in the lab fresh saliva samples should be **frozen at least overnight at -20 °C.**

Each saliva sample has to be frozen, thawed, and centrifuged in order to separate the mucins by centrifugation.

Storage: immediately at -20 °C

Then samples must be thawed and centrifuged for 5 to 10 minutes at 10 000 g.

Thereafter, the clear supernatant must be transferred into a fresh tube.

Only this clear supernatant can be used as sample for the ELISA.

If a set of multiple samples has to be tested, the lab has to **mix aliquots of the supernatant** of the 5 single samples in a separate sampling device and perform the testing from this mixture.

Supernatant

Storage: up to 7 days at 2 °C to 8 °C

at least 2 months at -20 °C, in aliquots

The supernatant should be frozen only once.

Thawed supernatant should be inverted several times prior to testing!

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more analyte than the highest standard, the specimen can be diluted with 0.9% NaCl and re-assayed as described in “Assay Procedure”.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

a) dilution 1:2: 50 µL sample + 50 µL 0.9% NaCl (mix thoroughly)

b) dilution 1:10: 10 µL sample + 90 µL 0.9% NaCl (mix thoroughly)

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Optical density is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **100 µL** of each **Standard, Control** and **sample** with new disposable tips into appropriate wells.
3. Incubate for **60 minutes** at room temperature.
4. Dispense **100 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
5. Incubate for **60 minutes** at room temperature.
6. Rinse the wells **3 times** with **400 µL** diluted *Wash Solution* per well, if a plate washer is used.
- OR -
Briskly shake out the contents of the wells.
Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.
Important note:
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
7. Add **100 µL** of **Substrate Solution** to each well.
8. Incubate for **30 minutes** at room temperature.
9. Stop the enzymatic reaction by adding **100 µL** of **Stop Solution** to each well.
10. Determine the optical density of the solution in each well at **450 nm (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average optical density (OD) values for each set of standards, controls and patient samples.
2. Using graph paper, construct a standard curve by plotting the mean OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean OD value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4-Parameter Rodbard or 4-Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 100 pg/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Density (450 nm)
<i>Standard 0</i> (0 pg/mL)	1.99
<i>Standard 1</i> (1 pg/mL)	1.80
<i>Standard 2</i> (5 pg/mL)	1.45
<i>Standard 3</i> (10 pg/mL)	1.20
<i>Standard 4</i> (50 pg/mL)	0.48
<i>Standard 5</i> (100 pg/mL)	0.24

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with apparently healthy individuals, using the DRG Salivary Estradiol HS ELISA the following data were observed:

Population	n	Mean (pg/mL)	Median (pg/mL)	2.5 th - 97.5 th Percentile (pg/mL)	Range (min. - max.) (pg/mL)
Males	58	3.74	3.45	1.26 – 9.12	0.74 – 11.46
Females					
Pre-menopausal					
Follicular phase	26	5.85	3.45	1.32 – 29.35	1.08 – 46.42
Luteal phase	26	4.28	3.12	1.46 – 14.97	1.19 – 22.94
Post-menopausal	26	5.56	3.03	1.21 – 28.43	1.02 – 41.48
Stages of pregnancy					
1. Trimester	24	8.27	5.08	2.33 – 25.29	1.90 – 28.98
2. Trimester	24	44.16	39.52	7.29 – 97.68	7.20 – 97.86
3. Trimester	24	75.35	75.69	23.25 – 126.49	21.38 – 149.68

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0.366 pg/mL – 100 pg/mL.

9.2 Specificity of Antibodies (Cross-Reactivity)

The following substances were tested for cross-reactivity of the assay:

Substance	Added Conc. (pg/mL)	Cross-reactivity (%)
Androstenedione	1000	0.03
DHEA	1440	0.00
Testosterone	1000	0.00
Cortisol	80000	0.02
17-OHP	1000	0.03
Estrone	300	1.77
Progesterone	1000	0.00
DHT	1000	0.00
Estriol	1000	0.61
17-a-Ethinylestradiol	20000	0.47

9.3 Sensitivity

The analytical sensitivity of the DRG ELISA was calculated by subtracting 2 standard deviations from the mean of 20 replicate analyses of *Standard 0* and was found to be 0.195 pg/mL.

The Limit of Blank (LoB) is 0.145 pg/mL.

The Limit of Detection (LoD) is 0.366 pg/mL.

The Limit of Quantification (LoQ) is 0.736 pg/mL.

9.4 Reproducibility

9.4.1 Intra-Assay

The within-assay variability was determined by measuring each sample 10 times per run (n = 10):

Sample	n	Mean (pg/mL)	CV (%)
1	10	2.50	9.4
2	10	4.95	9.4
3	10	26.45	6.3
4	10	92.54	1.4

9.4.2 Inter-Assay

The between-assay variability was determined by measuring each sample 10 times per run for 3 days (n = 30):

Sample	n	Mean (pg/mL)	CV (%)
1	30	2.51	14.9
2	30	5.49	12.3
3	30	29.48	9.4
4	30	93.20	1.7

9.4.3 Inter-Lot

The inter-assay (between-lots) variation was determined by measuring each sample 6 times with 3 different kit lots:

Sample	n	Mean (pg/mL)	CV (%)
1	18	2.67	13.5
2	18	6.82	6.2
3	18	26.35	5.9
4	18	85.69	5.5

9.5 Recovery

Samples have been spiked by adding estradiol solutions with known concentrations.

The recovery (%) was calculated by multiplying the ratio of measured and expected values with 100.

	Sample 1	Sample 2	Sample 3	Sample 4	
Concentration (pg/mL)	6.34	16.63	22.96	71.01	
Average Recovery (%)	100.8	87.4	88.8	90.0	
Range of Recovery (%)	from	98.8	86.8	87.4	88.9
	to	104.3	88.3	90.4	91.1

9.6 Linearity

Samples were measured undiluted and in serial dilutions with 0.9 % NaCl. The recovery (%) was calculated by multiplying the ratio of expected and measured values with 100.

	Sample 1	Sample 2	Sample 3	Sample 4	
Concentration (pg/mL)	32.02	63.24	65.23	80.89	
Average Recovery (%)	97.5	97.7	104.5	101.0	
Range of Recovery (%)	from	89.2	96.4	95.2	91.1
	to	109.9	100.3	109.8	111.0

9.7 Comparison Studies

A comparison of DRG Salivary Estradiol HS ELISA (SLV-6145) (y) and the reference method LC-MS/MS (x) using clinical samples gave the following correlation:

$$n = 9$$

$$r = 0.980$$

$$y = 0.775x + 7.536$$

A comparison of DRG Salivary Estradiol HS ELISA (SLV-6145) (y) and the IBL 17 β -Estradiol Saliva ELISA (REF 30121045) (x) using clinical samples gave the following correlation:

$$n = 48$$

$$r = 0.992$$

$$y = 1.691x - 3.462$$

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

The following substances will increase the estradiol values by more than 20 %: coffee, ethanol, coffee with milk, toothpaste, caffeine, chewing gum, nicotine.

In contrast, hemoglobin will decrease the estradiol values by more than 20%.

10.2 Drug Interferences

The Salivary Estradiol HS ELISA should not be used for patients being treated with the drugs Fulvestrant (Faslodex[®]) and Mifepristone (Mifegyne[®]), since this could lead to falsely elevated results.

Patients who are under treatment with L-Thyroxin may also show falsely elevated test results. To avoid interference with high dose L-Thyroxin treatment, the sampling must be done either directly in the morning before taking L-Thyroxin or at least 12 hours after taking L-Thyroxin.

10.3 High-Dose-Hook Effect

A High-Dose-Hook Effect is not known for competitive assays.

11 LEGAL ASPECTS

11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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

12 REFERENCES / LITERATURE

1. Tsang BK et al. (1980) Steroid biosynthesis by isolated human ovarian follicular cells in vitro. *J Clin Endocrinol Metab* 1980, 51, 1407 – 11.
2. Gore-Langton RE et al. Follicular steroidogenesis and its control. In: *The physiology of Reproduction*. 1988, Ed.: Knobil et al., 331-385, Raven press, New York.
3. Hall, PF. Steroid synthesis: Organization and Regulation. In: *The Physiology of Reproduction*. 1988, Ed.: Knobil et al., 975 – 988, Raven press, New York
4. Saldanha CJ, Remage-Healey L, and Schlinger BA. "Synaptocrine signaling: steroid synthesis and action at the synapse." *Endocrine reviews* 32.4, 2011, 532–549.
5. Wu CH et al. Free and protein-bound plasma estradiol-17 beta during the menstrual cycle". *The Journal of Clinical Endocrinology and Metabolism*. 1976, 43 (2), 436–45.
6. Siiteri PK et al. The serum transport of steroid hormones. *Rec Progr Horm Res* 1982, 38, 457 – 510
7. Martin B. et al. Binding of steroids by proteins in follicular fluid of the human ovary. *J Clin Endocrinol Metab* 1982, 35, 443 – 447.
8. Celec P et al. Salivary Sex Hormones during the Menstrual Cycle. *Endocrine Journal* 2009, 56 (3), 521-523
9. Gandara BK, Leresche L, and Mancl L. Patterns of Salivary Estradiol and Progesterone across the Menstrual Cycle *Ann N Y Acad Sci*. 2007, 1098, 446–450.
10. Baird DT Ovarian steroid secretion and metabolism in women. In: *The Endocrine function of the human ovary*. Eds.: James V.H.T., Serio M. and Guisti G., 1976, 125 – 33, Academic press, New York
11. McNtty KP et al. Concentration of estrogens and androgens in human ovarian venous plasma and follicular fluid throughout the menstrual cycle. *J. Endocrinol*. 1976, 71, 77 –85.
12. Gruhn JG, Kazer RR. *Hormonal Regulation of the Menstrual Cycle: The Evolution of Concepts*. Springer Science & Business Media. 2013, 69–73. ISBN 978-1-4899-3496-3.
13. March CM et al. Roles of estradiol and progesterone in eliciting mid-cycle luteinizing hormone and follicle stimulating hormone surges. *J. Clin. Endocrinol. Metab*. 1979, 49, 507 – 12.
14. Stricker R et al. Establishment of detailed reference values for luteinizing hormone, follicle stimulating hormone, estradiol, and progesterone during different phases of the menstrual cycle on the Abbott ARCHITECT analyzer". *Clinical Chemistry and Laboratory Medicine*. 2006, 44 (7), 883–7.
15. Simpson ER and McDonald PC. *Endocrinology of pregnancy*. In: *Textbook of Endocrinology*, Ed.: Williams R.H. 1981, 412 – 22, Saunders Company, Philadelphia.
16. Lipson, SF and Ellison, PT. Comparison of salivary steroid profiles in naturally occurring conception and non-conception cycles. *Hum Reprod*, 1996, 11(10), 2090-96.
17. Choe, JK, Khan-Dawood, FS, and Dawood, MY. Progesterone and estradiol in saliva and plasma during the menstrual cycle. *Am J Obstet Gynecol*, 1983, 147(5), 557-62.
18. Lewis JG. Steroid Analysis in Saliva: An overview. *Clin Biochem Rev* 2006, 27 August 139-146.
19. Gameiro, CM, Romão, F, and Castelo-Branco, C. Menopause and aging: Changes in the immune system –A review. *Maturitas*, 2010, 67(4), 316-20.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico in vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	No. di Cat.	No de catálogo	Référence
	Batch code *	Fertigungslosnummer, Charge *	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution *	Achtung *	Attenzione	Precaución	Attention
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement
<i>Volume/No.</i>	Volume / No.	Volumen / Anzahl	Volume / Quantità	Volumen / Número	Volume / Quantité
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Micropozzetti	Placas multipocillo	Microplaques
<i>Antiserum</i>	Antiserum	Antiserum	Antisiero	Antisuero	Antisérum
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complejo enzimático	Complexe enzymatique
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d'arresto	Solución de parada	Solution d'arrêt
<i>Zero Standard</i>	Zero Standard	Nullstandard	Standard zero	Estándar cero	Zero Standard
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)	1N NaOH	1N NaOH
<i>1 N HCl</i>	1 N HCl	1 N HCl		1 N HCl	1N HCl
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungs-medium	Diluyente dei campioni	Solución para dilución de la muestra	Solution pour dilution de l'échantillon
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungs-medium	Diluyente del tracciante	Solución para dilución del conjugado	Solution pour dilution du conjugué